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EN

1. MODELS	9
2. INTENDED USE	9
2.1 INTENDED USE AND CLINICAL BENEFITS	9
2.2 TARGET PATIENTS	9
2.3 PATIENT SELECTION CRITERIA	9
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	9
2.5 USERS AND INSTALLERS	9
3. REFERENCE STANDARDS	9
4. INTRODUCTION	9
4.1 DEVICE LABELLING AND TRACEABILITY	9
4.2 SIMBOLI	9
4.3 WARRANTY AND SERVICE	9
5. WARNINGS/DANGERS	9
6. SPECIFIC WARNINGS	10
7. RESIDUAL RISK	10
8. TECHNICAL DATA AND COMPONENTS	10
9. PROPER USE	10
9.1 CHOOSING THE SIZE	10
9.2 APPLYING THE EXTRICATION COLLAR	10
9.3 NEC UP ADJUSTMENT	10
10. CLEANING AND MAINTENANCE	10
10.1 CLEANING	10
10.2 MAINTENANCE	10
10.3 LIFE SPAN	11
11. TROUBLESHOOTING TABLE	11
12. ACCESSORIES	11
13. SPARE PARTS	11
14. DISPOSAL	11

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- NEC UP
- SPECTAR
- DALLAS
- NEC YORK
- YES

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Cervical collars are devices to be used to help maintain neutral alignment of the cervical spine. Depending on the model, they may or may not be equipped with a tracheotomy hole and be mono- or bi-valve.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group. The range of sizes for each model allows it to be utilised by an extensive group of patients, as long as the chosen size achieves the intended purpose of the device.

2.3 PATIENT SELECTION CRITERIA

The expected patients are typically those for whom cervical spine immobilisation is necessary. The use of the cervical collar is contraindicated in cases of embedded/penetrating bodies.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

2.5 USERS AND INSTALLERS

The intended users are rescue workers with in-depth knowledge related to the immobilisation and handling of individuals with neck trauma. The user must be able to evaluate the most suitable type of device for use on the specific patient. These devices are not intended for lay people.

Operators must be able to provide the necessary patient care.

The product must be used only by personnel trained in the use of this product and not on other similar products.

3. REFERENCE STANDARDS




REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT).

4.2 SIMBOLI

Simbolo	Significato	Simbolo	Significato
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the instructions for use
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Disposable

4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

Warranty and service conditions are available on the website <http://support.spencer.it>

5. WARNINGS/DANGERS

Product features

Use of the product for any purpose other than that described in the User Manual is prohibited.

- The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to + 50°C.
- Storage temperature: -10°C to +60°C.

General warnings for medical devices

- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.

With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

IT

EN

DE

FR

ES

PT

EL

RO

CS

SL

ET

6. SPECIFIC WARNINGS

To use the cervical collars, you must also have read, understood and carefully follow all the instructions in the user manual.

-  Perform immobilisation simulations with dummies before putting the device into service.
-  Improper use may result in injury or permanent disability. Always follow the procedures and protocols approved by the relevant Emergency Medical Service.
- To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- **Re-use of disposable devices may cause cross infections.**
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- Qualified personnel must be present during use of the device.
- At least two rescuers must be present for correct application of the device.
- Operators must be able to assess the patient's injuries and decide whether the use of the device is appropriate.
- Before using the collar, assess the need for its application according to the clinical condition of the patient and make sure that the most suitable type of device is used.
- Operators must be able to assess the most suitable type and size of cervical collar for use, taking into account the type of patient and the characteristics of the injuries.
- The device must not be exposed, much less come into contact with thermal sources of combustion or flammable agents.
- Always check the conditions of all parts before use.

On the Spectar model, the chin guard must be secured by inserting the red chin guard clip into the corresponding hole at the end. Do not pull on the clip but press it into the fastening slot. Improper use may cause the clip to break.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in this user manual, have been identified.

8. TECHNICAL DATA AND COMPONENTS

	SPECTAR	NEC UP	NEC YORK	YES	DALLAS
					
PARTS	DESCRIPTION/MATERIALS				
Main body	Constitutes the primary support element. Made of rigid polyethylene in Nec Up, Spectar, Nec York. Made of polyethylene foam in the remaining models.				
Lining	If present, also called padding; improves patient comfort. Made of polyethylene foam in Nec Up, Spectar, Nec York and Yes.				
Chin guard	If present, supports the patient's chin. Made from the main body and the lining				
Fastening strap	Closes around the patient's neck. Made of nylon fabric and is adhered by means of strap inserts.				

9. PROPER USE

Operators performing immobilisation must be able to select the most suitable type of cervical collar and the most appropriate collar size. Regardless of the type of collar used, make sure that there are no conditions that are not compatible with the use of this device (i.e. penetrating bodies, neck injuries, difficulty in maintaining airway patency).

9.1 CHOOSING THE SIZE

- Measure the distance between the insertion point of the clavicle and the lower edge of the chin.
- Identify the most suitable collar size by comparing the measurement taken above with the distance between the chin strap in its intended position of application and the part of the collar that will rest on the patient's shoulder.

9.2 APPLYING THE EXTRICATION COLLAR

- One rescuer must always maintain neutral alignment of the cervical spine while the second rescuer applies the collar.
- In the case of bi-valve collars, position the front valve, sliding it from the chest towards the chin.
- The rear valve should be slid under the nape of the neck, taking care to maintain neutral alignment at all times.
- Attach the two valves with the straps, tightening them to ensure correct immobilisation.
- In the case of mono-valve collars, apply by sliding the back of the collar behind the nape of the neck and then close the opposite end on the front side. Secure it with the strap bands to ensure correct immobilisation.

9.3 NEC UP ADJUSTMENT

Two sizes of Nec Up collars are available. After assessing the size as described in the previous paragraph, bring the chin strap up to the appropriate height. The selected size is shown by a red indicator near the chin guard. Then lock the chin guard, pressing on both locking tabs identified by the word "Lock".

10. CLEANING AND MAINTENANCE

10.1 CLEANING

Cervical collars are disposable devices and therefore there are no recurring cleaning procedures. If the devices become soiled before use, or between uses for the same patient, clean them with lukewarm water and mild soap. Allow the devices to dry thoroughly before use.

10.2 MAINTENANCE

The device does not require a routine maintenance program, but checks must be made before use to verify:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Fulfilment of the requirements of the user manual in section Warnings and Specific Warnings.

No periodic overhaul is foreseen for the device.

10.3 LIFE SPAN

If unused, the device has a life span of 5 years from the date of purchase.

11. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The straps do not adhere correctly.	The straps are dirty.	Remove any residue from the straps.
Broken fixing/locking clips	Improper use	Take the device out of service immediately and replace it with a similar one.
The collar does not provide adequate immobilisation.	The collar is not the correct size.	Check that the chosen size is suitable for the patient.
	Different valve sizes have been coupled.	Check that both collar shells have the same size indication and strap colouring.
The patient complains of discomfort from the collar.	The collar has not been applied correctly or the patient's clinical condition is not suitable for application of the collar	Check that the collar has been applied according to the instructions on the collar.

12. ACCESSORIES

QC10120	COLLAR BAG 3 - BAG WITH DIVIDERS MONO-VALVE COLLARS	QC19010	COLLAR BAG 1 - BAG FOR MONO-VALVE COLLARS
JM30028	SPECTAR FOAM - STERNUM PAD FOR COLLAR		

13. SPARE PARTS

No spare parts are available.

14. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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© Copyright Spenser Italia S.r.l.