BLUE SPLINT**FIXO SPLINT****URRA'****AIR SPLINT****RES-Q-SPLINT /
RES-Q-SPLINT PLUS****EMSOFT / EMSOFT PRO****Manuale d'uso e Manutenzione
STECCOBENDE E ACCESSORI**

IT

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EN

1.	MODELS	8
2.	INTENDED USE	8
2.1	INTENDED USE AND CLINICAL BENEFITS	8
2.2	TARGET PATIENTS	8
2.3	PATIENT SELECTION CRITERIA	8
2.4	CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	8
2.5	USERS AND INSTALLERS	8
3.	REFERENCE STANDARDS	8
4.	INTRODUCTION	8
4.1	DEVICE LABELLING AND TRACEABILITY	8
4.2	SIMBOLI	8
4.3	WARRANTY AND SERVICE	8
5.	WARNINGS/DANGERS	8
6.	SPECIFIC WARNINGS	9
7.	RESIDUAL RISK	9
8.	TECHNICAL DATA AND COMPONENTS	9
9.	PROPER USE	10
10.	CLEANING AND MAINTENANCE	10
10.1	CLEANING	10
10.2	MAINTENANCE	10
10.3	LIFE SPAN	10
11.	TROUBLESHOOTING TABLE	10
12.	ACCESSORIES	11
13.	SPARE PARTS	11
14.	DISPOSAL	11

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- BLUE SPLINT
- RES-Q-SPLINT
- FIXO SPLINT
- RES-Q-SPLINT PLUS
- URRRA'
- EMSOFT
- AIR SPLINT
- EMSOFT PRO

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group.

The design of the product allows it to be applied to any subject as long as the chosen size achieves the intended purpose of the device.

2.3 PATIENT SELECTION CRITERIA

The intended patients are typically those for whom it is necessary to immobilise a fractured or injured limb.

The use of the splints is not recommended in cases of embedded/penetrating bodies.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

2.5 USERS AND INSTALLERS

The intended users are rescue workers with in-depth knowledge related to the immobilisation and handling of individuals with fractures or suspected fractures.

These devices are not intended for lay people.

Operators must be able to provide the necessary patient care.

The product must be used only by personnel trained in the use of this product and not on other similar products.

3. REFERENCE STANDARDS

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

If the assigned Lot/SN cannot be traced, the device must be reconditioned, provided only under the responsibility of the manufacturer.

4.2 SIMBOLI

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death
	Medical device		See the user manual
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner (only for USA Market)

4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

Warranty and service conditions are available on the website <http://support.spencer.it>.

5. WARNINGS/DANGERS

Product features

Use of the product for any purpose other than that described in the User Manual is prohibited.

- The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to +50°C.
- Storage temperature: -10°C to +60°C.

General warnings for medical devices

- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.

With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

6. SPECIFIC WARNINGS


To use the splints, you must also have read, understood and carefully follow all the instructions in the user manual.

- ⚠ Perform immobilisation simulations with dummies before putting the device into service.
- ⚠ Improper use may result in injury or permanent disability. Always follow the procedures and protocols approved by the relevant Emergency Medical Service.
- To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- Qualified personnel must be present during use of the device.
- The splint must be applied by at least two trained responders with good dexterity and common sense.
- Before shaping the splint, evaluate the application according to the visibility of any wounds or the fracture site..
- The device must not be exposed, much less come into contact with thermal sources of combustion or flammable agents.
- Only use the device as described in this manual.
- Always check the conditions of all parts before use.
- Operators must be able to assess the patient's injuries and decide whether the use of the splint is appropriate.
- Professionals must be able to assess the most appropriate type and size of splint to use, considering the type of patient and the characteristics of the fracture.
- If pumps not supplied by Spencer are used, it is necessary to ensure that the pump being used is suitable for the correct immobilisation of the limb. The use of unsuitable pumps may prevent the splint from achieving its intended results. The suitability of this part must be checked before using the device.

7. RESIDUAL RISK

No residual risks are identified, i.e. risks that could arise despite compliance with all the warnings in this user manual.

8. TECHNICAL DATA AND COMPONENTS

BLUE SPLINT (Semi-rigid)		EMSOFT (Vacuum)	
	Flexible core 1 Made of aluminium, it allows the splint to be shaped as required.		Splint body 1 Made of PVC, it is filled with small polystyrene balls
	Splint body 2 Made of polyurethane, it wraps around the limb.		Valve 2 Made of PVC, allows the splint to be depressurised and released once use is complete.
	Closing phase Equipped with strap elements, they allow the splint to be applied to the patient's limb.		Closing phase 3 Equipped with strap elements, they allow the splint to be applied to the patient's limb.
FIXO SPLINT (Rigid)		AIR SPLINT (Inflatable)	
	Splint body 1 Made with a PVC-coated polyurea core		Splint valve 1 Made of PVC, allows the splint to be inflated
	Closing phase 2 Equipped with strap elements, they allow the splint to be applied to the patient's limb.		Splint body 2 Made of PVC, it stiffens around the patient's limb as it inflates.
Res-q-splint/Res-q-splint pro (Vacuum)		URrà (Semi-rigid)	
	Splint body 1 Made of PVC, it is filled with small polystyrene balls		Splint body 1 Made with an aluminium core coated with XPE
	Valve 2 Made of aluminium, allows the splint to be depressurised and released once use is complete.		
	Closing phase 3 Equipped with strap elements, they allow the splint to be applied to the patient's limb.		

Characteristics	Blue splint				
	Leg	Arm	Forearm	Wrist	Elbow-ankle
Size (mm)	580x450	590x280	390x300	310x240	540x280
Weight (grams)	360	300	206	145	273

Characteristics	Fixo splint		
	XS	S	L
Size (mm)	120x220x310	120x230x500	150x230x730
Weight (grams)	410	640	870

Characteristics	URRA'		
	S	M	L
Size (mm)	22x11x0,5	46x11x0,5	91x11x0,5
Weight (grams)	43	82	145

IT
EN
DE
FR
ES
PT
EL
HU
RO
SK

Res-q-splint			
Characteristics	S	M	L
Size (mm)	560x460	490X380X700	750x540x900
Weight (grams)	295 ± 50	630 ± 50	1050 ± 50

Res-q-Splint Plus			
Characteristics	S	M	L
Size (mm)	735x640	1100x695 ±	1300x720 ±
Weight (grams)	600 ± 10	850 ± 10	1050 ± 10

EMSoft			
Characteristics	S	M	L
Size (mm)	750x600	900x700	1100x800
Weight (grams)	585	830	1130

Sizes subject to ± 10 mm tolerances.

9. PROPER USE

Operators performing immobilisation must be able to select the most suitable type of splint and the most appropriate splint size. Before and after the application of the splint, it is necessary to check the circulation of the limb to be immobilised.

Applying rigid splints (Blue Splint , Fixo Splint, Urrá)

- Select the most suitable size according to the characteristics of the limb to be immobilised.
- Shape the splint to give it the most suitable shape (only Blue Splint).
- Apply the splint gently to the limb.
- Close the strap bands, ensuring that the tissues of the fractured limb are not disturbed. Urrá must be closed with gauze bandages.

Application of inflatable splints (Air splint)

- Select the most suitable size according to the characteristics of the limb to be immobilised.
- Unzip the splint
- Gently raise the limb and slide the open splint under the limb.
- Zip up the splint and inflate it with the valve nozzle, ensuring that the compression generated by the splint does not harm the tissue. Preferably use the QMX109 pump for inflation.
- Close the valve opening
- Check pulse and paleness, movement and sensitivity below the area of application as compression may cause reduced tissue circulation.

Application of vacuum splints (Res-Q-Splint, Res-Q-Splint Plus, EMSOft)

- Select the most suitable size according to the characteristics of the limb to be immobilised.
- Lay the splint on the ground and distribute the filling material evenly.
- Gently raise the limb and slide the open splint under the limb.
- Gently wrap the splint around the limb and close the strap bands
- Connect the suction pump to the valve, and begin depressurising while keeping the device adherent to the limb. With EMSOft, use the conical fitting supplied with the splint to connect the splint to the pump. Once the required rigidity has been achieved, close the valve and adjust the length of the closing straps if necessary.

After use, return the device to atmospheric pressure by opening the valve.

With EMSOft, close the protective cap

Res-Q-Splint and EMSOft splints are compatible with QMX109, QMX120, QMX125, QMX130 pumps.

EMSOft, requires the use of the conical fitting supplied with the splint.

10. CLEANING AND MAINTENANCE

10.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. Never use solvents or stain removers.

Rinse thoroughly with lukewarm water, making sure you have removed all traces of detergent, which may deteriorate or compromise conditions and durability. The use of high-pressure water should be avoided.

Allow to dry thoroughly before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources. If blood is visible, oxidise it before washing with water.

Blue splint can be machine washed at 40°C, after removal of the metal core.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

10.2 MAINTENANCE

The device does not require a routine maintenance program, but checks must be made to verify:

- General functionality of the device
 - Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
 - Fulfilment of the requirements of the user manual in section Warnings and Specific Warnings.
 - Fulfilment of the requirements of the manual in section on Proper use.
- No periodic overhaul is foreseen for the device.

10.3 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase.

11. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The pump cannot be connected to the splint	Loss of a fitting or improper use	Check that the pump type is suitable for the device in use
Fabric or belts are torn	Normal wear or improper use.	Immediately remove the device from service and contact the manufacturer

It is not possible to attach one or more buckles or straps.	Possible incorrect buckle orientation or buckle failure.	Check the correct orientation of the buckle and its conditions. If damaged, immediately remove the device from service and contact the manufacturer
	Dirty straps	Check the cleanliness of the strap elements

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

12. ACCESSORIES

QM22199	QMX 199 - VACUUM SYSTEM REPAIR KIT	JM90004	URRA' SPLINT CASE
JM80008	BLUE SPLINT - SPLINT BAG WITH SILKSCREEN PRINTING	QC70100	FIXO SPLINT - YELLOW PVC CARRY BAG
QS24201	ECO AIR SPLINT - NYLON CARRY BAG WITHOUT POCKETS	QM22541	RES-Q-SPLINT - YELLOW WELDED PVC CARRY BAG
QM22704	EMSOFT CONNECTOR	QS24102	AIR SPLINT - NYLON CARRY BAG WITH POCKETS

13. SPARE PARTS

RIQM004	RES Q SPLINT/RQS PLUS - SPARE VALVE	RIQM009	QMX120 – ALUMINIUM VACUUM PUMP
RIQM006	EMSOFT - SPLINT VALVE SPARE CAPS (5 pieces)	RIQM010	QMX125 – ALUMINIUM DOUBLE EFFECT PUMP
RIQM008	QMX109 – DOUBLE EFFECT ABS PUMP	RIQM011	QMX130 – VACUUM SPLINT PUMP

14. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

EL

HU

RO

SK

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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