1. GENERAL INFORMATION

1.1 AIM AND CONTENTS

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

1.2 CONSERVATION OF THE INSTRUCTION MANUAL

The instruction and maintenance manual must be kept with the product, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

1.3 SYMBOLS USED (See chart 1.3 on reverse)

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For any information regarding the use, maintenance and installation, please contact the Spencer Customer Care Service on tel. 0039 0521 541111, fax 0039 0521 541222, e-mail info@spencer.it or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY. In order to facilitate the assistance service, please always indicate or communicate the serial number (SN) or lot number (LOT) shown on the label applied on the box or on the device.

1.5 DEMOUNTING

When the devices are no more suitable for being used, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations about demolition.

1.6 LABELLING

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the Manufacturer, the product, CE mark, lot number (LOT). It must never be removed or covered. The terminal unit and the inlets for gas therapy are individually branded with the following information: "CE" mark followed by number "0123" which identifies the European notified body who approved the "CE" certification of the device. The manufacturer guarantees the traceability of the devices for gas therapy fill each direct client who on behalf of the manufacturer is obliged to guarantee the traceability to the end user.

2. WARNINGS

2.1 GENERAL WARNINGS

- * The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- * Before carrying out any kind of operation on the appliance, the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- * If the instructions belong to another device and not the device received, inform the Manufacturer immediately and avoid use of the device.
- * In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- * Regularly check the appliance. Carry out the prescribed maintenance in order to keep the appliance in good condition and to guarantee correct functioning and a long life.
- * In the case of any abnormalities or damage to the appliance, which could jeopardize the functioning and the safety, the appliance must be immediately removed from service.
- * Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- * The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- * Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- * Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments.
- * Attention: Laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
- * Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- * Use of the device in anyway other than described in this manual is forbidden.
- * Handle with care.
- * With reference to the D. Lgs. 24th February 1997, n. 46 amended by D. Lgs. 25/01/2010, n. 37 - Acknowledgement of Directive 93/42/CEE and 2007/47/CE, we remind both public and private operators that they are obliged to report any accident that involves any medical device to the Ministry of Health and to the Manufacturer as specified and within time given by the European regulations.

In addition, both public and private operators are obliged to inform the Manufacturer of any measures that should be adopted to make the steps necessary to guarantee the safety and the health of the patients and the users o any medical device.

2.2 SPECIFIC WARNINGS

- * Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the Competent Authorities and/or Manufacturer if requested.
- * If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- * Before each use of device the perfect operating state of the device must be checked as specified in the Instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the Manufacturer must be contacted.
- * When the device is being used, the assistance of qualified staff must be guaranteed.
- * Never leave an unassisted patient. The presence of at least one operator is essential at all times when the medical device is in use.
- * Before using the device, check that all connected devices are compatible, especially if they are supplied by other manufacturers.
- * Once the installation is completed, check the tightness, by making a conformity check to norm EN ISO 7396-1. No leak is admitted.
- * Verify that inlet and outlet are projected for supplying the same gas and with reference to the same norm.
- * Do not use the device with different gas or pressure from the ones for which it has been designed.
- * Do not bump the inlet, because its functionality could be damaged.
- * The device must be exposed to or come into contact with any source of combustion or inflammable agents.
- * Store in a cool, dry, dark place and do not expose to direct sun.
- * Do not store the device underneath any heavy objects which could cause structural damage.
- * Store and transport device in its original packaging.
- * Position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.

2.3 CONTRAINDICATIONS AND SIDE EFFECTS

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

3. DESCRIPTION OF PRODUCT

3.1 INTENDED USE

The Spencer rapid coupling attachment for gas therapy is an easy to use and safe device. It has been designed for distribution systems for medical gas, where the operator can make many connexions and disconnexions to use different gas as a source for anaesthesia, ventilation and medical appliances.

The device does not produce any noise (pressure acoustic level lower than 70 dBA) and dangerous vibrations for people.

The Spencer rapid coupling attachment for gas therapy is a Class II b medical device, "CE" branded and manufactured in conformity with the 93/42/CEE and the UNI 9507 normative regarding systems for distribution of non inflammable gas for medical use. End connections and UNI EN ISO 9170-1 Systems for the distribution of medical gas. End connectors for compressed and vacuumed medical gas.

- * is not for single use, nor supplied as sterile
- * can't contain drugs, biological or animal membrane, or human blood
- * does not contain phthalates, classified as a carcinogen, mutagen and toxic for reproduction agent

3.2 MAIN COMPONENTS

The main components of the terminal unit are the following:

- Base block
- Outlet

The main components of the inlet are the following:

- Inlet

3.3 MODELS (See chart 3.3 on reverse)

3.4 TECHNICAL DATA (See chart 3.4 on reverse)

3.5 ENVIRONMENTAL CONDITIONS

Transport and storage temperature: -15 + +50 °C
Operating temperature as for EN ISO 9170-1: -20 + +60 °C
Humidity: 10 + 90%
Atmospheric pressure: 700 + 1060 hPa

4. OPERATING INSTRUCTIONS

4.1 TRANSPORT AND STORAGE

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself. Keep the original packaging for use in case of any further transport. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

4.2 PREPARATION

On receipt of the product:

- * Remove the packaging and display the material so that all components are visible.
 - * Check that all the components/pieces on the accompanying list are present.
- The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage.

4.3 INSTALLATION Δ

The installation of the terminal unit and inlet for gas therapy should be carried out only by adequately trained personnel who must strictly abide to the manufacturer's instructions. It is absolutely forbidden to in any way alter any component during installation. Materials used during installation and any lubricants required, must be compatible for use with the designated gas. All precautions should be taken to guarantee and preserve the cleanliness of the components and if required the sterile specifications should be observed.

Before use it is essential to have correct connection to earth as prescribed by the norm.

Once the installation is completed, check the tightness, by making a conformity check to norm EN ISO 7396-1. No leak is admitted.

4.3.1 TERMINAL UNIT

Proceed as follows:

- Connect the terminal unit to the tube using the nozzle and with an adequate single use cable tie, in conformity to norm EN ISO 5359
- Check that the fixing is air tight: leaks are not permitted. Then make checks according to norm EN ISO 7396-1
- Fix the inlet to the panel

4.3.2 INLET

Proceed as follows:

- Fix the terminal nozzle to the tube using and with an adequate single use cable tie, in conformity to norm EN ISO 5359
- If required insert the indicator of the destined gas
- Check that the fixing is air tight: leaks are not permitted. Then make checks according to norm EN ISO 7396-1

4.4 FUNCTIONING CONNECTION

Insert the coupling in the relative inlet using the energy necessary to contrast the pressure until the complete locking of the plastic nut ring.

Use only inlets conform to UNI 9507. Do not swap components destined or used for other gases. Δ

UNHOOKING OF INLET

Keep the inlet to avoid a fast disconnection.

Press the nut ring and unhook from the inlet.

Before unhooking check that the ray of protection of the inlet does not interfere with vital and or functional parts of the operator or with other people exposed to it's use.

Do not bump the inlet, because its functionality could be damaged. Δ

4.5 SAFETY

The user should organise training for personnel regarding the risks of injuries, of safety devices for the operator and of the general safety rules as stipulated in current laws.

The operator responsible for use of the Rapid Coupling Attachment for gas therapy, (both frequent and random users), should have adequate operating knowledge and capability such as to ensure safety for both himself and for third parties and be able to keep the device from possible damage.

The use of the terminal unit and inlet for gas therapy is severely forbidden to those who have not read and studied carefully this manual; by those who are not competent; or by those who are not in a good state of health.

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Before using the device, check that all components are in perfect condition and that the device is correctly connected to the tubes. Before use practise should be made on the controls and their functions. The use of the device is forbidden in case of breakdown or incorrect functioning of the safety systems. If the rapid coupling attachment for gas therapy is out of order or undesign maintenance, disconnect the source and attach a clear sign indicating "OUT OF USE DO NOT USE". If adequate illumination is not available do not use the terminal unit and inlet for gas therapy. When use has terminated, disconnect from source. Do not use the terminal unit and inlet for gas therapy with the same gas at different pressure rates and/or concentration to those for which it has been designed. The medical device has been designed in conformity to UNI 9507 and UNI EN ISO 9170-1 norms and it has been made in a way that it excludes the possibility of an exchange between different gases.

5. MAINTENANCE AND CLEANING

5.1 CLEANING

The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures. For cleaning the device, use only distilled water. Wash using warm water and neutral detergent, do not use stain removers or solvents. In case of disinfecting procedures, make sure that the products used are not corrosive and are adapted for use on each material of the various components of the device. Rinse carefully with warm water to remove all traces of detergents. Failure to do so could compromise the product and its life span. The device must be left to dry completely before storage. To dry the product after washing, or if used in a humid atmosphere, do not use direct sources of heat or flame.

5.2 MAINTENANCE

5.2.1 Precautionary Maintenance

The person who carries out the precautionary maintenance of the appliance has to guarantee the following basic requirements:

- * Technical knowledge of the appliance and of the periodic maintenance procedures as described in these instructions.
- * Specific qualifications and training in the maintenance operations of the appliance in question.
- * The use of components/replacement parts/accessories that are either original or approved by the supplier, in such a way that each operation causes no alteration or modification to the appliance.
- * Possession of the checklist of operations carried out on the appliance.
- * Guarantee complete adherence to the instructions of the Directive 93/42/CEE which includes also the obligation to provide the Manufacturer to maintain post sales records and traceability of the appliance if requested.

A precautionary maintenance for the end terminals is necessary to guarantee safety and to preserve the correct functioning conditions.

The frequency of controls depends on different factors, such as the legal prescriptions, the kind of use, the frequency of use, the environmental conditions during functioning and storage. △

Periodical maintenance services by the Manufacturer or Authorised Centres are not planned, but we suggest to make the following verifications before and after every use of the device, and at least once every 3 months:

- ease of locking, unlocking e disconnection
- wearing or damage
- contamination
- labelling and markings
- air thickness according UNI EN ISO 7396-1 norm (every 3 months)
- supply according UNI EN ISO 7396-1 norm (every 3 months)

The maintenance programme should include an immediate procedure for signalling incorrect functioning to the Manufacturer and for the repairing and substitution. Any operations made on the end unit will be registered in a log book. Only the manufacturer is authorised to substitute any part of the rapid coupling attachment for gas therapy. If spare parts are required, only original ones can be used. Once completed the maintenance, make again the checks required by UNI EN ISO 7396-1 norm.

During all checking, maintenance and cleaning procedures, the operator must wear adequate personal protection such as gloves, mask, glasses etc.

If above mentioned revision is not carried out the conformity to the Directive 93/42/CE for medical devices will no longer be valid therefore, even though the CE mark is present, it is possible that the device no longer answers all the requirements as indicated by the Manufacturer at purchase.

Spencer Italia S.r.l. will take no responsibility the incorrect functioning or any damage caused by a device that has not undergone regular revision.

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

5.2.2 Special servicing

Only the Manufacturer or centres with written authorisation are authorised to complete any special servicing operations.

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device, if used as indicated in the following instruction manual, has an average life span of 5 years.

Spencer Italia S.r.l. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres.

6. ACCESSORIES AND SPARE PARTS

6.1 ACCESSORIES

- OX05103A White label for oxygen inlet
 - OX05106A Blue label for nitrous oxide inlet
 - OX05105A Black/white label for air inlet
 - OX05104A Yellow label for vacuum inlet
- There aren't any spare parts for these devices.

7. WARRANTY

Dear Customer, We would like to thank you for purchasing a Spencer product and we are certain that you will be satisfied with its quality. Spencer Italia S.r.l. guarantees this product from manufacturing or material faults for 12 (twelve) months (according to the law) starting from the original date of purchase.

If in the course of this warranty period you encounter any manufacturing or material faults, our authorised distributors or Spencer itself will repair or substitute the product or the defective component within the terms and conditions indicated below, without any charges for labour or spare parts.

Conditions.

1. This warranty does not cover costs, damages and faults caused by: improper use of the equipment (use which does not conform to the instructions given in the manual); use for means the product is not designed for; lack of maintenance, negligence, tampering with the equipment (e.g. interventions/repairs) over powering or under powering of energy caused by malfunctioning mains; system or equipment defects to which the product was connected; inobservance to the norms of technical/environmental use or to the correct way in which to use the equipment; incorrect installation and use of the product not conforming to the technical norms or safety regulations of the country in which the product is used; fortuitous events, lightning, floods, fires, gas leaks, accidents (knocks etc) or any other cause unrelated to Spencer Italia S.r.l.'s incorrect working/damage caused during transport.
2. This warranty is not valid if the indication of the model or serial number on the product has been modified, cancelled out, removed or made illegible.
3. The following are excluded from the warranty: any interventions for periodical maintenance and the repair or replacement of parts subjected to normal wear and tear; the intervention, repair of breakdowns or the time in repairs determined by the use of the equipment different to what the product is designed for; all of the labour costs of the technical assistance and transport from the client's premises to the servicing centre laboratory and vice versa, as well as any other relevant risks.
4. This warranty does not effect neither the rights of the purchaser established in the current national laws in force, nor does it effect the rights of the client towards the supplier deriving from the sales contract.
5. Any transport costs and risks to and from the cited servicing centres will be at the purchaser's expense.
6. The producer is not responsible of direct or indirect damage of any nature to people or things for the improper use of the product or for the impossibility of use during the time necessary for repairs.
7. For any controversy deriving from the interpretation of this warranty it is of exclusive competence of the judicial authority of the manufacturers place of origin to intervene.
8. The products must be repaired at the producers premises or authorised laboratory. Long with the faulty product the following must also be included: The warranty form correctly completed in all relevant sections for reseller or for end user; proof of purchase containing the lot number and matriculation (transport document for products sold in Italy and the invoice issued by Spencer Italia S.r.l. for products sold abroad); completed form stating the exact fault. Articles sent with the incorrect or incomplete documents will not be covered for warranty.
9. Spencer Italia S.r.l. reserves the right to contest the validity of the warranty certificate when, after objective checking, it is found that the product functioned correctly for a long time before the date indicated as date of purchase. It is recommended to read the instructions of use carefully before contacting Spencer Italia or its authorised distributors.

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All Rights reserved.
Spencer Italia S.r.l. Strada Cavi, 7 - 43044 Collecchio (Parma) ITALY
Copyright Spencer Italia S.r.l.