

# User's Manual

# Switch for medical gas distribution line







 $\mathbf{C} \mathbf{\epsilon}_{0123}$  This appliance conforms with the Directive 93/42/CEE "Medical Devices"

Guarantee of Quality system for the production and the final control of the products certified by the notifying body  $T\ddot{U}V$   $S\ddot{U}D$  Product Service GmbH

# **INDEX**

| General information | page 2 | Operating instructions      | page 4 |
|---------------------|--------|-----------------------------|--------|
| Warnings            | page 2 | Maintenance and cleaning    | page 6 |
| Product description | page 3 | Accessories and spare parts | page 6 |

First edition: 30/06/10 Rev. 2: 23/02/15

#### Thank you for choosing a Spencer product

#### 1. GENERAL INFORMATION

#### 1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

#### Conservation of the instruction manual

The instruction and maintenance manual must be kept together with the product, for the whole life of the device, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

#### Symbols used 1.3 **Symbol** Meaning

General or specific warnings

See instructions for use

Lot number

Product code

The product is compliant with the specifications of the Directive 93/42/CEE

#### 1.4 Servicing request

For any information regarding the correct interpretation of the instruction manual, the use, maintenance, installation and restore of the product, please contact the Spencer Customer Care Service tel. 0039 0521 541111, fax 0039 0521 541222, email service@spencer.it or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY. In order to facilitate the assistance service, please always indicate the lot number (LOT) shown on the label applied on the box or on the device.

#### 1.5 **Demolition**

When the devices are no more suitable for being used, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations about demolition.

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the manufacturer, the product, the CE mark, the lot number (LOT). It must never be removed or covered.



#### 2.1 **General warnings**

- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the Competent Authorities and/or Manufacturer if requested.
- Spencer Italia S.r.l. is always at your disposal to plan trainings on products.
- Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- If the instructions belong to another device and not the device received, inform the Manufacturer immediately and avoid use of the device.
- In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
- Regularly check the appliance, carry out the prescribed maintenance and respect the average life span, as indicated by the manufacturer in this user's manual.
- Before each use of device the perfect operating state of the device must be checked as specified in the Instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and or of the user are detected, the device must be immediately removed from service and the Manufacturer must be contacted.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- Use of the device in anyway other than described in this manual is forbidden.
- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.

- Handle with care.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments.
- When the device is being used, the assistance of qualified staff must be guaranteed.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store in a cool, dry, dark place and do not expose to direct sun.
- Store and transport device in its original packaging.
- The device not be exposed to or come into contact with any source of combustion or inflammable agents.
- Position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Attention: laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario
  for use. This means that in some cases the performance of the product could be notable different from results to date
  obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with
  adequate technical formation.
- With reference to the D. Lgs. 24th February 1997, n. 46 emended by D. Lgs. 25/01/2010, n. 37 Acknowledgement of Directive 93/42/CEE and 2007/47/CE, we remind both public and private operators that they are obliged to report any accident that involves any medical device to the Ministry of Health and to the Manufacture as specified and within time given by the European regulations.
- In addition, both public and private operators are obliged to inform the Manufacturer of any measures that should be
  adopted to make the steps necessary to guarantee the safety and the health of the patients and the users o any medical
  device.
- As a Distributor or End Users of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly
  required to have a basic knowledge of any legal requirements applying to the devices contained in this supply that are in
  power in the goods final destination Country (including laws and norms regarding technical specifications and/or safety
  requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects
  regarding the total conformity of the products to the regulations in the relevant territory.
- Promptly notify Spencer Italia S.r.l. regarding any revisions to be made by Manufacturer in order to guarantee the
  conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of
  other nature).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant User Manual.
- Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the Manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- You are aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully
  responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your
  non-compliance with the present regulatory provisions.



#### 2.2 Specific warnings

- Establish a maintenance program and periodic testing, identifying a reference employee. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in the user's manual.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
- Use only accessories/spare parts that are original or approved by Spencer Italia S.r.l., in order to carry out any operation without causing any alteration or modification to the device, otherwise we assume no responsibility for the proper functioning or damage resulting from device to the patient or the operator and warranty and will be considered void according to the compliance to the Medical Device Directive 93/42/CEE.
- Do not lubricate the product.
- Avoid contact with sharp objects.
- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer and the lack of warranty. The manufacturer won't be anymore responsible for the product..
- The product must be used by trained personnel only.
- Do not force the knob over the end position. It could cause damage.
- Do not use teflon tape or other materials and sealants liquid for input and output connections of the medical gas; their use could cause damage to the device.
- When not in use, the source of medical gas to the deviator must be closed.

#### 2.3 Contraindications and side effects

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

# 3. DESCRIPTION OF PRODUCT

## 3.1 Intended use

The manual deviator has been designed for use in low pressure medical gas distribution systems (oxygen, medical air, vacuum and nitrous oxide) so that in just one simple operation, the continuous supply of gas distribution on appliances is ensured.

#### 3.2 Main components

#### 3.2.1 Main components Pop

Deviation module

Hex keys and fixing nuts for attachment to Pop block

Pop block and caps with sliding block complete with springs

Knob with fixing hex key Turning connections Deviator label for Pop

#### 3.2.2 Main components Rek-Ox

Deviation module Rek-Ox block Knob with fixing hex key Turning connections Deviator label for Rek-Ox

#### 3.3 Models

These models could be modified, with reference to codes and/or descriptions without any previous notification.

OX04306A POP DEVIATION MODULE CYLINDER 1/CYLINDER 2/OFF

OX07200A REK-OX DEVIATION MODULE CYLINDER 1/CYLINDER 2/OFF WITHOUT MONITOX

OX07202A REK-OX DEVIATION MODULE CYLINDER 1/CYLINDER 2/OFF IN LINE

#### 3.4 Technical data

| 5.4 recilificat data       |                   |
|----------------------------|-------------------|
| Characteristics            | Technical data    |
| Dimensions Pop (mm)        | 100 x 102 x 102,5 |
| Dimensions Rek-Ox (mm)     | 97 x 93 x 110     |
| Dimensions Rek-Ox lin (mm) | 90 x 70 x 65      |
| Weight Pop (kg)            | 0,53              |
| Weight Rek-Ox (kg)         | 0,5               |
| Functioning pressure (bar) | 3,5 ± 0,5         |



Pressures higher than 5 bar (500 Kpa) might damage the device.

#### 3.5 Reference standards

| 5.5 Reference standards              |   |  |
|--------------------------------------|---|--|
| Reference                            | Title of document   |  |
| MDD 93/42/CEE                        | European Directive about Medical Devices  |  |
| MDD 2007/47/CEE                      | Modifications to 90/385/CEE Directive about active implants, Directive 93/42/CEE about medical devices and Directive 98/8/CE about the introduction of biocides onto the market |  |
| Legislative Decree 24/02/1997, n. 46 | Application of the 93/42/CEE Directive about Medical Devices  |  |
| Legislative Decree 25/01/2010, n. 35 | Modifications and additions to the 20/02/97 Decree n. 46  |  |
| UNI EN ISO 14971                     | Application of risks managing to medical devices  |  |
| UNI CEI EN ISO 15223-1               | Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements  |  |
| UNI CEI EN 1041                      | Information supplied by the medical devices manufacturer  |  |
| CEI EN 62366                         | Medical Devices - Application of the utilisation characteristics of engineering to medical devices  |  |
| MEDDEV 2.4/1a-b                      | Guideline for the classification of medical devices   |  |
| NB-MED 2.5.1 /Rec 5                  | Technical Documentation   |  |
| MEDDEV 2.7.1                         | Clinical Data   |  |
| MEDDEV 2.12/1                        | Medical Devices vigilance system  |  |
| UNI EN 14155                         | Clinical evaluation of the medical devices for human beings - Part 2:<br>Clinical evaluation plans  |  |

#### 3.6 Environmental conditions

Functioning temperature: from -10 to +50 °C Storage temperature: from -20 to +60 °C Relative humidity: from 5 to 95%

#### 4. OPERATING INSTRUCTIONS

#### 4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

Keep the original packaging for use in case of any further transport and for storage. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

## 4.2 Preparation

On receipt of the product:

- Remove the packaging and display the material so that all components are visible.
- Check that all the components/pieces on the accompanying list are present.
- Connect the two F ¼" inlets G of the turning connections to the non toxic oxygen tube crimped with M ¼" G leading from the two pressure reducers on the cylinders. Use fittings with O-ring seals already integrated.
- Connect the outlet on the opposite side of the knob to the plug block using non toxic oxygen tube crimped with M ¼" G. Use a sealant adapt for use with oxygen or tape for oxygen.



**For Pop module**: Make the calibrated lodgement for the POP module. The dimensions of the hole for a single module is: height 75 mm, length: 92 mm.

Make sure the dimensions of the lodging are as indicated above. Failure to do so would compromise the correct assembly of the module.

- Adjust the Pop module to the thickness of the wall to which it is assembled. Adjustment is done by cutting the plastic extensions on the base of the lodging. The extensions are respectively 3/5/7/10 mm.
- Follow the operations shown in the picture below:
  - Lower the moving part of the module as far as possible. Insert module in the hole made. First push lower part as far as possible into the hole and then the upper part.
  - Check that the part is level with the wall and position the other half on top. Screw up the attachments on the sides.
  - After checking that the module is level with the wall, screw up the attachments on the edges of the module until complete fixation to the wall/piece.



Do not over tighten these attachments, doing so could cause irreparable damages to the module. Screw only until it feels tight and then just give an extra half turn.

 Assemble the lid with its label on the module checking that it is correctly positioned on the structure. Place the red knob on the fixing hub, position and push on knob. Tighten the screws inside until firmly blocked.



Do not over tighten the screw, doing so could cause irreparable damages to the module. Screw only until it feels tight and then just give an extra half turn.

Before use check that the knob is firmly attached and is moving correctly to each position.

- The appliance must be checked before every use so as To reveal any working abnormalities and/or damage caused by transport and/or storage.
- **For Rek-Ox module**: Make the appropriate 4 holes, Ø 4 mm, and secure with countersunk M4 screws class 10.9. In the rear part of the wall to which the Rek-Ox module is fixed, is necessary that there is a counter or a structure to guarantee the resistance of the fastener.



With the gas closed, Position the knob vertically. Open the cylinder so that the gas pressure of inlets of deviator is at correct level for use.

Turn the knob anticlockwise as far as possible to activate cylinder no 1 and check the pressure in the terminal unit.

Turn the knob clockwise as far as possible to activate the cylinder n° 2 and check the pressure in the terminal unit. When no longer in use, reposition the knob vertically in the OFF position.



The knob does not automatically position itself, so the correct position of the knob for the application needed must be checked.

4.4 Troubleshooting

| PROBLEM   | CAUSE  | REMEDY   |
|---|--|--|
| Leakage in the distribution system of medical gas near to the deviator                    | Washers, sealing products or tube clamps are not adapt | Check that the joints are correctly tightened and control leakage in distribution system |
| The system has not closed or changed cylinder with the turning of the knob by 90° or 180° | Hole on Knob for connection to the deviator is damaged | Replace the adjustment knob  |

#### 5. MAINTENANCE AND CLEANING

#### 5.1 Cleaning

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.



The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures.

The exposed metal parts are usually treated and/or painted in order to increase their resistance.


The board has been made out of incontaminable material, in order to increase hygiene and easy cleaning.

Clean the exposed parts with water and delicate soap then dry with a soft cloth. In order to obtain a shine effect, it is possible to use car waxes and creams.

Do not clean with high pressure water; this will damage the joints and the lubricated parts.

If the stretcher is not cleaned regularly, this may cause risks in terms of cross-contamination.

We recommend the use of the polishing detergent Spencer STX 99.



Rinse thoroughly with warm water making sure that you have removed all traces of detergent, which could degrade or compromise the integrity and durability of the device. The use of high pressure water should be avoided. Water penetrates the joints and removes the oil, creating the risk of corrosion of components.

Allow to dry thoroughly before storing. Drying after washing or after use in wet environments must be natural and not forced, do not use flames or other sources of direct heat.

#### 5.2 Maintenance

#### 5.2.1 Precautionary maintenance

The person who carries out the precautionary maintenance of the appliance (user in person, Manufacturer/supplier or a third party) has to guarantee the following basic requirements:

- Technical knowledge of the appliance and of the periodic maintenance procedures as described in these instructions.
- Specific qualifications and training in the maintenance operations of the appliance in question.
- The use of components/replacement parts/accessories that are either original or approved by the supplier, in such a way that each operation causes no alteration or modification to the appliance.
- Possession of the checklist of operations carried out on the appliance.
- Guarantee complete adherence to the instructions of the Directive 93/42/CEE which includes also the obligation towards the Manufacturer to maintain post sales records and traceability of the appliance if requested.



# During all checking, maintenance and cleaning procedures, the operator must wear adequate personal protection such as gloves, mask, glasses etc.

Checks to be carried out before and after each use, and at least every month, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- · Absence of cuts, holes, tears on the structure
- Correct fixation of all nuts, bolts and screws
- Integrity of components

The inspection frequency is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage. Please note that you must do the cleaning as described in paragraph 5.1 and verify functionality before and after each use. Spencer Italia S.r.l. declines any responsibility for the proper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance warranty and will void the compliance to the Medical Device Directive 93/42/CEE. The person responsible for every day maintenance can substitute the spare parts indicated on paragraph 6.2 "Spare Parts", only if authorized by the manufacturer or by a centre authorised by Spencer.

Use only accessories/original spare parts approved by Spencer Italia S.r.l., otherwise we will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres. Warranty will be considered void in compliance with the Medical Device Directive 93/42/EEC.

#### 5.2.2 Periodic maintenance

The device does not require periodic review of the planned interventions by the Manufacturer or an authorized centre, but is prescribed to clean it and check it as described on the paragraphs 5.2.1.



## 5.2.3 Special servicing

# Only the Manufacturer or centres with written authorisation are authorised to complete any special servicing operations.

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device, if used as indicated in the following instruction manual, has an average life span of 5 years. The life span can be expanded only following a general revision of the product that must be carried out by the Manufacturer or by a centre authorised by the Manufacturer.

Spencer Italia S.r.I. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres, making void the guarantee and the conformity to the Medical Devices Directive 93/42/CEE.

#### 6. ACCESSORIES AND SPARE PARTS

#### 6.1 Accessories

The tube for medical gas may be supplied only on specific request with the ends already terminated with threaded connection. Before ordering, please make a specific request for the required length and connections.

#### 6.2 Spare parts

OX05014B Deviator label for caps of Pop (only for Pop models)

OX07031B Red regulation knob

# ATTACHMENT A - TRAINING REGISTER



The product must be used by trained personnel only, having attended specific training for this device and not for similar products.

| Operator's name | Training date  |                   | Training method (user's manual, during service, Trainer             |         |
|-----------------|----------------|-------------------|---|---------|
| Operator's name | Basic training | Advanced training | Training method (user's manual, during service, former class, etc.) | irainer |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   |   |         |
|                 |                |                   | <u> </u>  |         |

#### ATTACHMENT B – MAINTENANCE REGISTER



Keep this document at least 10 years from the end of life of the device.

|     | A)A |   |
|-----|-----|---|
| - 1 | m   | ٨ |
| - 4 | ж   | в |
|     |     |   |

Perform the required maintenance and to respect the life span of the device, as indicated by the Manufacturer in the User's Manual.

| Code and description of the device |  |
|------------------------------------|--|
| Purchase date                      |  |
| Lot (LOT) or serial number (SN)    |  |
| Bought by                          |  |

| SERVICE DATE | KIND OF SERVICE<br>(Maintenance/<br>check/<br>extension of life span) | OPERATIONS MADE ON THE<br>DEVICE | RESULT | PERSON IN CHARGE OF SERVICE (Operator/ Authorized centre/ Manufacturer) |
|--------------|---|----------------------------------|--------|---|
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |
|              |   |                                  |        |   |



#### Warning

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l. Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, lay-out or technical aspects that are herein described.

## © Copyright Spencer Italia S.r.l.

All rights reserved.

No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.