

## B-life / Air Cuffed Mask

Face masks for resuscitation



B-life



Air Cuffed Mask

**CE**<sub>0123</sub> This appliance conforms with the Directive 93/42/CEE "Medical Devices"

Guarantee of Quality system for the production and the final control of the products certified by the notifying body  
TÜV SÜD Product Service GmbH

### INDEX

General information	page 2	Operating instructions	page 5
Warnings	page 2	Maintenance and cleaning	page 6
Product description	page 4	Accessories and spare parts	page 7

First edition: 30/06/10  
Rev. 5: 04/03/14

## 1. GENERAL INFORMATION







### 1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

### 1.2 Conservation of the instruction and maintenance manual

The instruction and maintenance manual must be kept together with the product, for the whole life of the device, inside a dedicated container and above all, away from any substances or liquids which could compromise perfect legibility.

### 1.3 Symbols used

Symbol	Meaning
	General or specific warnings
	See instructions for use
	Lot number
	Product code
	The product is compliant with the specifications of the Directive 93/42/CEE
	Single use

### 1.4 Servicing request

For any information regarding the correct interpretation of the instruction manual, the use, maintenance, installation and restore of the product, please contact the Spencer Customer Care Service tel. 0039 0521 541111, fax 0039 0521 541222, e-mail [service@spencer.it](mailto:service@spencer.it) or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY. In order to facilitate the assistance service, please always indicate the lot number (LOT) shown on the label applied on the box or on the device.

### 1.5 Demolition

When the devices are no more suitable for being used, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations about demolition.

### 1.6 Labelling

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the manufacturer, the product, the CE mark, the lot number (LOT). It must never be removed or covered.

## 2. WARNINGS

### 2.1 General warnings

- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
- Spencer Italia S.r.l. is always at your disposal to plan trainings on products.
- Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- If the instructions belong to another device and not to the device received, inform the manufacturer immediately and avoid use of the device.
- In case of any doubts about the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
- Regularly check the appliance, carry out the prescribed maintenance and respect the average life span, as indicated by the manufacturer in this user's manual.
- Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the

patient and of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.

- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- Use of the device in anyway other than described in this manual is forbidden.
- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- Handle with care.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments.
- When the device is being used, the assistance of qualified staff must be guaranteed.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store in a cool, dry, dark place and do not expose to direct sun.
- Store and transport device in its original packaging.
- The device not be exposed to or come into contact with any source of combustion or inflammable agents.
- Position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Attention: laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
- With reference to the D. Lgs. 24th February 1997, n. 46 emended by D. Lgs. 25/01/2010, n. 37 – Acknowledgement of Directive 93/42/CEE and 2007/47/CE, we remind both public and private operators that they are obliged to report any accident that involves any medical device to the Ministry of Health and to the Manufacture as specified and within time given by the European regulations.
- In addition, both public and private operators are obliged to inform the manufacturer of any measures that should be adopted to make the steps necessary to guarantee the safety and the health of the patients and the users of any medical device.
- As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly required to have a basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the goods final destination Country (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
- Promptly notify Spencer Italia S.r.l. regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user's manual.
- Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- The distributor or final user is aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore Spencer Italia S.r.l. expressly disclaims any responsibility and/or liability for your non-compliance with the present regulatory provisions.



## 2.2 Specific warnings

- The face masks for resuscitation may only be used by personnel trained in cardio-pulmonary resuscitation techniques (CPR or ACLS).
- The resuscitation mask may not be exposed to combustion heat sources and inflammable agents.
- The oxygen supply to the patient in the presence of hydrocarbons creates explosive mixes.
- Do not use the face mask in a polluted environment.
- Before every use, check the general condition of the masks (integrity and cleanness).
- **Air Cuffed Mask resuscitation masks are devices in disposable version**, to be used for a single patient. Do not wash or sanitize after use. Reuse may cause cross-infection. Do not use the device after the expiry date printed on the package.

### 2.3 Contraindications and side effects

The use of the device, if used by personnel, trained in cardio-pulmonary resuscitation techniques (CPR or ACLS), does not present any contraindications and side effects.

## 3. DESCRIPTION OF PRODUCT

### 3.1 Intended use

The resuscitation masks are ideal for use in assisted and controlled ventilation; can be used with both manual resuscitators and active medical devices. The structure has been shaped anatomically and ergonomically in order to guarantee a perfect fit on the patient's face and an optimal grip for the operator.

### 3.2 Main components

- 1 Cushion
- 2 Adult shell
- 3 Mask (for available sizes check paragraph 3.3)
- 4 Valve (only for Air Cuffed Mask)




Fig. A

### 3.3 Models

These basic models could be modified, with reference to codes and/or descriptions without any previous notification.

#### 3.3.1 Models of B-life

RM20400B	B-life - Silicone mask size 0
RM20401B	B-life - Silicone mask size 1
RM20402B	B-life - Silicone mask size 2
RM20403B	B-life - Silicone mask size 3
RM20404B	B-life - Silicone mask size 4
RM20405B	B-life - Silicone mask size 5

#### 3.3.2 Models of Air Cuffed Mask

RM10860A	Air Cuffed Mask – Single patient mask without valve , PVC, size 0, 10 pcs
RM10861A	Air Cuffed Mask - Single patient mask without valve, PVC, size 1, 10 pcs
RM10862A	Air Cuffed Mask - Single patient mask without valve, PVC, size 2, 10 pcs
RM10863A	Air Cuffed Mask - Single patient mask without valve, PVC, size 3, 10 pcs
RM10864A	Air Cuffed Mask - Single patient mask without valve, PVC, size 4, 10 pcs
RM10865A	Air Cuffed Mask - Single patient mask without valve, PVC, size 5, 10 pcs
RM10870A	Air Cuffed Mask - Single patient mask with valve, PVC, size 0, 10 pcs
RM10871A	Air Cuffed Mask - Single patient mask with valve, PVC, size 1, 10 pcs
RM10872A	Air Cuffed Mask - Single patient mask with valve, PVC, size 2, 10 pcs
RM10873A	Air Cuffed Mask - Single patient mask with valve, PVC, size 3, 10 pcs
RM10874A	Air Cuffed Mask - Single patient mask with valve, PVC, size 4, 10 pcs
RM10875A	Air Cuffed Mask - Single patient mask with valve, PVC, size 5, 10 pcs
RM10889A	Air Cuffed Mask - Mask in PVC, colour coded without valve, size 0
RM10890A	Air Cuffed Mask - Mask in PVC, colour coded without valve, size 1
RM10891A	Air Cuffed Mask - Mask in PVC, colour coded without valve, size 2
RM10892A	Air Cuffed Mask - Mask in PVC, colour coded without valve, size 3
RM10893A	Air Cuffed Mask - Mask in PVC, colour coded without valve, size 4
RM10894A	Air Cuffed Mask - Mask in PVC, colour coded without valve, size 5
RM10900A	Air Cuffed Mask - Mask in PVC, colour coded with valve, size 0
RM10901A	Air Cuffed Mask - Mask in PVC, colour coded with valve, size 1
RM10902A	Air Cuffed Mask - Mask in PVC, colour coded with valve, size 2
RM10903A	Air Cuffed Mask - Mask in PVC, colour coded with valve, size 3
RM10904A	Air Cuffed Mask - Mask in PVC, colour coded with valve, size 4
RM10905A	Air Cuffed Mask - Mask in PVC, colour coded with valve, size 5

### 3.4 Technical data

Component	B-life	Air Cuffed Mask
Paediatric mask	Silicone	PVC
Cushion	Silicone	PP
Adult shell	Polycarbonate	-
Package	Non sterile	Non sterile



The masks contain high molecular weight phthalates (DINP type) which are not considered to be harmful to human health.

### 3.5 Reference standards

Reference	Title of document
MDD 93/42/CEE	European Directive about Medical Devices
MDD 2007/47/CEE	Modifications to 90/385/CEE Directive about active implants, Directive 93/42/CEE about medical devices and Directive 98/8/CE about the introduction of biocides onto the market
Legislative Decree 24/02/1997, n. 46	Application of the 93/42/CEE Directive about Medical Devices
Legislative Decree 25/01/2010, n. 35	Modifications and additions to the 20/02/97 Decree n. 46
UNI EN ISO 13485	Medical devices – Quality management systems – Requirements for regulation purposes
UNI EN ISO 14971	Application of risks managing to medical devices
UNI CEI EN ISO 15223-1	Medical devices - Symbols for use in the medical device labels, labelling and information to be provided. Part 1: general requirements
UNI CEI EN 1041	Information supplied by the medical devices manufacturer
CEI EN 62366	Medical Devices - Application of the utilisation characteristics of engineering to medical devices
MEDDEV 2.4/1a-b	Guideline for the classification of medical devices
NB-MED 2.5.1 /Rec 5	Technical Documentation
MEDDEV 2.7.1	Clinical Data
MEDDEV 2.12/1	Medical Devices vigilance system
UNI EN 14155	Clinical evaluation of the medical devices for human beings - Part 2: Clinical evaluation plans

### 3.6 Environmental conditions

Functioning temperature: from -20 to +60 °C

Storage temperature: from -20 to +60 °C

## 4. OPERATING INSTRUCTIONS

### 4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

Keep the original packaging for use in case of any further transport and for storage. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The device must be stored in a dry and cool area.

### 4.2 Preparation

Upon receipt, the mask requires a visual inspection, which must be repeated before each use to verify its integrity. The mask is used with the ventilation equipment available.

### 4.3 Functioning

1. Keep the patient lying down, holding his face facing upward.
2. Unleashing the mouth and throat of the patient from any foreign bodies and secreted.
3. It is suggested to introduce an oropharyngeal airway into the patient's mouth before using a manual resuscitator.
4. Stand behind the patient's head, holding it backward, then at the same time pull the chin upwards and backwards.
5. Place the mask over the patient's face at the nose and mouth, keeping it stable with the forefinger and thumb.
6. Observe the patient's face and lip colour.
7. Use a mask that completely covers both the nose and the mouth of the patient, covering the entire nasal pyramid and the lower lip, while excluding the chin, which must remain open.
8. Free the airways by removing, if necessary, vomiting, liquids, or foreign objects that obstruct them.
9. Hyperextend the patient's neck with due caution.

#### 4.4 Troubleshooting

PROBLEM	CAUSE	REMEDY
The mask let air despite of its correct position	The face of the patient has a moustache or beard or the skin is sweaty.	In the case of a moustache or beard try to make a greater pressure on the mask, in the case of clammy skin wipe the face of the patient.
	The measure of the mask used is not compatible with the morphology of the face of the patient.	Replace the mask using a size compatible with the and structure the patient's face.
The mask has structural or functional defects.	Wear of one or more components.	Immediately put out of service the device and contact the service centre.

### 5. MAINTENANCE AND CLEANING

#### 5.1 Cleaning

Cleaning methods	B-life (silicone)	Air Cuffed Mask (single patient)
Wash with detergent (avoid detergents containing phenol)	Yes	No
Disinfection with disinfectant	No	No
Sterilization with ethylene oxide through a validated cycle (UNI EN ISO 11135-1)	No	No
Autoclaving, through a validated cycle (UNI EN ISO 17665-1)	Yes	No

Make sure you have completely dried the mask in order to avoid the presence of residues of disinfectants that may compromise the integrity and durability of the mask.

The use of the autoclave for sterilization decreases the life time of the mask.



Failure to perform the cleaning operations can involve the risk of cross-infection due to the presence of secretions and/or residuals.

#### 5.2 Maintenance

##### 5.2.1 Precautionary maintenance

The person who carries out the precautionary maintenance of the appliance (user in person, Manufacturer/supplier or a third party) has to guarantee the following basic requirements:

- Technical knowledge of the appliance and of the periodic maintenance procedures as described in these instructions.
- Specific qualifications and training in the maintenance operations of the appliance in question.
- The use of components/replacement parts/accessories that are either original or approved by the supplier, in such a way that each operation causes no alteration or modification to the appliance.
- Possession of the checklist of operations carried out on the appliance.
- Guarantee complete adherence to the instructions of the Directive 93/42/CEE which includes also the obligation towards the Manufacturer to maintain post sales records and traceability of the appliance if requested.

**The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures.**

Checks to be carried out before and after each use, and at least every 3 months, are as follows:

- General functionality of the device.
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections).
- Absence of cuts, holes, tears on the structure.
- Status of wear and damage to the device due to the number of sterilization cycles undergone.
- Expiration date or average life-time, if any.

**Periodic review of planned interventions by the manufacturer or service center authorized by him are not required, but it is required to carry out cleaning and inspections indicated in the relevant sections "Cleaning" and "Precautionary Maintenance".** Spencer Italia S.r.l. declines any responsibility for the proper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance warranty and will void the compliance to the Medical Device Directive 93/42/CEE.



**For disposable devices it is necessary to check the expiration date/time of life, but there are no cleaning and maintenance activities, being non-reusable devices.**

### 5.2.2 Special servicing



**Only the Manufacturer or centres with written authorisation are authorised to complete any special servicing operations.**

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

**The device, if used as indicated in the following instruction manual, has an average life span of 5 years.**

Spencer Italia S.r.l. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres, making void both the guarantee and the conformity to the Medical Devices Directive 93/42/CEE.

## 6. ACCESSORIES AND SPARE PARTS

### 6.1 Accessories

There aren't accessories for these products.

### 6.2 Spare parts

There aren't spare parts for these products.



#### Warning

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l. Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, lay-out or technical aspects that are herein described.

© Copyright Spencer Italia S.r.l.

All rights reserved.

No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.