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IT

Instructions for use
PICK-UP STRETCHERS

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EN

1. MODELS	12
2. INTENDED USE	12
2.1 INTENDED USE AND CLINICAL BENEFITS	12
2.2 TARGET PATIENTS	12
2.3 PATIENT SELECTION CRITERIA	12
2.4 CONTRAINDICATIONS AND SIDE EFFECTS	12
2.5 USERS AND INSTALLERS	12
2.5.1 USER TRAINING	12
2.5.2 INSTALLER TRAINING	12
3. REFERENCE STANDARDS	12
4. INTRODUCTION	12
4.1 USING THE INSTRUCTIONS FOR USE	12
4.2 DEVICE LABELLING AND TRACEABILITY	12
4.3 SYMBOLS	13
4.4 WARRANTY AND SERVICE	13
5. WARNINGS/DANGERS	13
6. SPECIFIC WARNINGS	13
7. RESIDUAL RISK	14
8. TECHNICAL DATA AND COMPONENTS	15
9. COMMISSIONING	15
10. PROPER USE	16
10.1 APPLYING IMMOBILIZATION ACCESSORIES	16
10.2 PREPARATION AND POSITIONING	16
10.3 STORING THE STRETCHER	17
10.4 USING FIXING SYSTEMS	17
11. CLEANING AND MAINTENANCE	17
11.1 CLEANING	18
11.2 ROUTINE MAINTENANCE	18
11.3 PERIODIC OVERHAUL	18
11.4 SPECIAL MAINTENANCE	18
11.5 LIFE SPAN	18

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- SX RED KEY 2 EDGES
- SX RED KEY 3 EDGES

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Pick-up stretchers are lifting devices designed to move severely injured patients to other transport devices.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity and within the limits of the size of the device.

2.3 PATIENT SELECTION CRITERIA

The expected patients are typically subjects for whom recovery (collection) is necessary with as little travel as possible in addition to transfer to another transport device.

It is common practice to use a head restraint with the device to prevent any aggravation of neck injuries.

Operators must also be trained to perform device application manoeuvres.

2.4 CONTRAINDICATIONS AND SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device.

2.5 USERS AND INSTALLERS

The intended users are rescue workers with in-depth knowledge related to the immobilisation and handling of people suffering from road traffic injuries, spinal injuries and fractures.

- Personnel trained for use of the device must also have training in managing lifting and handling suspended loads with people.

These devices are not intended for lay people.

Pick-up stretchers are devices intended for professional use only. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs, should it be necessary to lift and/or support the device and the patient. Operators' ability must be assessed before the definition of roles in use of the stretcher.

Operators must be able to provide the necessary patient care.

2.5.1 USER TRAINING

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of these instructions for use before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.

- The product must be used only by personnel trained in the use of this product and not on other similar products.

- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**

- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses.

2.5.2 INSTALLER TRAINING

The installer must have the appropriate skills and qualifications to ensure correct attachment of any ambulance restraints. The device in itself does not require installation.

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices
UNI EN 1865-1	Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment

4. INTRODUCTION

4.1 USING THE INSTRUCTIONS FOR USE

The purpose of these instructions for use is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: The instructions for use are an integral part of the device and therefore must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer product instructions for use can be downloaded from the site <http://support.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to read and understand the contents of these instructions before installing, operating, or servicing this product.


4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). **This must never be removed or covered.**

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it>, or inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the instructions for use
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal law restricts the sale of this device by or on the order of a licensed professional (US market only)
		Production identification	
		Alphanumeric code that identifies the production units of the device, composed of:	
(01)8057711230006 (11) 200626 (10) 1234567890		(01)805771123	company prefix
		000	progressive GS1
		6	control number
		(11)200626	date of production (YYMMDD)
		(10)1234567890	lot number

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at <http://support.spencer.it>.

Note: record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal (<http://service.spencer.it/>) available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

5. WARNINGS/DANGERS

 Warnings, dangers, notes, and other important safety information are provided in this section and are clearly visible throughout the instructions for use.

Product features

 Use of the product for any purpose other than that described in the instructions for use is prohibited.

- Before each use, always check the conditions of the product, as specified in the instructions for use. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
-  If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The product must not be tampered with or modified without the manufacturer's authorisation (modification, tweaking, additions, repair, use of non-approved accessories), as they may constitute imminent danger of injury to persons and material damage. Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE marking and the product warranty shall be null and void.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -10°C to +60°C.

Regulatory requirements

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the instructions for use.
- **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

The user must carefully read the following in addition to the general warnings.

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-  Application of the device is not expected to last longer than the time required for recovery and transfer to the device to be used for ambulance transport (e.g. Main stretcher).
-  Qualified personnel and at least two operators must be present during use of the device.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Follow the internal procedures and protocols approved by your organisation.
-  Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.
- Do not use drying machines to dry the device.

6. SPECIFIC WARNINGS

To use the device, you must also have read, understood and carefully follow all the instructions in the instructions for use.

-  Always comply with the maximum capacity, if any, indicated in the instructions for use. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. In addition, the operator should assess whether the overall size of the patient reduces the functionality of the product.
- The maximum weight, which weighs on each operator, must comply with local health and safety requirements.
- Make sure that operators are in good physical condition before lifting, as listed in the instructions for use.
- Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with routine maintenance of the device must ensure the basic requirements envisaged by the manufacturer within these operating instructions.
- All maintenance activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the manufacturer when requested.
- Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.
- Never leave the patient on the device unsupervised, as they could get injured.
- Avoid contact with sharp objects.
- Follow approved Emergency Medical Service procedures for patient immobilization and transportation.
-  Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Before lifting, make sure that operators have a secure grip on the supporting structure of the device.
- **Perform rescue simulations with a stretcher and a patient simulating load and accessories before putting the device into service.**
-  At least 2 operators in suitable physical conditions are required for use of the device. They must be endowed with strength, balance, coordination, common sense and must be trained on the correct operation of the device.
- For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of more than two operators is recommended in addition to the minimum foreseen operators.
- Before each use, always check the conditions of the device and its components, as specified in the instructions for use. In case of faults or damage that may compromise the functionality and safety of the device, and therefore the patient and the operator, remove the device from service or replace the components that are not intact.
- Do not lift the stretcher if the weight is not properly distributed.
- Use the devices only as described in these instructions for use.
-  Do not alter or modify the device to adapt it to unforeseen conditions of use: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- **The device is expected to come into contact through the patient's clothing. In case of direct contact with the skin, place a protective surgical drape over it to avoid contamination from substances that may have contaminated the device.**
-  During lifting, the stretcher must be kept horizontal with respect to the ground. Any abnormal inclination can cause serious damage to the patient, the device and the operator.
-  To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- Always observe the maximum applicable static safety load indicated in these instructions for use. Maximum static load is intended as a mechanical force applied slowly, not quickly, beyond which the device may not be safe. This value does not take into account the dynamic forces to be added to the static load, such as shocks, vibrations and possible weather and climate conditions during use of the device.
- Never leave the patient unsupervised when the device is in use, as they could get injured.
- After washing, the device and all its components must be left to dry completely before storage away from sunlight and direct heat sources.
- Avoid contact with sharp objects.
- Never use solvents or stain removers.
- Do not use the device if there are cracks, abrasions or burns on the structure or seams or fraying on the belts.
- Do not drag the device on rough surfaces.
- Always check the conditions of all parts of the straps and buckles before each use.
- Immediately replace devices with worn or damaged straps or buckles.
- Position and adjust the straps in such a way that they do not hinder rescuer operations or use of the use of rescue equipment.
- Keep the appropriate documentation for a period of ten years from the date of transfer to the final consumer and, therefore, show, where required, to trace the origin of the products.
- Use of the device without proper patient immobilization can result in serious damage. Always make sure that the patient is properly immobilized before using the device.
- The rescuer must assess the actual need for use of this type of device in accordance with local guidelines.
- If the device is used with its own dedicated Spencer fixing system, make sure that this system has been correctly installed. After placing the device on its fixing, make sure that it is firmly anchored to its stop.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in these instructions for use, have been identified.

IT

EN

DE

FR


ES

PT


EL

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

	PARTS	DESCRIPTION AND MATERIALS
 <p>SX 2 LAE/SX 3 EDGES</p>	1	Foot-side locking hook Made of aluminium, with a quick-release system, this allows the two halves of the stretcher to be disconnected.
	2	Leg edge Made of aluminium, provides support for the patient's legs
	3	Removable profile Made of aluminium, equipped with holes for length adjustment
	4	Stops Made of aluminium, allow sliding of the removable profile to be locked in one of the pre-set positions
	5	Chest/head edges Made of aluminium, there are two in the three-blade models and one in the two-blade models. They allow the patient to be picked up and supported
	6	Handles Recessed into the aluminium structure of the device, they allow operators to grasp for lifting.
	7	Head-side locking hook Made of aluminium, with a quick-release system, this allows the two halves of the stretcher to be disconnected.

FIXING ACCESSORIES

	PARTS	No.	COMPONENT DESCRIPTION	MATERIALS
	1	1	Spine supports	Steel
	2	2	Rod for wall installation	Steel
	3	3	Stops	Nylon/Brass

Characteristics	Spencer SX 2 edges	Spencer SX 3 edges
Shortened length	1660 mm	1660 mm
Intermediate length 1	1780 mm	1780 mm
Intermediate length 2	1900 mm	1900 mm
Maximum length	2050 mm	2050 mm
Length when folded	1190 mm	1190 mm
Thickness when folded	70 mm	70 mm
Thickness when extended	60 mm	60 mm
Width	440 mm	440 mm
Materials	Aluminium	Aluminium
Load capacity	170 kg	170 kg
Weight	7,6 ± 0,4 kg	7,9 ± 0,4 kg

Length measurements are subject to tolerances of ± 10 mm.

Characteristics	SXO	SXV
Dimensioni	510x25x60 ± 5mm	/
Dimensioni elemento lungo	/	440x25x95 ± 5mm
Dimensioni elemento corto	/	180x60x60 ± 5mm
Materiali	Steel, Nylon, brass	
Peso	1,5 ± 0,2 kg	790 ± 10 g

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- Product cleanliness
- Absence of cuts, holes, tears or abrasions on the entire device

Do not modify the device or its parts for any reason as this could cause damage to the patient and/or rescuers.

 Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself. For subsequent use, perform the operations specified in paragraph 11.

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If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

10. PROPER USE

Primary medical evaluations must be carried out before intervening on the patient.

Before using the pick-up stretcher, you should carefully read the instructions for use for all accessories that you plan to use with the device, such as belts, cervical collars and spine boards.

Strictly follow the local emergency service guidelines before placing the patient on the pick-up stretcher.

10.1 APPLYING IMMOBILIZATION ACCESSORIES


Before placing the patient on the stretcher, it is generally advisable to immobilise the cervical spine using a cervical collar.

If necessary, prepare the head restraint, placing it on the pick-up stretcher.



Super SX is the Spencer head restraint dedicated to this type of stretcher.

To apply this head restraint, carefully clean the surface of the edge in correspondence with the patient's head. Then apply the fabric portion using the adhesive side. Coupling with each half of the base will be made by means of the "buckle-loop" surfaces.

Make sure that the base of the head restraint is firmly attached to the stretcher.



10.2 PREPARATION AND POSITIONING



Open the pick-up stretcher, placing it on the floor so that the two halves are perfectly aligned.

When the stretcher is level, push the foot end towards the head end until you feel the locking mechanism engage.

Then place the stretcher at the patient's side.

Unlock the locking systems by turning them outwards, then adjust the length of the stretcher according to the patient's height, closing the locks when one of the available positions is set. Make sure that the two halves of the stretcher are securely locked together.

Separate the two halves of the pick-up stretcher by pressing the buttons at the ends.

For one-piece straps, lay them flat on the floor at the shoulders, abdomen and pelvis. For two-piece belts, apply them by means of a loop knot at the handles, then spread them out to the sides.

While two operators lightly lift one side of the patient, the third operator carefully inserts one half of the stretcher underneath the patient. The operation must then be carried out on the opposite side as well, taking care to match the attachment points between the two halves of the stretcher both to ensure that the connection is secure and to

IT
EN

DE

FR

ES

PT

EL

prevent further movement of the patient.
 Make sure that all locking systems are properly engaged.
 Fasten and tighten the straps to ensure patient safety.


You can now lift and transfer the patient. If a head restraint is not used, an operator should hold the head stabilised during the transfer.
 Never use cranes, slings or other systems to lift the device.
 Lifting is only permitted by hand and by an appropriate number of operators (minimum 2).
 Operators must be positioned symmetrically and in such a way to allow safe, stable lifting of both the head and foot sides of the stretcher.
 The maximum load on each operator must never exceed the maximum load permitted by work safety regulations.

After reaching the primary stretcher for transport, remove the pick-up stretcher, as this is not suitable for ambulance transport.

10.3 STORING THE STRETCHER

After use, the stretcher can be folded and stored in dedicated compartments in the ambulance or where the stretcher is normally stored.
 To fold the stretcher, unlock the hooks described above and extend the stretcher until the tubes are fully extended, i.e., the joints are fully exposed.
 Once this position has been reached, the foot end can be folded over the head end, thus reducing the space required.


Caution: If the joints are not fully extended and an attempt is made to close the stretcher, the stretcher may become irreparably damaged.


10.4 USING FIXING SYSTEMS

To use the SXV fixings, close the pick-up stretcher completely making sure that it is locked (not folded) and orient it so that the head end is facing up and the patient end is facing the wall.

Place the head-end frame on the top stop, guide the stretcher against the wall and, after bringing the stretcher up against the bottom stop, turn the locking handle to a horizontal position. Make sure that the stretcher is properly secured by making a few repeated pulling movements.

To use the SXO fixings, close the pick-up stretcher completely making sure that it is locked (not folded) and orient it so that the patient side faces the wall.
 Place the right or left side of the frame on the supports at the top of the stop. Slowly lower the stretcher towards the wall until it comes to rest against it. Rotate the locking handles to a vertical position. Make sure that the stretcher is properly secured by making a few repeated pulling movements.

11. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the instructions for use, identifying a reference person who meets the basic requirements set forth in the instructions for use.
- **The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.**
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l.

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declines any responsibility for any direct or indirect damage which is a consequence of improper use of spare parts and/or any repair work carried out by unauthorised parties.

- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the instructions for use in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- If required, the product and all its components must be washed and left to dry completely before storage.

11.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

After complete drying, proceed with lubrication as described below.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices.

Carefully follow the instructions of the manufacturer of the product used with regards to the application method and contact time.

Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

11.2 ROUTINE MAINTENANCE

You must establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged within these instructions for use.

All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

The maintenance schedule must comply with the following table:

	DISINFECTION	AT EACH USE	IF NECESSARY	MONTHLY
Disinfection		•		
Cleaning			•	
Inspection		•	•	•

The inspection to be carried out after each use involves:

- Check that all components are present
- Check device conditions – There must be no cracks, fissures, holes or cuts.
- Check conditions of wear – There must not be any level of abrasion that would compromise the safety of the product, for example the thinning of one or more of its parts.
- Check that the moving parts slide properly.
- General check of the state of wear of each component
- Check that all the accessories provided are present, functional and intact.
- Disinfection

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

11.3 PERIODIC OVERHAUL

No periodic overhaul is foreseen for the device.

11.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

11.5 LIFE SPAN

The device and its fasteners, if used as described in the following instructions, has an average life span of 5 years from the date of purchase.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have exceeded the maximum allowable life span.

12. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The device does not carry the load correctly, showing excessive bending	Damaged parts or excessive load	Immediately remove the device from service and contact the manufacturer
The locks do not open/close	Dirty or spring broken	Clean the unit thoroughly, removing any debris. If the problem persists, immediately remove the device from service and contact the manufacturer.
Lower limb extension does not unlock for storage	The joint has not fully extended	Stretch the stretcher further until the joint is fully extended
The support profiles rotate around the edge	Broken fixing rivets	Remove the device from service and contact customer service.

If a problem or fault is detected that is not among those indicated above, please contact Spencer Italia srl customer care service.

13. ACCESSORIES

ST05055	SX BAG - YELLOW PICK-UP STRETCHER CARRY BAG	ST05018	SVX – Vertical fixing system for 20G pick-up stretchers
ST05019	SXO – Horizontal fixing system for 20G pick-up stretchers	SH00500	SUPER SX YELLOW HEAD RESTRAINT FOR PICK-UP STRETCHER

14. SPARE PARTS

ST00598	STX 598 - Two piece strap with plastic hook	ST00592	STX 592 - Two piece strap with metal hook
ST00523	K BELT 2 - SET OF 3 BELTS STX 592 W/BAG	ST00594	K BELT 1 - SET OF 3 BELTS STX 598 W/BAG
ST00593	K BELT 3 - OF 3 BELTS STX 597 W/BAG	ST00597	STX 597 - 1PC ORANGE BELT PLASTIC H.

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

EL

Notice

The information in these instructions for use is subject to change without notice.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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