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EN

1. MODELS	7
2. INTENDED USE	7
2.1 INTENDED USE AND CLINICAL BENEFITS	7
2.2 TARGET PATIENTS	7
2.3 PATIENT SELECTION CRITERIA	7
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	7
2.5 USERS AND INSTALLERS	7
3. REFERENCE STANDARDS	7
4. INTRODUCTION	7
4.1 DEVICE LABELLING AND TRACEABILITY	7
4.2 SYMBOLS	7
4.3 WARRANTY AND SERVICE	7
5. WARNINGS/DANGERS	7
6. SPECIFIC WARNINGS	8
7. RESIDUAL RISK	8
8. TECHNICAL DATA AND COMPONENTS	8
9. PROPER USE	8
10. CLEANING AND MAINTENANCE	8
10.1 CLEANING	8
10.2 MAINTENANCE	8
10.3 LIFE SPAN	8
11. TROUBLESHOOTING TABLE	9
12. ACCESSORIES	9
13. SPARE PARTS	9
14. DISPOSAL	9

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- NRT WHITE

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Recovery stretchers are devices to be used when the morphological characteristics of the intervention area require particularly compact equipment and adherence to the patient.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity and within the limits of the size of the device.

2.3 PATIENT SELECTION CRITERIA

The patients foreseen are those with injuries that prevent them from walking in a given rescue situation, or who are in a state of unconsciousness.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

2.5 USERS AND INSTALLERS

The intended users are rescue workers, typically technical, carrying out operations related to the use of harnessing systems.

- Personnel trained for use of the device must also have training in managing lifting and handling suspended loads with people.
- Personnel who carry out interventions in situations classified as high risk or which are purely technical must be suitably trained and experienced in rescue. These devices are not intended for lay people.

Operators must be able to provide the necessary patient care.

- The product must be used only by personnel trained in the use of this product and not on other similar products.
- Training in the use of the product must be recorded and supported by practical tests on use, the result of which must be documented and certified by a training officer. It is advisable to repeat training annually. This documentation must be kept for at least 10 years after the end of the device's life. .

3. REFERENCE STANDARDS

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

4. INTRODUCTION

4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

If the assigned Lot/SN cannot be traced, the device must be reconditioned, provided only under the responsibility of the manufacturer.

4.2 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the user manual.
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner (only for USA Market)

4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it

Warranty and service conditions are available on the website <http://support.spencer.it>

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

5. WARNINGS/DANGERS

Product features

Use of the product for any purpose other than that described in the User Manual is prohibited.

- The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to + 50°C.
- Storage temperature: -10°C to +60°C.

General warnings for medical devices

- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities

IT

EN

DE

FR

ES

PT

EL

DA

NL

for their respective actions.

- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

6. SPECIFIC WARNINGS


To use NRT, you must also have read, understood and carefully follow all the instructions in the user manual.

-  **Perform simulations with dummies before putting the device into service.**
- To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The device must only be used by trained personnel.
- Qualified personnel must be present during use of the device.
- The device must not be exposed, much less come into contact with thermal sources of combustion or flammable agents.
- Only use the device as described in this manual.
- Always check the conditions of all parts before use.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in this user manual, have been identified.

8. TECHNICAL DATA AND COMPONENTS

1	Head side handle
2	Head stabilisation strap
3	Adjustable belts
4	Side handles
5	Foot side handle
6	Transport bag
Dimensions	1080x1830x20 ± 5mm
Weight	7 ± 0,5 kg
Capacity	220 kg (Using all handles simultaneously)

9. PROPER USE

Use of the device requires a minimum of two operators, who must always operate in compliance with work safety regulations, particularly regarding the maximum weight that weigh on each operator.

To use the device:

- Remove the stretcher from the transport bag and lay it on the floor with the head straps facing up.
- Trained personnel then position the patient on the stretcher using recognised and approved patient handling techniques, taking into account the patient's clinical condition.
- Stabilise the patient's head using the head strap.
- Wrap the stretcher around the patient and secure it with the straps, making sure to tighten them properly so that the patient is stable.
- Then proceed with recovery manoeuvres.

10. CLEANING AND MAINTENANCE

10.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

The use of high-pressure water should be avoided.

If the device gets wet, let it dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

Once dry, clean the device using a brush with stiff bristles.

Do not use chemicals.

Do not immerse in water unless the device has been immersed in salt water. If so, immerse in cold water for about 2 hours. Rinse several times until the salt is completely removed and let it dry naturally.

If blood is visible, oxidise it before washing with water.

During all operations, be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

10.2 MAINTENANCE

The device does not require a routine maintenance program, but checks should be carried out every 3 months and after each use to verify:

- General functionality of the device
 - Wear of any parts of the fabric
 - Conditions of internal rods
 - Conditions of belt grommets
 - Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- No periodic overhaul is foreseen for the device.

10.3 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase.

11. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
Excessive patient mobility	Incorrect immobilisation	Re-position the straps in order to stabilise the patient in the device.
Internal rod breakage	Improper use	Immediately remove the device from service and contact the manufacturer
Damage to grommets	Improper use	Immediately remove the device from service and contact the manufacturer
Fabric wear/rips	Normal wear or improper use	Immediately remove the device from service and contact the manufacturer

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

12. ACCESSORIES

ST30153 NRT EXTENSIONS KIT

13. SPARE PARTS

RIST129 NRT TRANSPORT BAG

14. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

EL

DA

NL

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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