**Istruzioni d'uso
TELI PORTAFERITI**

IT

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EN

1. MODELS	12
2. INTENDED USE	12
2.1 INTENDED USE AND CLINICAL BENEFITS	12
2.2 TARGET PATIENTS	12
2.3 PATIENT SELECTION CRITERIA	12
2.4 CONTRAINDICATIONS AND SIDE EFFECTS	12
2.5 USERS AND INSTALLERS	12
2.5.1 USER TRAINING	12
2.5.2 INSTALLER TRAINING	12
3. REFERENCE STANDARDS	12
4. INTRODUCTION	12
4.1 USING THE INSTRUCTIONS FOR USE	12
4.2 DEVICE LABELLING AND TRACEABILITY	12
4.3 SYMBOLS	13
4.4 WARRANTY AND SERVICE	13
5. WARNINGS	13
6. SPECIFIC WARNINGS	14
7. RESIDUAL RISK	14
8. TECHNICAL DATA AND COMPONENTS	14
9. COMMISSIONING	15
10. PROPER USE	15
10.1 PATIENT POSITIONING	16
10.2 ROD INSERTION (WOW ONLY)	16
10.3 APPLYING AND ADJUSTING SIDE STRAPS (WOW ONLY)	16
10.4 USING THE FOOTREST (IF PRESENT)	16
10.5 LIFTING THE PATIENT IN A SUPINE POSITION	16
10.6 LIFTING THE PATIENT IN A SEATED POSITION	17
10.7 ACCESSORY POCKET (WOW ONLY)	17
11. CLEANING AND MAINTENANCE	17
11.1 CLEANING	17
11.2 ROUTINE MAINTENANCE	17
11.3 PERIODIC OVERHAUL	17
11.4 SPECIAL MAINTENANCE	17
11.5 LIFE SPAN	17
12. TROUBLESHOOTING TABLE	18
13. ACCESSORIES	18
14. SPARE PARTS	18
15. DISPOSAL	18

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SK

1. MODELS

The following basic models may be subject to implementation or change without notice.

- HIRVIN 6 HANDLES
- HIRVIN 8 HANDLES
- HIRVIN 10 HANDLES
- HIRVIN PLUS
- WOW
- WOW WITH FOOTREST

The above models are produced in different colour variations.

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Injury transportation sheets are devices to be used to move and transport patients in a lying or seated position. It is not foreseen that patients be able to intervene on devices.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity and within the limits of the size of the device. If paediatric subjects must be transported, it will be in the role of the rescuer to determine whether the devices are suitable for their transport or if it will be necessary to use another device.

2.3 PATIENT SELECTION CRITERIA

The patients foreseen are those with injuries that prevent them from walking in a given rescue situation.

2.4 CONTRAINDICATIONS AND SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the instructions for use.

2.5 USERS AND INSTALLERS

The intended users are persons trained in first aid procedures and in use of medical equipment in an emergency medical service (EMS) environment. These devices are not intended for lay people.

Injury transportation sheets are devices intended for professional use only. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs, should it be necessary to lift and/or support the device and the patient. Operators' ability must be assessed before the definition of roles in use of the stretcher.

Operators must be able to provide the necessary patient care.

2.5.1 USER TRAINING

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of these instructions for use before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses.

2.5.2 INSTALLER TRAINING

The installer must be able to ensure safe storage inside the ambulance. The device in itself does not require installation.

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices
UNI EN 1865-1	Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment

4. INTRODUCTION

4.1 USING THE INSTRUCTIONS FOR USE

The purpose of these instructions for use is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: The Instructions for use are an integral part of the device and therefore must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer products User Manuals can be downloaded from the site <http://support.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to read and understand the contents of these instructions for use before installing, operating, or servicing this product.

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking or lot number (LOT). This must never be removed or covered.

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it>, or inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning						
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.						
	Medical device		See the instructions for use						
	Manufacturer		Lot number						
	Date of manufacture		Product code						
	Unique Device Identifier		Caution: Federal law restricts the sale of this device by or on the order of a licensed professional (US market only)						
 <p>(01)8057711230006 (11) 200626 (10) 1234567890</p>		Production identification Alphanumeric code that identifies the production units of the device, composed of:							
		(01)805771123	company prefix	000	progressive GS1	6	control number	{11}200626	date of production (YYMMDD)

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at <http://support.spencer.it>.

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal (<http://service.spencer.it>) available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

5. ⚠️ WARNINGS

Warnings, notes, and other important safety information are provided in this section and are clearly visible throughout the instructions for use.

⚠️ Product features

Use of the product for any purpose other than that described in the instructions for use is prohibited.

- Before each use, always check the conditions of the product, as specified in the instructions for use. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations.
- The product must not be tampered with (modification, tweaking, additions, repair). Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE certification (when required by law) and the product warranty shall be null and void.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the structure of the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -10°C to +60°C.

Regulatory requirements

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfillments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the instructions for use.
- **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfillment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

⚠️ The user must carefully read the following in addition to the general warnings.

- The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point.
- Qualified personnel and at least four operators must be present during use of the device.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Follow the internal procedures and protocols approved by your organisation.

IT

EN

DE

FR

ES

PT

RO

CS

SK

- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.
- Do not use drying machines to dry the device.

6. SPECIFIC WARNINGS

To use the device, you must also have read, understood and carefully follow all the instructions in the instructions for use.


- Always comply with the maximum capacity, if any, indicated in the instructions for use. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. Moreover, the operator should assess whether the overall size of the patient reduces the functionality of the product.
- If foreseen for the device, make sure that operators are in good physical condition before lifting, as listed in the instructions for use.
- Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with routine maintenance of the device must ensure the basic requirements envisaged by the manufacturer within these instructions for use.
- To protect patients from possible contamination and/or reactions resulting from contact with the materials, as well as the device itself, a surgical drape should be placed between the patient and the device to separate it from the skin, in accordance with biocompatibility requirements.
- Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Before any type of handling, make sure that operators have a secure grip on the device.
- Do not drag the device over surfaces that could tear it.
- Do not lift with a crane or other mechanical lift.
- Do not use drying machines.
- Do not use the device in conjunction with other equipment not expressly approved by the manufacturer.
- Perform rescue simulations with a device and a patient simulating load and accessories before putting the device into service.**
- For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 4 as expected in standard conditions).
- Before handling the sheet with the patient on it, coordinate in order to perform synchronised lifting and handling.
- Before each use, always check the conditions of the device, as specified in the instructions for use. In case of faults or damage that may compromise the functionality and safety of the device, and therefore the patient and the operator, replace the damaged parts.
- The drape should not be used in the event of respiratory and/or cardiac problems.** Ensure that the clinician has carried out the necessary primary assessments.
- The device is not suitable for use in cases of spinal trauma and where the use of rigid aids is required.
- Make sure that the sheet does not interfere with any handling/holding/lifting system of the main stretcher.
- Do not move the device if the weight is not properly distributed.
- Before handling patients with the device, inform them of the manoeuvres to be performed.
- Before lifting, stand with your feet wide apart to provide a more stable base of support.
- No equipment, device or operator should ever pass over the patient.
- Avoid excessive force when lifting the device: unnecessary force can cause damage to the patient and the product itself and can adversely affect rescue operations.
- When transporting on stairs, make sure that the patient is facing the direction of travel.
- Hold the device firmly when moving patients.**
- Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise rescue operations. If you cannot clear the way, choose an alternative route.
- Do not alter or modify the product arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Before starting loading procedures on the sheet, remove any glasses or unstable dentures from the patient.
- Lift the device without excessive force but firmly. Any uncertainties during lifting could cause the patient's head to hit the ground.
- During all movements, take care to maintain adequate distance from the floor as the device does not protect the patient from potential impact.
- Do not use the device in or near MRI machines.**
- Even in its version with footrest, the device must not be used to transport the patient in an upright position.**

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in these instructions for use, have been identified.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

	HIRVIN		DESCRIPTION AND MATERIALS	
			Sheet	
1			Made of PVC, it is the main component of the device, providing support to the patient.	
			Side handles	
2			Made of nylon straps, these are an extension of the transverse reinforcement straps and are intended to be grasped by the operators when lifting patients.	
	Hirvin 6 handles	Hirvin 8 handles	Hirvin 10 handles	Hirvin Plus
Length	1850 ± 10 mm	1850 ± 20 mm	1850 ± 10 mm	1850 ± 20 mm
Width	760 ± 5 mm	760 ± 5 mm	760 ± 5 mm	760 ± 5 mm
Wrapped dimensions (1)	150x80x760 mm	150x80x760 mm	150x80x760 mm	150x80x760 mm
# handles	6	8	10	8
Strap width	50 mm	50 mm	50 mm	50 mm
Length of handles	170 mm	170 mm	170 mm	170 mm
Capacity	150 kg	150 kg	300 kg	150 kg
Materials	Nylon, PVC	Nylon, PVC	Nylon, PVC	Nylon, PVC
Weight	1,3 ± 0,2 kg	1,4 ± 0,2 kg	1,5 kg	1,5 ± 0,2 kg

(1) The measurement refers to wrapped sheets. This is purely indicative, as the device can be folded or wrapped in numerous ways

PARTS		DESCRIPTION AND MATERIALS	
1	Upper side strap attachment loop Made of nylon; the anchor point (head side) of the side strap	9	Side strap length adjustment loops Made of nylon, these allow quick adjustment of the length of the side strap
2	Patient handles Made of nylon, these offer the patient a gripping point, increasing stability during transport	10	Side strap attachment spring catch Made of aluminium, hooks the side strap to the sheet
3	Lower side strap attachment loop Made of nylon; the anchor point (foot side) of the side strap	11	Side strap length adjustment spring catch Made of nickel-plated brass, this allows the length of the side strap to be adjusted, attaching it to the respective loops.
4	Upper footrest attachment loop (if present) Made of nylon, is one of the attachment points of the footrest	12	Side strap belt Made of nylon, it distributes the load better on the operator during lifting instead of concentrating it exclusively on the lifting handles
5	Lower footrest attachment loop (if present) Made of nylon, is one of the attachment points of the footrest	13	Accessory pocket Pocket made in the sheet for reinforcement rod storage when they are not in use
6	Footrest attachment spring catch (if present) Made of aluminium, hooks the footrests to the respective loops.	14	Rod insertion point Seams in the sheet into which the reinforcement rods are to be inserted
7	Sheet fabric Made of PVC, it is the main component of the device, providing support to the patient.	15	Reinforcement rods Made of aluminium, they are jointed and are used to provide rigidity to the sheet
8	Handles Made of nylon, used to lift the patient	16	Elastic strap Made of reinforced elastic material, this allows the jointed portions of the rods to remain joined together.

SPECIFICATIONS		VALUE	
Device	WOW		WOW W/FOOTREST
Length	1870 mm		2560 mm
Upper part width		960 mm	
Lower part width		700 mm	
Dimensions when stored (1)		500x600x150 mm	
Dimensions of single folded rod		74x450x30 mm	
# Handles		10	
Strap width		50 mm	
Length of upper handles		160 mm	
Length of lower handles		260 mm	
Rod length		1210 mm	
Number of side strap adjustments		6	
Capacity		150 kg	
Weight	4,2 Kg		4,7 Kg

(1) The measurement refers to wrapped sheets. This is purely indicative, as the device can be folded or wrapped in numerous ways.

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- Product cleanliness
- Absence of cuts, holes, tears or abrasions on the entire device

Do not modify the device or any of its parts for any reason as this could cause damage to the patient and/or rescuers.

⚠ Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.

Practice with devices without patients before putting the device into service.

For subsequent use, perform the operations specified in paragraph 12.

If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

10. PROPER USE

Primary medical evaluations must be carried out before the patient is moved, lifted or transported. Once the diagnosis has been made, it is preferable to advise the patient to actively contribute to the transfer procedures, while at the same time informing them of the risks they may be running. Before loading the patient, bring the device as close as possible to the patient.

IT
EN
DE
FR
ES
PT
RO
CS
SK

10.1 PATIENT POSITIONING

Note: The procedures described below are intended to illustrate one possible way of using the product. Follow the guidelines approved by your emergency medical service when using the device.

In the case of the Wow injury transportation sheets, remove the rods and the side straps, opening the accessory pocket.

Lay the sheet next to the patient with the side with the yellow and black pattern facing up and the wider side facing the patient's head.

For Hirvin models, the side that should be facing up is the side with the smoothest PVC surface treatment.

Fold the sheet in half lengthwise.

The upper half of the sheet must be folded back to expose the yellow and black pattern.

Bring the sheet as close as possible to the patient lying on their back.

The operator on the folded side of the sheet helps the patient to turn on his/her side, thus turning his/her back to the sheet.


The operator on the opposite side pushes the sheet under the patient, having its midline coincide with the patient's spine.

After the patient has been returned to the supine position, he/she is turned onto the other side so that the free operator can grasp the handles of the sheet and spread it out. Return the patient to the supine position.

When using the Wow sheet, the log roll manoeuvre can also be performed with the rods already in place, if necessary.

10.2 ROD INSERTION (WOW ONLY)

 Open each rod, straightening it and inserting the narrowest part of each section into the next section until they are completely flush with each other.

Be careful when opening the rods: **Crushing hazard.**


Insert the rods into the slots at the top of the sheet. To lift the patient in a lying position, the rods should be inserted along their entire length.

In the case of tight spaces, the rods can also be inserted through the holes immediately below the transverse strap located in the middle of the sheet and visible from the patient's side.

IT

10.3 APPLYING AND ADJUSTING SIDE STRAPS (WOW ONLY)

The side straps supplied help to reduce fatigue for the rescuers by distributing the load on the shoulders as well.

1. The shoulder strap has one end with an integrated small spring catch (no. 11 par. 8), and one free end. Attach one of the supplied spring catches to the lower loop for attaching the side strap (no. 3 par. 8 and no. 1 in the picture opposite).
2. Insert the strap into the second spring catch and attach it to the upper loop about 85 cm from the lower loop (no. 2 picture opposite).
3. Then use the spring catch built into the side strap to secure it in the closed position, choosing the loop most suited to your height (no. 3 picture opposite).

The side strap can be adjusted to 6 positions by changing the attachment loop on the spring catch with loop.

In the picture, letter "A" indicates the loop for shorter operators. Letter "B" identifies the loop for taller operators. Each intermediate position corresponds to a variation in height of the operator of approximately 5cm.

10.4 USING THE FOOTREST (IF PRESENT)

The footrest is an element of the device that limits the longitudinal movement of the patient.

After checking that the patient is correctly positioned, fold the lower end of the sheet over the patient's legs.

Then hook the footrest spring catches to the lower or upper loop on the sides of the sheet (no. 4 and 5 par. 8). Choose the loop that best suits the patient's height.

10.5 LIFTING THE PATIENT IN A SUPINE POSITION

After the operations described above have been performed, the patient can be lifted.


The number of operators must be proportionate to the weight to be lifted in accordance with occupational safety regulations.

For Hirvin models, 4 operators should be placed on either side of the sheet so that they can safely grasp the handles and lift the sheet with the patient on it.

When using the Wow sheet, proceed as follows:

 Two operators will stand on either side of the sheet looking down at the patient's feet and will place the side strap on the shoulder opposite to the patient's side, thus distributing the weight over the entire spine.

Caution: Do not place the side strap in such a way as to concentrate the load on the neck.

SK

Ask the patient to grasp the red handles that will be at shoulder height by crossing their arms. In this way, he/she will be protected from the movements of the operators, will be more stable during transport and will feel safer.

The remaining operators then grasp the handles on the sheet that are appropriate for the type of transport and, after synchronising with each other, lift the patient.

10.6 LIFTING THE PATIENT IN A SEATED POSITION


In addition to transporting the patient in a supine position, the Wow injury transportation sheet also allows for transport in a seated position. This is particularly suitable in confined spaces and makes the device more manageable.

For this type of transport, only 2/3 of the internal rods must be inserted inside the sheet, while the remaining part must be kept outside the sheet and folded on the side opposite the patient, securing it with the red straps on the back of the sheet.

In order to ensure better load distribution and to prevent injury from strain, it is essential to support the weight also with the help of the handles (not only the side straps).

The operator can choose which handles to use according to his/her height and with the aim of minimising fatigue.

It is advisable to always grasp the middle handle of the sheet (located at the point where the sheet narrows). It is this handle that carries the most weight during transportation in a seated position.

10.7 ACCESSORY POCKET (WOW ONLY)

The device has a bag with access from the top of the sheet, which is kept closed with the help of two magnets.

A pillow or blanket can be placed inside during transport to increase patient comfort. After use, rods and shoulder straps can be inserted, ensuring easy retrieval of accessories and reduced storage space.

11. CLEANING AND MAINTENANCE

11.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

After complete drying, proceed with lubrication as described below.

If disinfecting, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices.

Carefully follow the instructions of the manufacturer of the product used with regards to the application method and contact time.

Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

11.2 ROUTINE MAINTENANCE

You must establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged within the instructions for use.

All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

The maintenance schedule must comply with the following table:

MINIMUM MAINTENANCE INTERVALS	AT EACH USE	IF NECESSARY	MONTHLY
Disinfection	•		
Cleaning		•	
Inspection	•	•	•

The inspection to be carried out after each use involves:

- Check that all components are present
- Check device conditions – There must be no cracks, holes or cuts
- Check conditions of wear – There must not be any level of abrasion that would compromise the safety of the product
- Check that the moving parts slide properly.
- General check of the state of wear of each component
- Check that all the accessories provided are present, functional and intact.
- Correct operation of the spring catches and rods
- Absence of cuts, holes, tears or abrasions on the entire product
- Intact conditions of side straps and handles

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

11.3 PERIODIC OVERHAUL

No periodic overhaul is foreseen for the device.

11.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

11.5 LIFE SPAN

The device, if used as described in the following instructions, has an average life span of 5 years from the date of purchase, beyond which it must be disposed of.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have exceeded the maximum allowable life span.

IT

EN

DE

FR

ES

PT

RO

CS

SK

12. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
Tears in the fabric	Wear or improper use	Remove the device from service and replace it
Spring catches do not work properly	Deformed spring catch	Immediately take the device out of service and replace the spring catches with those listed in the spare parts section
	Return spring breakage	
The sections of the internal rods cannot be joined together	The elastic strap is damaged	Immediately take the device out of service and replace the rods
	One or more sections are deformed	
The sheet cannot support the weight of the patient	Sheet damaged	Immediately take the device out of service and contact the service centre
Handles are cracked	Wear or improper use	Immediately take the device out of service and contact the service centre

13. ACCESSORIES

There are no accessories available for this device.

14. SPARE PARTS

RIQM001C	PAIR OF BLACK ALUMINIUM SPRING CATCHES	RIQM002C	PAIR OF WOW INTERNAL RODS
RIQM003C	PAIR OF WOW COMPLETE SIDE STRAPS		

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by organic agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

RO

CS

SK

Notice

The information in these instructions is subject to change without notice.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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