
Manuale d'uso utente

JET HOMECARE

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
Rappresentante Europeo:

DeVilbiss Healthcare GmbH
Kamenzerstraße 3
68309 Mannheim
Germany

User manual
JET HOMECARE
 Suction device for home use

CE 0044 Class IIa Medical Device according to Directive. 93/42/CEE

Warning

The information contained in this manual is subject to change without notice.

The images are inserted only for reference and may vary slightly from the actual device.

Spencer Italia S.r.l. assumes no responsibility for any errors contained herein or for damage, accidents or consequences connected with the supply, performance or use of this manual.

First edition: 2014
 Rev. 0: 10/12/2014

INDEX

1.	MODELS	13
2.	INTENDED USE	13
3.	REFERENCE STANDARDS	13
4.	INTRODUCTION	13
4.1	Symbols	13
4.2	Warranty	13
5.	WARNINGS	14
6.	TECHNICAL DATA AND COMPONENTS	17
6.1	Technical data	17
6.2	Components	17
7.	INSTRUCTIONS FOR USE	18
8.	CLEANING INSTRUCTIONS	19
9.	MAINTENANCE	20
10.	TROUBLESHOOTING	20
11.	ACCESSORIES/SPARE PARTS	20

Warning

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l. Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, lay-out or technical aspects that are herein described.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.

1. MODELS

JET COMPACT HOMECARE
JET WIRE HOMECARE

2. INTENDED USE

The device is to be used to remove fluids from the airway or respiratory support system and infectious materials from wounds. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection canister. The fluids are trapped in the collection canister for proper disposal. It is for use on the order of a physician only. To maximize product life and performance, follow recommended operating and maintenance procedures.


3. REFERENCE STANDARDS

Approvals

IEC 601-1; IEC 68; CAN/CSA-C22.2 No. 601.1-M90; UL 2601-1, CE EN 60601-1-2, ISO10079-1:1999 Meets RTCA/DO-160E DO-160E - section 21 Category M (for battery operation only; commercial aircraft, airborne equipment)

4. INTRODUCTION

4.1 Symbols

	General or specific warnings, consult instruction guide	~	Alternating current
	Standby - On/Off	 12 V DC	Direct current
	IP12: Vertically falling water drops shall have no harmful effects when the enclosure is tilted at an angle up to 15° on either side of the vertical.		Center positive polarity indicator
	Type BF equipment-applied part		This device contains electrical and/or electronic equipment that must be recycled per EC Directive 2012/19/UE- Waste Electrical & Electronic Equipment
REF	7310PD-S	Model number of the suction device main body	

4.2 Warranty

Two-Year Limited Warranty

The compressor portion of the Suction Unit (excluding internal rechargeable battery and collection canister) is warranted to be free from defective workmanship and materials for a period of two years from date of purchase. Internal rechargeable batteries are warranted for 6 months. Any defective part(s) will be repaired or replaced at Spencer's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your Spencer Provider or Spencer for instructions.

NOTE - Be sure to retain a dated proof of purchase document to verify unit is within 2-year warranty period.

NOTE - This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW. ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

5. WARNINGS

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

DANGER—Urgent safety information for hazards that will cause serious injury or death.

WARNING—Important safety information for hazards that might cause serious injury.

CAUTION—Information for preventing damage to the product.

NOTE—Information to which you should pay special attention.

READ ALL INSTRUCTIONS BEFORE USING.

SAVE THESE INSTRUCTIONS

DANGER

This suction unit is a vacuum suction device designed for the collection of nonflammable fluid materials in medical applications only.

Improper use during medical applications can cause injury or death. For all medical applications:

1. All suctioning should be done in strict accordance with appropriate procedures that have been established by a licensed medical authority.
2. Some attachments or accessories may not fit the tubing supplied. All attachments or accessories should be checked prior to use to assure proper fit.

To reduce the risk of electrocution:

1. Do not use while bathing.
2. Do not place or store product where it can fall or be pulled into a tub or sink.
3. Do not place in or drop into water or other liquid.
4. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING

To reduce the risk of burns, electrocution, fire or injury to persons:

1. Close supervision is necessary when this product is used by, on, or near children or physically incapacitated individuals.
2. Use this product only for its intended use as described in this guide.
3. Never operate this product if:
 - a. It has a damaged power cord or plug.
 - b. It is not working properly.
 - c. It has been dropped or damaged.
 - d. It has been dropped into water. Return the product to an authorized Spencer service center for examination and repair.
4. Keep the power cord away from heated surfaces.
5. Never use while drowsy or asleep.

International Travel

The suction unit is equipped with an AC to DC Adapter/Charger allowing operation on any AC voltage (100-240 VAC, 47-63 Hz). However the correct power cord must be used to connect to adaptable wall power. Your unit is approved for use on commercial aircraft; specific standards are listed in the Specifications in this guide. Certificate is available upon request; if needed, contact your provider for a copy of approval letter before travel.

NOTE—Check power cord for adaptability before using.

WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment. The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTE– The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems			
This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.			
Emissions Test	Compliance	Electromagnetic Enforcement – Guidance	
RF Emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class B Radiated and Conducted Emissions	This device is suitable for use in all establishments including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonics IEC 61000-3-2	Class A		
Flicker IEC 61000-3-3	Complies		
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/burst IEC 61000-4-4	±2kV on AC Mains	±2kV on AC Mains	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms from 150 kHz to 80 MHz	V1 = 3 Vrms	Portable and mobile RF communications equipment should be separated from the device by no less than the recommended separation distances calculated/listed below: $D=(3.5V1)\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	E1 = 3V/m	$D=(3.5/E1)\sqrt{P}$ 80 to 800 MHz $D=(7/E1)\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and this device. This device and system are NOT Life-Supporting

This device is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Maximum Output Power (Watts)	Recommended Separation Distances for the device (meters)		
	150 kHz to 80 MHz $D=(1.1667)\sqrt{P}$	80 to 800MHz $D=(1.1667)\sqrt{P}$	800 MHz to 2.5 GHz $D=(2.3333)\sqrt{P}$
0.01	0.11667	0.11667	0.23333
0.1	0.36894	0.36894	0.73785
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3785
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.


Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

6. TECHNICAL DATA AND COMPONENTS

6.1 Technical data

	JET COMPACT HOMECARE	JET WIRE HOMECARE
Dimensions (WxDxH)	245x70x185 mm	185x205x220
Weight	1,54 kg	1,98 kg
Power supply weight	520 g	
Electrical Requirements	100-240 VAC 50/60 Hz 0,75 A max; 12 VDC, 33 W max	
Vacuum Range	Da 50 a 550 mmHg	
Air Flow at pump inlet	27 LPM (free flow) typical (may be less when running from internal battery)	
Collection Canister Capacity	300 ml	1000 ml
Operating Temperature Range	0°C – 40 °C	
Operating Relative Humidity	0-95%	
Operating Atmospheric Pressure	10,2 Psi (70 kPa) – 15,4 Psi (106 kPa)	
Storage & Transport Temperature Range	-40 °C – 70 °C	
Storage & Transport Relative Humidity	0-95%	
Storage & Transport Atmospheric Pressure	7,3 Psi (50 kPa) – 15,4 Psi (106 kPa)	
Equipment Classifications		
With respect to protection from electric shock	Classe I e alimentazione interna	
Degree of protection against electric shock	Parti applicate di tipo BF	
Degree of protection against ingress of liquids	IP12 e alimentazione ordinaria	
Mode of Operation	Funzionamento intermittente: 30 minuti acceso, 30 minuti spento	
ISO Classification	7310 Series - Electrically powered medical suction equipment for field and transport use according to ISO 10079-1 : 1999	
	High Flow/High Vacuum	

6.2 Components

N°	Descrizione	N°	Descrizione
1	On/Off button	8	Connection for Catheter
2	Vacuum Regulator Knob	9	Disposable canister for Jet Compact
3	12V DC Power Input	10	Antibacterial filter
4	Display Panel	11	Connection for Catheter
5	Unit Carry Handle/Catheter Holder	12	Canister for Jet Wire
6	Vacuum Inlet Port	13	Wire Canister Bracket Jet Wire
7	Battery Door		

Set-Up

Battery Connection

Allestimento

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- Using a coin or straight-blade screwdriver rotate latch to unlocked position.
- Remove door by pulling up on latch.
- Remove battery from compartment and plug connector into circuit board.
- Replace battery and door; rotate latch to locked position.
- Fully charge battery for 5 hours before using.

Jet Wire

Place canister into wire canister bracket; ensure inlet port marked "Patient" is accessible.

- Attach one end of the connection tubing to connection elbow.
- Insert unit connection elbow into side port.
- Attach other end of connection tubing to side of filter marked "Out".
- The 90° canister elbow connects to the clear side of filter marked "In" and to the top of canister lid marked "Vacuum".

NOTE-Verify clear side of filter marked "In" faces canister.

- Connect patient tubing to canister lid at inlet port labeled "Patient".
- Ensure all connections are secure to prevent leakage in the canister/tubing system.

Jet Compact

- This single-use sealed disposable canister includes an internal filter and fluid shut-off that automatically stops suctioning when it becomes wet.
- Firmly attach canister by pushing connection fitting straight into the open side port of unit. (Both the top connection fitting and the bottom support tab secure the canister).
- Securely attach appropriate suction tip to the tubing.

7. INSTRUCTIONS FOR USE

Control Panel Symbols

	On/Off
	External power: Supplied from AC or DC power cord. Illuminates in GREEN when power is supplied.
	Battery charging: Illuminates in YELLOW. The light will go out when the battery is fully charged.
	Low battery: illuminates in RED when battery reaches a discharge state.
mmHg	Suction level setting: Scale/strength illuminates in GREEN. This scale shows the level of suction strength in mmHg.

Power Source Options

AC OPERATION - Plug the 90° power connector of AC to DC adapter/charger into DC power input and attach the line cord. Plug the other end of AC adapter into a grounded AC supply.

BATTERY OPERATION - Your unit is equipped with a high capacity rechargeable battery. For initial charge on new unit, fully charge the battery for a minimum of 5 hours (see Battery Charging). To operate unit from the rechargeable battery, ensure that no external power sources are plugged into the suction unit.

NOTE- During charging or operating, the power supply may become warm to touch; this is normal.

WARNING

If the unit does not receive external power or the battery has not been charged, the low battery indicator light will remain on and the performance of the unit will drop rapidly. If Low Battery symbol lights, switch to another power source to avoid an interrupted suction procedure.

How to adjust the vacuum level

1. Once power source is selected, turn the unit on by pressing the "On" button. The GREEN light, indicating external power, will remain lit when external power is connected.
2. Occlude (block) the patient end of the tubing, then adjust vacuum level from 50 to 550 mmHg by turning the vacuum regulator knob clockwise to increase and counter-clockwise to decrease the vacuum. Release and occlude once more to confirm setting. The desired level of vacuum can be viewed on the LED display.

NOTE-The LEDs have two brightness levels. As the vacuum level is adjusted, the LEDs will illuminate in progression. When an LED is at half brightness, it indicates that the vacuum level is halfway between the previous fully lit LED and the half brightness LED. EXAMPLE: If the 150

mmHg LED is fully illuminated and the 200 mmHg LED is at half brightness, the suction level is 175 mmHg. When the 200 mmHg LED illuminates at full brightness, the unit has reached 200 mmHg.

3. Connect suction tip or catheter as appropriate.

NOTE—If the unit does not maintain vacuum, refer to *Troubleshooting*.

NOTE—Suction ceases if liquid level reaches float shut-off valve located on underside of canister lid.

CAUTION—Further suctioning may cause damage to the vacuum pump and voids warranty. Equipment service is required if fluid content is aspirated back into the unit.

Battery Charging

The device is equipped with a factory-installed high capacity rechargeable battery. Connect battery and charge before first use (refer to *Battery Replacement*).

1. Connect the unit to either an AC or DC power source.

2. The green external power light should illuminate; the yellow charge indicator will remain lit while the battery is charging.

3. Ensure that the yellow charging light is illuminated when charging begins. As the battery nears a full charge, the yellow battery charging light may flash for several minutes. This is normal. If your unit does not hold a charge, check that the yellow light turns on when external power is applied with the power button "Off." If problems persist, contact your medical equipment provider.

NOTE—A fully discharged battery will require up to 5 hours (depending on depth of discharge) of charging to reach full capacity. If unit is not in use for extended periods, the battery should be recharged a minimum of every 6 months. A fully charged battery will provide approximately 45-60 minutes of continuous operation at zero vacuum level (free flow).

CAUTION—Completely discharging the battery will shorten the battery life. Do not operate the unit for more than a few minutes if the low battery indicator light is lit. Recharge the battery as soon as possible.

Battery replacement

1. Using a coin or straight-blade screwdriver rotate latch to unlocked position.

2. Remove door by pulling up on latch.

3. Remove battery from compartment and unplug connector from circuit board.

4. Install new battery by reversing the above steps.

5. Fully charge battery before using.

8. CLEANING INSTRUCTIONS

Collection Canister (not applicable for Jet Compact):

NOTE—Elbow is permanently attached to the lid and cannot be removed.

1. Push power button to turn unit off. Wait for vacuum level to drop.

2. Disconnect external power source from input receptacle on unit (if applicable).

3. Remove canister from unit by disconnecting the patient tubing and filter.

NOTE—The reusable collection canister and lid are single-patient use and should be replaced every 1 to 2 months, depending on use and cleaning method.

4. Collection canisters and lids should be emptied and thoroughly cleaned after each use. Wash in a hot water/dishwashing detergent solution and rinse with clean, hot tap water. Then rinse with a commercial disinfectant or one part vinegar to three parts hot water solution. Thoroughly rinse with hot tap water and air dry. **Note:** If using a commercial disinfectant, follow manufacturer's recommended instructions and dilution ratios carefully.

Suction Unit:

1. With the power "Off," disconnect the unit from all external power sources.

2. Wipe the outside housing with a clean cloth dampened with commercial disinfectant/detergent (bacterial-germicidal).

CAUTION—Do not submerge suction unit in water; this will result in damage to vacuum pump.

Tubing:

1. Disconnect tubing from unit.

2. Tubing should be rinsed thoroughly after each use by running tap water through it, followed by a solution of one part vinegar to three parts hot water.

3. Rinse again with hot tap water and air dry.

4. Keep the outside surface of the tubing clean by wiping with a clean damp cloth.

9. MAINTENANCE

Changing Filter and Replacing Canister

1. Change the bacteria filter every 1 to 2 months OR immediately if overflow occurs.
2. Remove the bacteria filter by disconnecting it from the suction unit and lid assembly.
3. Replace it with a new bacteria filter.

NOTE – Verify clear side of filter marked “In” faces canister.

NOTE – Use only the bacteria filter provided by Spencer or one of its Distributors. Substitution may lead to contamination of the unit and/or poor performance and will void warranty.

4. Replace Collection Canister every 1 to 2 months, depending on use and cleaning method.

10. TROUBLESHOOTING

NOTE– Before you contact your equipment provider, follow Troubleshooting below:

DANGER

Electric shock hazard. Do not attempt to open or remove cabinet, there are no user-serviceable internal components. If service is required, return the suction unit to a qualified Spencer provider or authorized service center. Opening or tampering with the unit will void the warranty.

NOTE– If problem is not resolved, contact your authorized equipment provider.

Problem	Action
Unit does not power on. (Green external power indicator should be illuminated when power is applied.)	<ol style="list-style-type: none">1. Check power sources and connections.2. Ensure wall outlet is live by plugging in a lamp.3. Check that battery is fully charged.
Pump runs, but no vacuum.	<ol style="list-style-type: none">1. Check that all tubing is connected properly.2. Check tubing connections for breaks, leaks, or occlusions.3. Ensure that float shut-off is not activated due to full bottle.4. Check for leaks or cracks in bottle assembly.
Low vacuum.	<ol style="list-style-type: none">1. Use vacuum adjustment knob to increase vacuum level.2. Check system for leaks.
Battery will not hold a charge.	<ol style="list-style-type: none">1. Confirm battery is connected (refer to Set-Up).2. Verify that charge light turns on.3. Check electrical connections during charging.4. Ensure wall outlet is live by plugging in a lamp. If problem continues, contact your equipment provider.

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DeVilbiss Healthcare LLC