## SED - Estricatore / SED - Extrication device / SED - Rettungskorsetts / SED - Attelle d'extraction / SED – Dispositivo de extricación / SED - Dispositivo de extração / SED - Γιλέκο απεγκλωβισμού / SED - Гърбодържател / SED - Vyprošťovací korzet

1

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**EN** Release the groin straps from the back by letting them drop to the sides and wrap the thoracic flaps around the patient's torso, passing them under the arms. Adjust the position of the device by lifting it so that the chest flaps become a support for the armpits. During the application, operators should avoid moving the

patient's head to avoid aggravating any injury. Check that the device has been correctly positioned, i.e. as close as possible to the patient's back and well aligned with the spine.

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**EN** Fasten the chest belts in the following sequence. Generally accepted as correct:

- 1 – Intermediary (Yellow)
- 2 – Lower (Red)
- 3 – Upper (Green)

This sequence may not be indicated in all rescue situations. Follow the instructions of your local health authority. When applying it to pregnant patients, some of the thoracic flaps have to be folded to leave the abdomen uncovered. A possible closure configuration, which may

vary depending on the size of the abdomen and breasts, is as follows:

- Intermediary belt in the respective buckle
- Lower belt in upper buckle
- Upper belt in the lower buckle

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**EN** Close the groin straps by sliding them under the patient's legs using any free space left by the seat.

Check that the belts are straight in relation to the anchorage point and as close as possible to the symmetry axis of the device. Hook the belts into their buckles and tighten them until they come into contact with the patient. A further tightening must be carried out before extrication.

The use of crotch belts requires particular care in the case of femur or pelvis fractures. The device handles can help obtain better vertical positioning and alignment of the device.

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**EN** Do not immobilise the head until checking the correct position of the extrication device in relation to the size of the patient. Take care not to cause head and neck movements that could compromise the patient's condition. If there is a gap between the extrication device and the patient's head, collar or shoulders, fill the nape space with the appropriate support cushion, making sure not to hyper extend or flex the patient's neck.

While one operator holds the patient's head in position and positions the flaps, the other operator must secure the headband by tilting the back slightly downwards. This reduces the risk of the headband slipping over the head.

Then apply the chin strap, taking care to support the cervical collar and tilt it towards the patient's ears. This avoids immobilising the jaw and preventing the mouth from opening. Ensure that the straps are correctly secured and symmetrically positioned.

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**EN** When the stretcher or spine board has been brought to the extrication site, the final steps of the application can be carried out, i.e. further tightening of the extrication straps, therefore:

- 1 – Tighten the upper belt (green)
- 2 – Check and tighten the lower belt (red)
- 3 – Check and tighten the groin belts.

After checking the correct application, the two operators can grasp the handles of the device to rotate the patient so that his or her back faces the inside of the vehicle. Rotation must take place without twisting the spinal cord or the cervical spine.

The two operators can then grasp the nearest handle of the extrication device with one hand and the back of the patient's knees with the other. After coordinating the operation correctly, they can lift the patient.

If needed, a third operator could apply a spine board. The upper belt (green) must be loosened immediately after extrication to allow for

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