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IT

Instructions for use
VACUUM MATTRESSES

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EN

1. MODELS	11
2. INTENDED USE	11
2.1 INTENDED USE AND CLINICAL BENEFITS	11
2.2 TARGET PATIENTS	11
2.3 PATIENT SELECTION CRITERIA	11
2.4 CONTRAINDICATIONS AND SIDE EFFECTS	11
2.5 USERS AND INSTALLERS	11
3. REFERENCE STANDARDS	11
4. INTRODUCTION	11
4.1 USING THEINSTRUCTIONS FOR USE	11
4.2 DEVICE LABELLING AND TRACEABILITY	11
4.3 SYMBOLS	11
4.4 WARRANTY AND SERVICE	12
5. WARNINGS/DANGERS	12
6. SPECIFIC WARNINGS	12
7. RESIDUAL RISK	13
8. TECHNICAL DATA AND COMPONENTS	13
9. COMMISSIONING	14
10. PROPER USE	14
11. CLEANING AND MAINTENANCE	15
11.1 CLEANING	16
11.2 ROUTINE MAINTENANCE	16
11.3 PERIODIC OVERHAUL	16
11.4 SPECIAL MAINTENANCE	16
11.5 LIFE SPAN	16
12. TROUBLESHOOTING TABLE	16
13. ACCESSORIES	17
14. SPARE PARTS	17
15. DISPOSAL	17

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- COMBIMATT
- NEXUS
- RES-QMATT
- CDK
- RES-Q-MATT PLUS
- MINI MATT

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport. It is not foreseen that the patient be able to intervene on the devices.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group. The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity and within the limits of the size of the device. If paediatric subjects must be transported, it will be in the role of the rescuer to determine whether the devices are suitable for immobilization or if it will be necessary to use another device.

2.3 PATIENT SELECTION CRITERIA

The patients foreseen are those with injuries that prevent them from walking in a given rescue situation and requiring immobilization.

2.4 CONTRAINDICATIONS AND SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the instructions for use.

2.5 USERS AND INSTALLERS

The intended users are persons trained in first aid procedures and in use of medical equipment in an emergency medical service (EMS) environment. These devices are not intended for lay people.

Vacuum mattresses are devices intended for professional use only. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others. Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly. The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs, should it be necessary to lift and/or support the device and the patient. Operators' ability must be assessed before the definition of roles in use of the stretcher.

Operators must be able to provide the necessary patient care.

2.5.1 USER TRAINING

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of these instructions for use before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses.

2.5.2 INSTALLER TRAINING

The installer must be able to ensure safe storage inside the ambulance. The device in itself does not require installation.

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 USING THE INSTRUCTIONS FOR USE

The purpose of these instructions is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: The instructions for use are an integral part of the device and therefore must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer product instructions for use can be downloaded from the site <http://support.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of these instructions before operating or servicing this product.

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it>, or inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		See the instructions for use
	Medical device		Lot number

 Manufacturer	 Product code
 Date of manufacture	 Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
 Unique Device Identifier	 Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner (only for USA Market)
	Production identification Alphanumeric code that identifies the production units of the device, composed of:
(01)8057711230006 (11) 200626 (10) 1234567890	(01)805771123 company prefix 000 progressive GS1 6 control number (11)200626 date of production (YYMMDD) (10)1234567890 lot/SN

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it. To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself. Warranty and service conditions are available at <http://support.spencer.it>

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal (<http://service.spencer.it>) available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

5. WARNINGS/DANGERS

 Warnings, dangers, notes, and other important safety information are provided in this section and are clearly visible throughout the instructions for use. .

Product features

 Use of the product for any purpose other than that described in the instructions for use is prohibited.

- Before each use, always check the conditions of the product, as specified in the instructions for use. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
-  If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The product must not be tampered with or modified without the manufacturer's authorisation (modification, tweaking, additions, repair, use of non-approved accessories), as they may constitute imminent danger of injury to persons and material damage. Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE marking and the product warranty shall be null and void.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -10°C to +60°C.

Regulatory requirements

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the instructions for use.
- **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

The user must carefully read the following in addition to the general warnings.

-  It is not foreseen that application of the device lasts longer than the time required for first aid operations and subsequent transport to the nearest rescue point.
-  Qualified personnel and at least four operators must be present during use of the device.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Follow the internal procedures and protocols approved by your organisation.
-  Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.
- Do not use drying machines to dry the device.

6. SPECIFIC WARNINGS

To use the device, you must also have read, understood and carefully follow all the instructions in the instructions for use.


-  Always comply with the maximum capacity, if any, indicated in the instructions for use. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. In addition, the operator should assess whether the overall size of the patient reduces the functionality of the product.
- If foreseen for the device, make sure that operators are in good physical condition before lifting, as listed in the instructions for use.
- Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with routine maintenance of the device must ensure the basic requirements envisaged by the manufacturer within these operating instructions.
- All maintenance activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the manufacturer when requested.
- Never leave the patient on the device unsupervised, as they could get injured.
- Avoid contact with sharp objects.
- Follow approved Emergency Medical Service procedures for patient immobilization and transportation.
-  Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Before lifting, make sure that operators have a secure grip on the device.
- **Perform rescue simulations with a mattress and a patient simulating load and accessories before putting the device into service.**
-  At least 4 operators in suitable physical conditions are required for use of the device. They must be endowed with strength, balance, coordination, common sense and must be trained on the correct operation of the device.
- For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of more than two operators is recommended in addition to the minimum foreseen operators.
- Before each use, always check the conditions of the device and its components, as specified in the instructions for use. In case of faults or damage that may compromise the functionality and safety of the device, and therefore the patient and the operator, remove the device from service or replace the components that are not intact.
- Do not lift if the weight is not properly distributed.
- Use the devices only as described in these instructions for use.
-  Do not alter or modify the device to adapt it to unforeseen conditions of use: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- The device is expected to come into contact through the patient's clothing. In case of direct contact with the skin, place a protective surgical drape over it to avoid contamination from substances that may have contaminated the device.
-  During lifting, the mattress must be kept horizontal with respect to the ground. Any abnormal inclination can cause serious damage to the patient, the device and the operator.
-  To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- Always observe the maximum applicable static safety load indicated in these instructions for use. Maximum static load is intended as a mechanical force applied slowly, not quickly, beyond which the device may not be safe. This value does not take into account the dynamic forces to be added to the static load, such as shocks, vibrations and possible weather and climate conditions during use of the device.
- Never leave the patient unsupervised when the device is in use, as they could get injured.
- After washing, the device and all its components must be left to dry completely before storage away from sunlight and direct heat sources.
- Do not machine wash the device.
- Avoid contact with sharp objects.
- Never use solvents or stain removers.
- Do not use the device if cuts, burns, abrasions, open seams or fraying are present.
- Do not drag the device on rough surfaces.
- Always check the conditions of all parts of the straps and buckles before each use.
- Immediately replace devices with worn or damaged straps or buckles.
- Position and adjust the straps and buckles in such a way that they do not hinder rescuer operations or use of the use of rescue equipment.
- Keep the appropriate documentation for a period of ten years from the date of transfer to the final consumer and, therefore, show, where required, to trace the origin of the products.
- Use of the device without proper patient immobilization can result in serious damage. Always make sure that the patient is properly immobilized before handling the device.
- The rescuer must assess the actual need for use of this type of device in accordance with local guidelines, ensuring that immobilization with vacuum devices is appropriate in relation to the clinical condition of the patient.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in these instructions for use, have been identified.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

DESCRIPTION AND MATERIALS

Strap handles		Mattress	
1	Made of Nylon, allow lifting of the mattress and the patient.	6	Made of PVC with polystyrene filling, this is the part of the device subjected to vacuum.

2	Valve Equipped with a closing cap (Res-q-Matt closing cap), this is the element through which the vacuum is carried out and which must be closed in order to maintain the level of pressurisation obtained.	7	Belt slots Where present, these are slots cut into the PVC fabric that allow the belts to pass through, improving stability by limiting unwanted movement.
3	Wooden inserts Made of wood, they stiffen the support and limit the shrinkage of the mattress when in a vacuum.	8	Removable mattress in PE Made of polyethylene foam, this is only available in Combimatt and is inserted between the vacuum mattress and the additional PVC fabric on which the carrying handles are made.
4	Handles Made from additional PVC fabric under the mattress, they are made from a cord around which the PVC fabric is wrapped.	9	Rubber sleeves These cover some sectors of the perimeter cord allowing them to be used as lifting handles.
5	"Hook and loop" coupling closure Made with Hook and loop inserts, it joins the two parts of the Combimatt mattress. The mattress (No.8) is inserted between the two coupled parts.	10	Perimeter cord Made of polypropylene, this allows the mattress to be grasped for lifting at the sleeves and allows better application of the belts in the remaining portions.

	Res Q Matt	Res Q Matt Plus	CDK	Nexus	Combi Matt	Mini Matt
1 Strap handles	•					•
2 Valve	•	•	•	•	•	•
3 Wooden inserts				•		
4 Handles		•		•	•	
5 "Hook and loop" coupling closure	•				•	•
6 Longitudinal channel mattress				•		
7 Belt slots		•		•	•	
8 Removable mattress in Polyethylene					•	
9 Rubber sleeves			•			
10 Perimeter cord			•			

	Combimatt	Nexus	Res-q-matt
Length	2030 ± 10mm	2025 ± 10mm	2060 ± 20mm
Width	990 ± 10mm	855 ± 10mm	800 ± 10mm
# Handles	8	8	8
Materials	PVC, PE, Nylon, Polystyrene, Brass	PVC, Wood, Nylon, Brass, Polystyrene	PVC, Nylon, Polystyrene
Weight (not including pump)	6,4 ± 0,5 kg	5,6 ± 0,5 kg	7,8 ± 0,5 kg
Load capacity	150 kg	150 kg	150 kg


	Res-q-matt plus	CDK	Mini Matt
Length	2140 ± 20mm	2150 ± 30 mm	1300 ± 30 mm
Width	940 ± 10mm	900 ± 30 mm	600 ± 30 mm
# Handles	8	9	4
Materials	PVC, Nylon, Brass, Polystyrene	PP, PVC, Nylon, Brass, Polystyrene	PVC, Nylon
Weight (not including pump)	5 ± 0,5 kg	4,1 ± 0,3 kg	2,2 ± 0,3 kg
Load capacity	150 kg	150 kg	30 kg

PUMPS	Materials	Dimensions	Weight	TUBE LENGTH
QMX109	ABS, PE	215x100x445 ± 5 mm	0,6 ± 0,1 kg	100 ± 5 cm
QMX120	Al/PVC/STEEL	500x140x80 mm	1,14 ± 0,1 kg	100 ± 5 cm
QMX125	Al/PVC/STEEL	500x140x80 mm	1,2 ± 0,1 kg	100 ± 5 cm

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- Product cleanliness
- Absence of cuts, holes, tears or abrasions on the entire device
- If the mattress model in use is the Res-q-matt, the adapter supplied with the mattress must be fitted to the pump tube prior to operation. This is done by squeezing the end of the pump tube as much as possible and then inserting the end of the pump tube into the adapter tube fitting. The adapter tube will be the one that will need to be connected to the mattress valve during use. Do not operate the mattress if this has not been done or if the connection is not stable.

See paragraph 10 for how to carry out the above-mentioned checks.

Do not modify the device or its parts for any reason as this could cause damage to the patient and/or rescuers.

⚠ Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself. For subsequent use, perform the operations specified in paragraph 11.



If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.


Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.



10. PROPER USE


Primary medical evaluations must be carried out before intervening on the patient.

Before using the vacuum mattress, you should carefully read the instructions for use for all accessories that you plan to use with the device, such as belts, cervical collars and spine boards. Strictly follow the local emergency service guidelines before placing the patient on the mattress. Carefully assess the need for a cervical collar or other immobilization aids in addition to the necessary restraint belts. Place the belts flat on the floor, spaced as far apart as you plan to apply them to the patient. Place the mattress firmly on top of the belts, on a flat surface, away from sharp objects such as glass, sheet metal or anything else that could damage the fabric of the mattress. If you have the Combimatt model, you can insert a spine board into the dedicated compartment to further stiffen the system.

Metal valve	Res-q-matt valve	
		<p>Open the mattress valve to allow air to pass through, thus making the mattress mouldable. Evenly distribute the filler material, thus avoiding localised accumulations which could reduce the performance of the device. The patient, who usually undergoes recovery procedures with a spinal board or pick-up stretcher before being placed on the mattress, should be positioned in the centre of the mattress to allow adequate immobilization. Once on the mattress, the spine board or pick-up stretcher used for recovery can be removed. During all operations, always maintain immobilization of the cervical spine, also using a cervical collar if appropriate.</p>

	<p>Connect one side of the tube to the suction pump. In QMX120, there is only one hose connector at the bottom of the pump. The QMX125 has two hose connectors. The one on the bottom side should be used as it is intended for vacuums and therefore suitable for use with vacuum mattresses. The hose connector on the side of the pump is for inflation and therefore should not be used with vacuum mattresses. The QMX109 model, which can be used for both inflating and creating vacuums, has the connections clearly marked on the pump. On this model, the tube must be connected to the vacuum connection port.</p>
----------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------


Metal valve	Res-q-matt valve	
		<p>The other end of the tube should be connected to the previously opened mattress valve. Ensure that the tube is correctly attached to the valve fitting as well as the Res-q-matt adapter.</p>

	<p>If using the QMX120/125 pumps, place a foot inside the bottom pump support. If using the QMX109 pump, place two feet on either side of the pump, above the dedicated supports. Two operators mould the mattress to the patient's body and fasten the belts. While the two operators keep the mattress snug against the patient, the third operator begins to create a vacuum in the mattress by repeatedly raising and lowering the pump plunger until the mattress is rigid enough to properly immobilize the patient. The mattress valve should be closed as soon as the necessary vacuum is achieved. Delayed closing of the valve will result in rapid loss of vacuum in the mattress. Re-tighten the belts around the patient in order to maintain an adequate level of immobilization.</p>
-----------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

APPLYING IMMOBILIZATION ACCESSORIES

It is essential to fit belts on the patient.
If the patient's condition requires and/or allows it, apply a cervical collar until an adequate degree of immobilization is achieved.
Make sure that the patient is correctly immobilized before moving him/her with the mattress.

LIFTING THE PATIENT

	<p>Never use cranes, slings or other systems to lift the device. Lifting is only permitted by hand and by an appropriate number of operators (minimum 4). Operators must be positioned symmetrically and in such a way to allow safe, stable lifting of both the head and foot sides of the spine board mattress. If the 4 minimum operators foreseen (illustrated in red) is not sufficient to ensure the safety of operations, two additional operators (illustrated in blue) must hold the device in the central area in order to distribute the load as best as possible. The maximum load on each operator must never exceed the maximum load permitted by work safety regulations.</p>
-------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

11. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the instructions for use, identifying a reference person who meets the basic requirements set forth in the instructions for use.
- **The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.**
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l.

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declines any responsibility for any direct or indirect damage which is a consequence of improper use of spare parts and/or any repair work carried out by unauthorised parties.

- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the instructions for use in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- If required, the product and all its components must be washed and left to dry completely before storage.

11.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. Never use solvents or stain removers.

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

After complete drying, proceed with lubrication as described below.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices.

Carefully follow the instructions of the manufacturer of the product used with regards to the application method and contact time.

Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

11.2 ROUTINE MAINTENANCE

You must establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged within the instructions for use.

All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

The maintenance schedule must comply with the following table:

Minimum maintenance intervals	At each use	If necessary	Monthly
Disinfection	•		
Cleaning		•	
Inspection	•	•	•
Check vacuum tightness	•		

The inspection to be carried out after each use involves:

- Check that all components are present
- Check device conditions – There must be no cracks, holes or cuts
- Check conditions of wear – There must not be any level of abrasion that would compromise the safety of the product
- Check that the moving parts slide properly
- General check of the state of wear of each component
- Check that all the accessories provided are present, functional and intact.

Vacuum tightness test consists of bringing the mattress to a vacuum and then simulating application. The mattress must remain rigid enough to allow immobilization for at least 30 minutes.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

11.3 PERIODIC OVERHAUL

No periodic overhaul is foreseen for the device.

11.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

11.5 LIFE SPAN

The device, if used as described in the following instructions, has an average life span of 5 years from the date of purchase, beyond which it must be disposed of.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have exceeded the maximum allowable life span.

12. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The device stiffens during a drop in altitude	Change in relative pressure	Check the actual pressure status of the device by opening the valve or depressurising the pump. If the problem persists, remove the device from service and contact customer service.
The device tends to lose the set vacuum and regenerate a pressure level equal to the ambient value inside it.	Valve deteriorated.	Immediately remove the device from service and contact customer service.
	Possible holes in the internal chamber	
No vacuum in the device	Inefficiency of the pump or wear/breakage of the connection system between the valve and the surface of the device	Immediately remove the device from service and contact customer service.
Device does not adapt to the patient's shape	Non-conforming device material Limited flexibility	Immediately remove the device from service and contact customer service.
Damage (holes, cuts and abrasions) to the external containment structure	Improper use. Incorrect storage.	Immediately remove the device from service and contact customer service.

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

13. ACCESSORIES

QM22803	TRANSPORT BAG FOR COMBIMATT	QM22902	TRANSPORT BAG FOR NEXUS
QM22170	QMX170 – BAG FOR ORANGE PVC VACUUM MATTRESSES	ST00597	STX597 – 1PC BELT WITH PLASTIC HOOK
ST00593	K BELT 3 – SET OF 3 BELTS STX 597 W/BAG		

14. SPARE PARTS

RIQM005	SPARE VALVE FOR VACUUM MATTRESSES	RIQM008	QMX109 – DOUBLE EFFECT ABS PUMP
RIQM009	QMX120 – ALUMINIUM VACUUM PUMP	RIQM010	QMX125 – ALUMINIUM DOUBLE EFFECT PUMP

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

EL

Notice

The information in these instructions is subject to change without notice.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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