Istruzioni d'uso
Fixo Kid

IT

Instructions for use
Fixo Kid

EN

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DE

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EN

1. MODELS	10
2. INTENDED USE	10
2.1 INTENDED USE AND CLINICAL BENEFITS	10
2.2 TARGET PATIENTS	10
2.3 PATIENT SELECTION CRITERIA	10
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	10
2.5 USERS AND INSTALLERS	10
2.5.1 FORMAZIONE UTILIZZATORI	10
2.5.2 INSTALLER TRAINING	10
3. REFERENCE STANDARDS	10
4. INTRODUCTION	10
4.1 USING THE INSTRUCTIONS FOR USE	10
4.2 DEVICE LABELLING AND TRACEABILITY	10
4.3 SYMBOLS	11
4.4 WARRANTY AND SERVICE	11
5. WARNINGS/DANGERS	11
6. SPECIFIC WARNINGS	12
7. RESIDUAL RISK	12
8. TECHNICAL DATA AND COMPONENTS	12
9. COMMISSIONING	12
10. PROPER USE	13
11. CLEANING AND MAINTENANCE	14
11.1 CLEANING	14
11.2 ROUTINE MAINTENANCE	14
11.3 PERIODIC OVERHAUL	14
11.4 SPECIAL MAINTENANCE	14
11.5 LIFE SPAN	14

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RO

1. MODELS

Only one Fixo Kid model is available.

The device may be subject to implementation or modification without prior notice.

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Fixo Kid is a restraint system for paediatric patients to be used to obtain a level of restraint that allows the patient to be transported safely when used in accordance with the product in use.

2.2 TARGET PATIENTS

The target patients are paediatric patients weighing between 4 and 20 kg.

Further limitations of use may be related to the device with which the belt is used. For this purpose, please refer to the instructions of the transport device in use.

2.3 PATIENT SELECTION CRITERIA

The target patient selection criteria are those applicable to the device with which the paediatric restraint system is used.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the instructions for use.

2.5 USERS AND INSTALLERS

The intended users are rescue workers and experts in patient handling, immobilisation and transport.

These devices are not intended for lay people.

Fixo Kid is a device intended exclusively for professional use. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

Operators using the device should be physically able to use the device and have good muscle coordination. Operators' ability must be assessed before the definition of roles in use of the stretcher.

Operators must be able to provide the necessary patient care.

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- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of these instructions for use before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: *Spencer Italia S.r.l. is always available for training courses.*

2.5.2 INSTALLER TRAINING

Installation is not required for this device.

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 USING THE INSTRUCTIONS FOR USE

The purpose of these instructions for use is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: *The instructions for use are an integral part of the device and therefore must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.*

Spencer products User Manuals can be downloaded from the site <http://support.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of these instructions before installing, operating, or servicing this product.

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). **This must never be removed or covered.**

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

If the assigned Lot/SN cannot be traced, the device must be reconditioned, provided only under the responsibility of the manufacturer.

EU Regulation 2017/754 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it>, or inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the instructions for use
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal law restricts the sale of this device by or on the order of a licensed professional (US market only)
		Production identification Alphanumeric code that identifies the production units of the device, composed of: (01)805771123 company prefix 000 progressive GS1 6 control number (11)200626 date of production (YYMMDD) (10)1234567890 lot number	
(01)8057711230006 (11) 200626 (10) 1234567890			

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.
 For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.
 To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.
 Warranty and service conditions are available at <http://support-spencer.it>
Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

5. WARNINGS/DANGERS

Warnings, dangers, notes, and other important safety information are provided in this section and are clearly visible throughout the instructions for use.

At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website www.spencer.it on the specific product page.

Product features

Use of the product for any purpose other than that described in the instructions for use is prohibited.

- Before each use, always check the conditions of the product, as specified in the instructions for use. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The product must not be tampered with or modified without the manufacturer's authorisation (modification, tweaking, additions, repair), as they may constitute imminent danger of injury to persons and material damage. Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE marking and the product warranty shall be null and void.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to +50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the device.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -10°C to +60°C.

Regulatory requirements

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the instructions for use.
- **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

The user must carefully read the following in addition to the general warnings.

- The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point.

IT
EN
DE
FR
ES
PT
RO

- Qualified personnel and at least two operators must be present during use of the device.
- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Follow the internal procedures and protocols approved by your organisation.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Do not use drying machines to dry the device.
- In case of exposed and/or injured skin, cover the surfaces in contact with the patient with a surgical sheet that respects bio-compatibility regulations to protect the patient's health.

6. SPECIFIC WARNINGS


- To use the device, you must also have read, understood and carefully follow all the instructions in the instructions for use.
- Always comply with the positioning specifications determined by the device with which the Fixo Kid is used. Moreover, the operator should assess whether the overall size of the patient reduces the functionality of the product.
 -  Carry out rescue simulations with the device applied to the product with which it is intended to be used and a load simulating the patient and accessories, before putting the device into service.
 - Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with routine maintenance of the device must ensure the basic requirements envisaged by the manufacturer within these operating instructions.
 - All maintenance activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the manufacturer when requested.
 - Avoid contact with sharp objects.
 - Follow approved Emergency Medical Service procedures for patient positioning, restraint and transport.
 - Before each use, always check the conditions of the device and its components, as specified in the instructions for use. In case of faults or damage that may compromise the functionality and safety of the device, and therefore the patient and the operator, remove the device from service. .
 -  Do not arbitrarily alter or modify the device to adapt it to unforeseen conditions of use: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
 -  To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
 - Make sure that the buckles are properly anchored to the strap and that they are intact and functional.
 - Never leave the patient unsupervised when the device is in use, as they could get injured.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in these instructions for use, have been identified.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

	Id	DESCRIPTION
	1	Main sheet
	2	Strap for attachment to stretcher backrest
	3	Shoulder belt
	4	Shoulder belts adjustment system
	5	Chest belts
	6	Strap for attachment to stretcher frame
	7	Adjustable padded crotch belt
	8	Groin belt adjustment system
	9	Stretcher frame buckle
	10	Belt buckle closure
	11	Half buckle central hook group
	12	Chest belt tension adjustment system
	13	Shoulder belt padding

	FIXO KID
Length when extended (mm)	800 ± 10
Width when extended (mm)	485 ± 5
Soft padding thickness	5 mm
Dimensions when placed in bag (mm)	430x330x80
Weight	2 ± 0.2 kg
Patient retention belts	5
Device fastening straps	3
Minimum load capacity (kg)	4
Maximum load capacity (kg)	20
Main cloth material	PVC
Belts material	Nylon
Hooks material	Stretcher attachment hooks: steel Patient fastening hook: POM

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device

- Product cleanliness
- Absence of cuts, holes, tears or abrasions on the entire structure
- Conditions of wear of the device.
- Check that the buckles are working properly.

Do not modify the device or its parts for any reason as this could cause damage to the patient and/or rescuers.

⚠ Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself. For subsequent use, perform the operations specified in paragraph 11.


If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

10. PROPER USE



Primary medical evaluations must be carried out before intervening on the patient.

Before using Fixo Kid, carefully read the instructions for the device with which it is to be used.



Before using the device, perform all the preliminary checks set out in these operating instructions.

1. Remove any standard straps from the stretcher.
2. Raise the stretcher backrest to an angle between 15 and 45 degrees.
3. Place Fixo Kid on the mattress with all straps extended.
4. Centre the device with respect to the mattress edges.
5. Encircle the stretcher's backrest with the stretcher strap and fasten the buckle. Leave a minimum amount of slack during this phase
6. Position the patient on the device

7. Reposition the strap around the stretcher back in such a way that it is perpendicular to the stretcher's perimeter frame.
8. Fasten the lateral straps (no. 6 section 8) to the stretcher frame over the side rail and adjust the tension so that they are of equal intensity on both sides and to ensure the device is not moving.
9. Fold the groin strap over the patient's abdomen, ensuring that the padded part is well positioned.
10. Pass the patient's arm through the gap between the right chest belt and the right shoulder belt (No. 5 and No. 3 section 8 respectively), then insert the half-buckle (No. 11 section 8) into the central buckle (No. 10 section 3.2).
11. Repeat the same operation for the left side.
12. Ensure that the two half-buckles are correctly coupled and inserted.
13. After checking that the padded parts are properly positioned to ensure the patient's comfort, adjust the tension of all the belts (shoulder, chest and groin) to achieve an effective retention system.

14. To remove the patient from the device, use the two red buttons on the half-buckles and remove the groin, chest and shoulder straps from the patient's body. Then lift the patient out of the device according to emergency procedures.

IT

EN

DE

FR

ES

PT

RO

11. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the instructions for use, identifying a reference person who meets the basic requirements set forth in the instructions for use.
- **The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.**
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l. declines any responsibility for any direct or indirect damage resulting from improper use of spare parts and/or any repair work carried out by unauthorised parties.
- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the instructions for use in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.

11.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Any metal parts exposed to external agents undergo surface treatments and/or coating in order to obtain better resistance. Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

In the case of disposable products, no cleaning is required except that the product be properly stored and packaged according to the manufacturer's specifications.

11.2 ROUTINE MAINTENANCE

You must establish a program for periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged within these instructions for use. All maintenance activities must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

The device does not require a routine maintenance program, but checks must be made before and after use to verify:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Absence of cuts, holes, tears on the structure, including the straps
- Correct tightening of belts
- Wear conditions
- Integrity of the stitching and cloths
- Integrity of the parts

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

11.3 PERIODIC OVERHAUL

No periodic overhaul is foreseen for the device.

11.4 SPECIAL MAINTENANCE

No special maintenance is foreseen for these devices.

11.5 LIFE SPAN

If used as described in the following instructions, the device has a lifetime of 2 years from the date of purchase, at the end of which the device must be disposed of.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have exceeded the maximum allowable life span.

12. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
It is not possible to hook the plastic buckle.	The buckle has not been correctly inserted.	Insert the buckle, ensuring that the two halves of the male part are correctly aligned with the female part. The release button must be able to fit into the seat of the female buckle.
The metal buckle does not hook.	Internal components damaged.	Put the device out of service immediately and replace it with a similar one.
Fabric damage	Improper use or wear	Take the device out of service immediately and replace it with a similar one.

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

13. ACCESSORIES

There are no accessories for these devices.

14. SPARE PARTS

There are no spare parts available for these devices.

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. Spencer products are exported to many countries where the same rules do not always apply. For this reason, there may be differences between what is described herein and the products delivered. Spencer is constantly working on improving all types and models of the products sold. We therefore rely on your understanding if we should reserve the right to make changes to the scope of delivery at any time in terms of form, equipment, set-up and technology.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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