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IT

**Instructions for use**  
**Davis - Traction system**

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**EN**

<b>1. MODELS</b>	<b>8</b>
<b>2. INTENDED USE</b>	<b>8</b>
2.1 INTENDED USE AND CLINICAL BENEFITS	8
2.2 TARGET PATIENTS	8
2.3 PATIENT SELECTION CRITERIA	8
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	8
2.5 USERS AND INSTALLERS	8
<b>3. REFERENCE STANDARDS</b>	<b>8</b>
<b>4. INTRODUCTION</b>	<b>8</b>
4.1 DEVICE LABELLING AND TRACEABILITY	8
4.2 SYMBOLS	8
4.3 WARRANTY AND SERVICE	8
<b>5. WARNINGS/DANGERS</b>	<b>8</b>
<b>6. SPECIFIC WARNINGS</b>	<b>9</b>
<b>7. RESIDUAL RISK</b>	<b>9</b>
<b>8. TECHNICAL DATA AND COMPONENTS</b>	<b>9</b>
<b>9. PROPER USE</b>	<b>10</b>
<b>10. CLEANING AND MAINTENANCE</b>	<b>10</b>
10.1 CLEANING	10
10.2 MAINTENANCE	11
10.3 LIFE SPAN	11
<b>11. TROUBLESHOOTING TABLE</b>	<b>11</b>
<b>12. ACCESSORIES</b>	<b>11</b>
<b>13. SPARE PARTS</b>	<b>11</b>
<b>14. DISPOSAL</b>	<b>11</b>

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## 1. MODELS

The following basic models may be subject to implementation or change without notice.

- DAVIS ADULT
- DAVIS PAEDIATRIC

## 2. INTENDED USE

### 2.1 INTENDED USE AND CLINICAL BENEFITS

Traction systems are devices that limit tissue damage caused by possible bone rubbing by spacing the two halves of a fracture.

### 2.2 TARGET PATIENTS

There are no particular indications related to the patient group.

The design of the product allows it to be applied to any subject as long as the chosen size of the limb allows for the intended purpose of the device to be achieved.

### 2.3 PATIENT SELECTION CRITERIA

The expected patients are those for whom traction is necessary.

### 2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

### 2.5 USERS AND INSTALLERS

The intended users are rescue workers with in-depth knowledge related to the immobilisation and handling of individuals with fractures or suspected fractures. These devices are not intended for lay people.

**Operators must be able to provide the necessary patient care.**

The product must be used only by personnel trained in the use of this product and not on other similar products.

The device does not require installation.

## 3. REFERENCE STANDARDS

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices


## 4. INTRODUCTION

### 4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

Regulation 2017/745/EU requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it>, or inform Customer Service (§ 4.4).

### 4.2 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the instructions for use
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal law restricts the sale of this device by or on the order of a licensed professional (US market only)
			Production identification Alphanumeric code that identifies the production units of the device, composed of:
(01)8057711230006 (11) 200626 (10) 1234567890		(01)805771123 000 6 (11)200626 (10)1234567890	company prefix progressive GS1 control number date of production (YYMMDD) lot number

### 4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase.**

Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, [service@spencer.it](mailto:service@spencer.it)

Warranty and service conditions are available on the website <http://support.spencer.it>

## 5. WARNINGS/DANGERS

### Product features

**Use of the product for any purpose other than that described in the User Manual is prohibited.**

- The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to +50°C.
- Storage temperature: -10°C to +60°C.

### General warnings for medical devices

- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities

for their respective actions.

- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

## 6. SPECIFIC WARNINGS

To use the traction system, you must also have read, understood and carefully follow all the instructions in the user manual.


- Perform traction simulations with dummies before putting the device into service.
- Improper use or use by untrained operators may result in injury or permanent disability. Always follow the procedures and protocols approved by the relevant Emergency Medical Service.
- To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- Qualified personnel must be present during use of the device.
- The traction system should be applied by at least two trained rescuers who have good manual dexterity and common sense and are trained in emergency rescue procedures.
- Before applying the traction system, appropriate primary medical evaluations should be performed to assess both the appropriateness of applying the device and any additional actions that should be taken in relation to the patient's injuries.
- The device must not be exposed, much less come into contact with thermal sources of combustion or flammable agents.
- Only use the device as described in this manual.
- Always check the conditions of all parts before use.

## 7. RESIDUAL RISK

No residual risks are identified, i.e. risks that could arise despite compliance with all the warnings in this user manual.

## 8. TECHNICAL DATA AND COMPONENTS

### DAVIS/DAVIS PAEDIATRIC

PARTS	DESCRIPTION AND MATERIALS
1	Gear wheel Made of aluminium, allows locking of the traction mechanism
2	Traction mechanism support Made of nylon, this is the structure on which the gear wheel, the pin around which the belt is wound, and the knob with which traction is performed are installed.
3	Traction belt Made of nylon and equipped with an eyelet to which a hook is attached; this is the element that transmits traction to the limb.
4	Locking ring Made of plastic material, these allow the removable portion of the frame to be locked in the desired position
5	Ankle brace Made of padded PVC and nylon, this is applied to the patient's ankle allowing traction of the limb
6	Support bands Made of nylon, these support the limb and, thanks to the closing straps with quick-release buckles, allow the device to be attached to the limb
7	Ischial support Made of polyurethane, this is the supporting part of the end of the femur
8	Fixed frame Made of aluminium, this is the fixed part of the telescopic frame
9	Mobile Made of aluminium, this is the telescopic part of the frame that allows the device to be extended to suit the height of the specific patient
10	Lock pin Made of aluminium, this is a pin that locks the support in the open position
11	Support Made of aluminium, when open it keeps the device slightly inclined and allows the limb to be moved away from the ground so that it does not touch the heel
12	Traction knob Made of plastic, when rotated this turns the pin on which the belt is wound, generating traction

Characteristics	David adult	Davis paediatric
Minimum length(1)	890 ± 10 mm	810 ± 10 mm
Maximum length(1)	1350 ± 10 mm	1170 ± 10 mm
Width	210 mm	185 mm
Base width	230 mm	225 mm
Width of leg support area	165 mm	135 mm
Traction belt length	From 0 to 500 mm	From 0 to 500 mm
Inclination	10° (variable depending on the extension)	10° (variable depending on the extension)
Materials	Steel,Al, Nylon	Steel,Al, Nylon
Weight	1,78 ± 0,1 kg	1,58 ± 0,1 kg
Bag dimensions	990x330 mm	920x300 mm
Bag weight	520 g	320 g

(1) The limb must be approximately 15 cm shorter than the indicated measurements  
Sizes subject to ± 10 mm tolerances.

IT

EN


DE

FR

ES

PT

## 9. PROPER USE

<p><b>1</b></p> <p>Place the device next to the leg to be immobilised, positioning the ischial support at the patient's hip. Loosen the locking rings by turning them counter-clockwise and adjust the length to the patient's limb. It is suggested that the length be adjusted so that the traction mechanism is approximately 20cm away from the heel of the healthy limb.</p> <p>Raise the device, open the support and, while holding it at an angle, check that it has locked properly.</p>		<p><b>2</b></p> <p>Once the required length has been reached, lock the locking rings by turning them clockwise and check that the telescopic frame has locked.</p> <p>Open all the bands, distributing them evenly along the frame, and prepare the patient by removing bulky clothing and shoes. Keep the limb aligned manually.</p>	
<p><b>3</b></p> <p>While a colleague gently lifts the fractured limb by a few centimetres, the second operator applies the ankle brace to the patient (snug but not tight) around the limb.</p> <p>Adjust the traction belt of the ankle brace as short as possible, thus providing the most space for traction available.</p>		<p><b>4</b></p> <p>The second operator grips the ankle brace and applies gentle traction while the colleague prepares by manually supporting the limb. The two rescuers lift the end of the limb until it is at an angle to the ground at which the device can be inserted.</p>	
<p><b>5</b></p> <p>Place the device under the limb, limiting limb movement as much as possible. Apply gentle traction to the ankle brace and insert the device under the limb so that the ischial padding is in contact with the patient's ischium.</p> <p>While maintaining light traction on the ankle brace, attach its ankle strap to the traction mechanism using the "S" hook.</p>		<p><b>6</b></p> <p>Fasten the ischial belt, passing the long end through the ring connected to the short end and then lock it in place.</p> <p>Tighten the belt, taking care not to obstruct circulation.</p>	
<p><b>7</b></p> <p>Once you have checked that the ankle brace is secure, turn the traction mechanism knob.</p> <p>The amount of traction should be carefully assessed according to the patient's condition.</p>		<p><b>8</b></p> <p>After checking the alignment of the limb, close all bands from the top to the ankle, taking care not to constrict circulation or compress the fracture site. Close the bands with the quick-release buckles and adjust the tension.</p> <p>For transport, it is best to load the patient on a spine board and cover with an insulated drape.</p>	

When the clinician considers it appropriate to remove the device from the patient, gently open the bands while supporting and holding the limb. Grasp the strap of the ankle brace in traction and pull the release tab of the traction mechanism. Gradually release the traction on the ankle brace. While continuing to support the limb, gently slide the device under the patient's leg until it is completely removed.

## 10. CLEANING AND MAINTENANCE

### 10.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**

Rinse thoroughly with lukewarm water, making sure you have removed all traces of detergent, which may deteriorate or compromise conditions and durability. **The use of high-pressure water should be avoided.**

Allow to dry thoroughly before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources. If blood is visible, oxidise it before washing with water.

The support bands and ankle brace can be removed and machine washed at 40°C.

**When disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

## 10.2 MAINTENANCE

The device does not require a routine maintenance program, but checks should be carried out before and after each use to verify:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Fulfilment of the requirements of the user manual in section Warnings and Specific Warnings.
- Fulfilment of the requirements of the manual in section on Proper use.

No periodic overhaul is foreseen for the device.

## 10.3 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase. Beyond this time, the device must be disposed of. Spencer Italia s.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have exceeded the maximum allowable life span.

## 11. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
The device does not maintain traction on the limb	Shortening of the Davis	Tighten the telescopic closures of the ankle brace and ischial belt
Device does not provide sufficient traction	Loose belts	Check tightening of straps on the ankle brace and ischial belt

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

## 12. ACCESSORIES

There are no accessories for these devices.

## 13. SPARE PARTS

RISR009	SET OF STRAPS AND BELTS FOR DAVIS ADULT
RISR010	SET OF STRAPS AND BELTS FOR DAVIS PAEDIATRIC

## 14. DISPOSAL

When devices are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

### Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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