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
### **Self Switch**

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
## User's Manual

### Self Switch

Exchange valve for cylinders containing medical gases

**CE**<sub>0123</sub> This appliance conforms with the Directive 93/42/CEE "Medical Devices"

Guarantee of Quality system for the production and the final control of the products certified by the notifying body TÜV SÜD Product Service GmbH

### INDEX

General information	page 2	Operating instructions	page 4
Warnings	page 2	Maintenance and cleaning	page 5
Description of product	page 3	Accessories and spare parts	page 6

First edition: 01/07/05  
Rev. 4: 23/02/15

## Thank you for choosing a Spencer product

### 1. GENERAL INFORMATION

#### 1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

#### 1.2 Conservation of the instruction manual

The instruction and maintenance manual must be kept together with the product, for the whole life of the device, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

#### 1.3 Symbols used

##### Symbol

##### Meaning

General or specific warning

See instructions for use

Serial Number

Product code

The product is compliant with the specifications of the Directive 93/42/CEE

#### 1.4 Servicing requests

For any information regarding the use, maintenance and installation, please contact the Spencer Customer Care Service on tel. 0039 0521 541111, fax 0039 0521 541222, e-mail [service@spencer.it](mailto:service@spencer.it) or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY. In order to facilitate the assistance service, please always indicate or communicate the serial number (SN) shown on the label applied on the box or on the device.

#### 1.5 Demolition

When the devices are no more suitable for being used, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations about demolition.

#### 1.6 Labelling

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the Manufacturer, the product, CE mark, serial number (SN). It must never be removed or covered.

### 2. WARNINGS

#### 2.1 General warnings

- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the Competent Authorities and/or Manufacturer if requested.
- Spencer Italia S.r.l. is always at your disposal to plan trainings on products.
- Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- If the instructions belong to another device and not the device received, inform the Manufacturer immediately and avoid use of the device.
- In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
- Regularly check the appliance, carry out the prescribed maintenance and respect the average life span, as indicated by the Manufacturer in this User's Manual.
- Before each use of device the perfect operating state of the device must be checked as specified in the Instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and or of the user are detected, the device must be immediately removed from service and the Manufacturer must be contacted.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- Use of the device in anyway other than described in this manual is forbidden.
- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- Handle with care.

- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids. Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments. When the device is being used, the assistance of qualified staff must be guaranteed.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store in a cool, dry, dark place and do not expose to direct sun.
- Store and transport device in its original packaging.
- The device not be exposed to or come into contact with any source of combustion or inflammable agents.
- Position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Attention: laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
- With reference to the D. Lgs. 24th February 1997, n. 46 emended by D. Lgs. 25/01/2010, n. 37 – Acknowledgement of Directive 93/42/CEE and 2007/47/CE, we remind both public and private operators that they are obliged to report any accident that involves any medical device to the Ministry of Health and to the Manufacture as specified and within time given by the European regulations.  
In addition, both public and private operators are obliged to inform the Manufacturer of any measures that should be adopted to make the steps necessary to guarantee the safety and the health of the patients and the users of any medical device.
- As a Distributor or End Users of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly required to have a basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the goods final destination Country (including laws and norms regarding technical specifications and / or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
- Promptly notify Spencer Italia S.r.l. regarding any revisions to be made by Manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant User Manual.
- Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the Manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- You are aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your non-compliance with the present "Regulatory provisions".

## 2.2 Specific warnings

- Establish a maintenance program and periodic testing, identifying a reference employee. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the Manufacturer in the user's manual.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the Competent Authorities and/or Manufacturer if requested.
- Use only accessories/spare parts that are original or approved by Spencer Italia S.r.l., in order to carry out any operation without causing any alteration or modification to the device, otherwise we assume no responsibility for the proper functioning or damage resulting from device to the patient or the operator and warranty and will be considered void according to the compliance to the Medical Device Directive 93/42/CEE.
- Avoid pulling the device on rough surfaces.
- The product may be used by trained personnel only.
- Do not switch the "self-switch" with gas pressure of 4.5 bar and over, in order to avoid damage to the device
- Keep the distribution pressure under control
- Clean the reducer filters regularly
- Check regularly the quick connecting joints
- Do not use countersunk head screw for fixation
- Do not open both connected oxygen bottles at the same time
- Protect the device from dust in order to avoid blocking of the device's movements.
- For lubrication of device, use only products compatible with oxygen.
- Do not use teflon tape or other materials and sealants liquid for input and output connections of the medical gas; their use could cause damage to the device.

## 2.3 Contraindications and side effects

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

## 3. DESCRIPTION OF PRODUCT

### 3.1 Intended use

Self-switch is an exchange for compressed medical gas, guaranteeing automatic exchange from an active bottle, tending to finish, and the spare bottle, correctly loaded with gas to be distributed to the patient. It has been designed in order to offer an ideal solution in case there is a need for continuous and safe medical gas distribution. Self-switch can be used with medical gasses like oxygen and medical air.

### 3.2 Main components

The device is assembled in a single body.

### 3.3 Models

This model could be modified, with reference to codes and/or descriptions without any previous notification.

OX07020 Self Switch medical gas bottle exchange system

### 3.4 Technical data

Characteristics	Technical data
Dimensions (mm)	17 x 50 x 30
Operation pressure (bar)	3,5 ± 0,5

### 3.5 Reference standards

Reference	Title of document
MDD 93/42/CEE	European Directive about Medical Devices
MDD 2007/47/CEE	Modifications to 90/385/CEE Directive about active implants, Directive 93/42/CEE about medical devices and Directive 98/8/CE about the introduction of biocides onto the market
Legislative Decree 24/02/1997, n. 46	Application of the 93/42/CEE Directive about Medical Devices
Legislative Decree 25/01/2010, n. 35	Modifications and additions to the 20/02/97 Decree n. 46
UNI EN ISO 14971	Application of risks managing to medical devices
UNI CEI EN 980	Graphic symbols used for medical devices labelling
UNI CEI EN 1041	Information supplied by the medical devices manufacturer
CEI EN 62366	Medical Devices - Application of the utilisation characteristics of engineering to medical devices
MEDDEV 2.4/1a-b	Guideline for the classification of medical devices
NB-MED 2.5.1 /Rec 5	Technical Documentation
MEDDEV 2.7.1	Clinical Data
MEDDEV 2.12/1	Medical Devices vigilance system
UNI EN 14155	Clinical evaluation of the medical devices for human beings - Part 2: Clinical evaluation plans

### 3.6 Environmental conditions

Functioning and storage temperature: from -5 to +40 °C

Relative humidity: from 30 to 80%

## 4. OPERATING INSTRUCTIONS

### 4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

Keep the original packaging for use in case of any further transport and for storage. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

### 4.2 Preparation

On receipt of the product:

- Remove the packaging and display the material so that all components are visible.
- Check that all the components/pieces on the accompanying list are present.

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage. In particular, check:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Absence of cuts, holes, tears on the structure, including the straps
- Integrity of components

### 4.3 Installation

**It is essential that the device is installed in an area inside the vehicle and immediately accessible by the operator.**

For fixing the Self Switch on a panel or another support use the 3 holes Ø 5,5 mm gained from the Self Switch structure.

Do not use countersunk head screw for fixation.

#### Connection for use

Connect the entrances of the Self Switch to the pressure reducers of the oxygen bottles, using an a-toxic tube compliant to the gas used and according to the reference norm UNI EN 5359.

Do not easy fast connection inlets because they do not grant a correct tightness using the device.

Use fittings with O-ring seals already integrated.

Make sure the medical gas distribution system's output is 3.5 ± 0.5 bar.

Verify the held of the system before use, in the following way:


- open the oxygen bottle in closed system mode
- measure the system's pressure value with a manometer
- close the oxygen bottle in closed system mode
- make sure that during a 2 h test period the loss of pressure must be less than 0,025%/h (for example using a pressure of 3,5 bar the loss of pressure must not be over 87 mbar/h).

During the test period a user must be connected to the exit of the Self Switch, in order to avoid both sockets automatically stabilize.

### 4.4 Functioning

The Self Switch is equipped with two entries (thread F1/4G) for connection to the bottles, one exit (thread F1/4G) for the use, red and green buttons for manual switch and two threads for the connection of pressure switches. Two visuals (red and green

button) allow to know which bottle is in use (button not pushed). The internal mechanism allows complete emptying of the used bottle. On the rear side of the Self Switch two threads 1/8G are to be found for connection of n°2 oxygen pressure switches. Pressure switches can be connected to visual or acoustic control systems in order to show which of the two bottles is in use or if the bottle pressure goes under 3 bar. Once connected the Self Switch, open one bottle at a time. In order to guarantee the outflow in case the mechanisms are altered, it is possible to use the manual switch, pushing one of the two buttons. During the test period a user must be connected to the exit of the Self Switch, in order to avoid both sockets automatically stabilize.

#### 4.5 Troubleshooting

PROBLEM	CAUSE	REMEDY
The automatic switch does not work	Insufficient lubrication; the adherence blocks the internal mechanism	Use the manual switch and send to Assistance Centre
	The pressure reducers of the bottles are broken	Use the manual switch and verify the correct functioning of the feeding system. If the problem continues to exist contact the Assistance Centre
Both the red and the green button are stabilized in central position	Both bottles have been opened contemporarily	Make sure the system has been connected correctly as indicated in the User's Manual.
	The system is not correctly installed at the user's end.	If the problem continues to exist contact the Assistance Centre

### 5. MAINTENANCE AND CLEANING

#### 5.1 Cleaning

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.

**The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures.**

#### 5.2 Maintenance

##### 5.2.1 Precautionary maintenance

The person who carries out the precautionary maintenance of the appliance (user in person, Manufacturer/supplier or a third party) has to guarantee the following basic requirements:

- Technical knowledge of the appliance and of the periodic maintenance procedures as described in these instructions.
- Specific qualifications and training in the maintenance operations of the appliance in question.
- The use of components/replacement parts/accessories that are either original or approved by the supplier, in such a way that each operation causes no alteration or modification to the appliance.
- Possession of the checklist of operations carried out on the appliance.
- Guarantee complete adherence to the instructions of the Directive 93/42/CEE which includes also the obligation towards the Manufacturer to maintain post sales records and traceability of the appliance if requested.

**During all checking, maintenance and cleaning procedures, the operator must wear adequate personal protection such as gloves, mask, glasses etc.**

Checks to be carried out before and after each use, and at least every month, are as follows:

- Check on the outflow pressure
- Cleaning of the pressure reducer's filters: the gas, if not filtered correctly, may damage the internal mechanisms of the system.

Checks to be carried out before and after each use, and at least every 3 months, are as follows:

The inspection frequency is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage. Please note that you must do the cleaning as described in paragraph 5.1

and verify functionality before and after each use. Spencer Italia S.r.l. declines any responsibility for the proper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance warranty and will void the compliance to the Medical Device Directive 93/42/CEE.

The person responsible for routine maintenance can identify damaged/worn parts, but the replacement or restoration of the them can only be done by the manufacturer or or by an authorized service centre.

Use only accessories/original spare parts approved by Spencer Italia S.r.l., otherwise we will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres. Warranty will be considered void in compliance with the Medical Device Directive 93/42/EEC.

### **5.2.2 Periodic maintenance**

The device must be serviced by the manufacturer or by an authorised centre **every year**.

If the correct revision is not carried out, the CE branding will no longer be considered valid as it will no longer be compliant with the 93/42/CE Directive for Medical Devices and therefore there is the possibility that it is no longer compliant with the safety standards declared by the Manufacturer at time of purchase. Spencer Italia S.r.l. will take no responsibility the incorrect functioning or any damage caused by a device that has not undergone regular revision.

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

### **5.2.3 Special servicing**

**Only the Manufacturer or centres with written authorisation are authorised to complete any special servicing operations.**

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device, if used as indicated in the following instruction manual, has an average life span of 5 years. The life span can be expanded only following a general revision of the product that must be carried out by the Manufacturer or by a centre authorised by the Manufacturer.

Spencer Italia S.r.l. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres, making void the guarantee and the conformity to the Medical Devices Directive 93/42/CEE.

## **6. ACCESSORIES AND SPARE PARTS**

### **6.1 Accessories**

OX10104B        ATOXIC TUBE FOR MEDICAL OXYGEN EN 5359  
OX10109B        ATOXIC TUBE FOR MEDICAL MEDICAL AIR EN 5359

The tube for medical gas may be supplied only on specific request with the ends already terminated with threaded connection. Before ordering, please make a specific request for the required length and connections.

### **6.2 Spare parts**

There aren't spare parts for this product.