TOTAL**TOTAL KON**

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IT

**Use and Maintenance Manual
TOTAL/TOTAL KON**

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EN

1. MODELS	11
2. INTENDED USE	11
2.1 INTENDED USE AND CLINICAL BENEFITS	11
2.2 TARGET PATIENTS	11
2.3 PATIENT SELECTION CRITERIA	11
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	11
2.5 USERS AND INSTALLERS	11
2.5.1 USER TRAINING	11
2.5.2 INSTALLER TRAINING	11
3. REFERENCE STANDARDS	11
4. INTRODUCTION	11
4.1 USING THE MANUAL	11
4.2 DEVICE LABELLING AND TRACEABILITY	11
4.3 SYMBOLS	12
4.4 WARRANTY AND SERVICE	12
5. WARNINGS/DANGERS	12
6. SPECIFIC WARNINGS	13
7. RESIDUAL RISK	13
8. TECHNICAL DATA AND COMPONENTS	13
9. COMMISSIONING	13
10. OPERATING CHARACTERISTICS	14
11. PROPER USE	14
11.1 LOADING THE PATIENT ONTO THE TOTAL STRETCHER	14
11.2 TRANSPORTING THE STRETCHER	15
11.3 STORING THE DEVICE	15

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- TOTAL
- TOTAL KON

2. INTENDED USE

2.1 INTENDED USE AND CLINICAL BENEFITS

Recovery stretchers are devices to be used when the morphological characteristics of the intervention area require particularly compact equipment and adherence to the patient.

Total can also be lifted with the aid of hoisting systems fixed to the ground, using the dedicated harness system.

It is not foreseen that the patient be able to intervene on the devices.

2.2 TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity and within the limits of the size of the device .

2.3 PATIENT SELECTION CRITERIA

The patients foreseen are those with injuries that prevent them from walking in a given rescue situation, or who are in a state of unconsciousness.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

2.5 USERS AND INSTALLERS

The intended users are rescue workers, typically technical, carrying out operations related to the use of harnessing systems.

- Personnel trained for use of the device must also have training in managing lifting and handling suspended loads with people.
 - Personnel who carry out interventions in situations classified as high risk or which are purely technical must be suitably trained and experienced in rescue.
- These devices are not intended for lay people.

Recovery stretchers are devices intended for professional use only. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs, should it be necessary to lift and/or support the device and the patient. Operators' ability must be assessed before the definition of roles in use of the stretcher.

Operators must be able to provide the necessary patient care.

■ 2.5.1 USER TRAINING

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note: Spencer Italia S.r.l. is always available for training courses.

■ 2.5.2 INSTALLER TRAINING

Installation is not required.

3. REFERENCE STANDARDS

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

4. INTRODUCTION

4.1 USING THE MANUAL

The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer products User Manuals can be downloaded from the site <http://support.spencer.it> or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). **This must never be removed or covered.**

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

If the assigned Lot/SN cannot be traced, the device must be reconditioned, provided only under the responsibility of the manufacturer.

EU Regulation 2017/754 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at <http://service.spencer.it> or inform Customer Service (see § 4.4).

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4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the user manual.
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner (only for USA Market)

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at <http://support.spencer.it>.

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal (<http://service.spencer.it>) available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

5. WARNINGS/DANGERS

 Warnings, dangers, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website www.spencer.it on the specific product page.

Product features

 Use of the product for any purpose other than that described in the User Manual is prohibited.

- Before each use, always check the conditions of the product, as specified in the User Manual. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
- If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations. Non-compliant devices must be taken out of service.
- The product must not be tampered with or modified without the manufacturer's authorisation (modification, tweaking, additions, repair, use of non-approved accessories), as they may constitute imminent danger of injury to persons and material damage. Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE marking and the product warranty shall be null and void.
- When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
- Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
- Avoid contact with sharp or abrasive objects.
- Device installation, if necessary, must be carried out by qualified personnel trained and enabled by Spencer Italia S.r.l. The times and methods of such shall be agreed upon between the customer and our Sales Offices.
- Operating temperature: from -5°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the structure of the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -10°C to +60°C.

Regulatory requirements

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the User Manual.
- **Participate in safety checks on products** placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

 The user must carefully read the following in addition to the general warnings.

- The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point.
- Qualified personnel and at least two operators must be present during use of the device.
-  Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Follow the internal procedures and protocols approved by your organisation.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's

warranty and release the manufacturer from all liability.

- Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.
- Do not use drying machines to dry the device.
- In case of exposed and/or injured skin, cover the surfaces in contact with the patient with a surgical sheet that respects bio-compatibility regulations to protect the patient's health.

6. SPECIFIC WARNINGS

⚠ To use the Total device, you must also have read, understood and carefully follow all the instructions in the user manual.

- Always comply with the maximum capacity, if any, indicated in the User Manual. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. In addition, the operator should assess whether the overall size of the patient reduces the functionality of the product.
- Make sure that operators are in good physical condition before lifting, as listed in the User Manual.
- **The maximum weight, which weighs on each operator, must comply with local health and safety requirements.**
- **The warranty seals, if present on the product, must not be removed;** otherwise, the Manufacturer shall no longer recognise the product warranty and shall decline all responsibility for incorrect operation or any damage caused by the product itself.
- Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with routine maintenance of the device must ensure the basic requirements envisaged by the manufacturer within these operating instructions.
- All maintenance activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the manufacturer when requested.
- Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.
- Never leave the patient on the device unsupervised, as they could get injured.
- Avoid contact with sharp objects.
- Follow approved Emergency Medical Service procedures for patient immobilization and transportation.
- Follow approved Emergency Medical Service procedures for patient positioning and transportation.
- Before lifting, make sure that operators have a secure grip on the device.
- The device is a patient transport stretcher and cannot be used as a stationing device.

- **⚠** The device is not intended for spinal immobilisation and if such injuries are suspected or evident, suitable devices should be used according to the guidelines of the relevant emergency medical service.
- Practice with an empty stretcher with a simulated patient load to make sure you are familiar with the manoeuvres.
- At least two operators for Total Kon and four operators for Total, in suitable physical conditions, are required for use of the device. They must therefore be endowed with strength, balance, coordination, common sense and must be trained in the correct operation of the Spencer stretcher device.
- For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of more than two operators is recommended in addition to the minimum foreseen operators.

- **⚠** Before each use, always check the conditions of the device and its components, as specified in the user manual. In case of faults or damage that may compromise the functionality and safety of the device, and therefore the patient and the operator, remove the device from service or replace the components that are not intact.
- Make sure that the straps are properly installed.
- Always immobilise the patient using the equipment supplied with the device. A failure to do so can cause serious damage to the patient.
- Do not move the stretcher if the weight is not properly distributed.
- Use the stretcher only as described in this user manual.
- Do not alter or modify the stretcher to adapt it to unforeseen conditions of use: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise proper functioning of the device. If you cannot clear the route, choose an alternative route.

- **⚠** During hoisting, the stretcher must be kept horizontal with respect to the ground. Any abnormal inclination can cause serious damage to the patient, the device and the operator.
- If it is necessary to use ropes, winches, ladders, lifting straps or other special equipment to move the stretcher or in the presence of a rescue that is classified as high risk or of a purely technical nature, these interventions must be carried out solely by personnel adequately trained and with proven experience in rescue.
- The Total stretcher can only be hoisted from a fixed position using the spring catch supplied with the device. Total Kon cannot be hoisted.
- Dragging the stretcher on any type of surface leads to premature deterioration of the stretcher, which reduces its useful life and initial safety conditions.
- **⚠** To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions.
- Recovery stretchers are not approved for use with aircraft.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in this user manual, have been identified.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

	DESCRIPTION	MATERIAL
1	Support bed	Polyethylene
2	Patient containment belts	Polyester
3	Rope	Polyester
4	Extension	Polyester
5	Removable handles	Polyester
6	Harnesses	Nylon
7	Grommets	Brass
8	Handles	Polyester
9	Transport bag	PVC

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device

IT

EN

DE

FR

ES

PT

EL

DA

NL

- Product cleanliness
- There are no cuts, holes, lacerations or abrasions on the entire structure of the device and its parts, including detachable parts.
- Correct fixing and holding of straps
- Check that the grommets for anchoring the handles, where provided, are intact and correctly attached to the device.
- Conditions of wear of the device and its standard parts.

See paragraph 11 for how to carry out the above-mentioned checks.

Do not modify the device or its parts for any reason as this could cause damage to the patient and/or rescuers.

⚠ Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself. For subsequent use, perform the operations specified in paragraph 12.

If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

10. OPERATING CHARACTERISTICS

See paragraph 11 - Proper use for operating characteristics.

11. PROPER USE

Primary medical evaluations must be carried out before intervening on the patient.

11.1 LOADING THE PATIENT ONTO THE TOTAL STRETCHER

Before moving the patient, a medical evaluation of the patient's condition must be performed to stabilize the patient's clinical condition, verify possible surrounding hazardous situations, assess how to move the patient and whether it may be necessary to use alternative devices. Once these priorities have been met, you may proceed to the next steps in using the device.

Unrolling and stretching out (Total and Total Kon)

- Remove the device from the case.
- Unlock all the straps and unroll the stretcher completely, while holding one end of the support bed (Fig. A - Appendix) firmly on the ground with one foot.
- Then roll up the upper half of the stretcher in the opposite direction to that in which it was taken out of the bag in order to eliminate the "memory" effect of the material. Carry out the same operation on the opposite side in order to make the stretcher perfectly flat. (Fig. B - Appendix)
- Place the stretcher on the ground according to the indications on the "HEAD SIDE" and "EXTERNAL SIDE" stickers.

Patient positioning - LOG ROLL method

- Place the device next to the patient. Make sure that the side of the device on which the patient's head is to be placed is adjacent to the head itself. Position the transverse straps under the stretcher.
- Roll the patient to one side and slide the device as far as possible under the patient. Gently lay the patient on the device (Fig.C - Appendix).

Patient positioning – PULL method (head)

- Fold the foot end of the stretcher to form a kind of ramp and place the belts under or to the side of the stretcher.
- Using forearms to stabilise the head, grasp the patient in the underarm area and then pull the patient onto the stretcher until the stretcher is in the proper position for anchoring (Fig.D - Appendix).

Patient positioning – PULL method (feet)

- Fold the head end of the stretcher to form a kind of ramp and place the belts under or to the side of the stretcher.
- Grasp the patient's legs and pull him or her onto the stretcher until he or she is in the proper position for anchoring.

Patient anchorage

- Make sure that you have correctly centred the patient in relation to the stretcher.
- Raise the sides of the stretcher and fasten the straps, wrapping the stretcher around the patient. Do not cross the straps.
- Grasp the belts on the foot end, pass them through the last grommet in the straight portion of the stretcher and fasten them with the buckle (Fig. E - Appendix). Adjust so that it is tensioned enough to ensure the patient's stability.
- Slightly roll up the head end of the stretcher to prevent the patient from sliding and adjust the belts so that this position is maintained.

Total Kon patient anchorage

- With Total Kon, the upper strap must be passed below the patient's armpits.
- To attach Total Kon to the patient, the groin belts must also be carefully positioned and anchored (Fig. F - Appendix).
- These are located at the bottom of the device. After passing the belts through the patient's groin, close the quick release buckles, making sure that the belts are securely fastened and able to support the patient's weight.

Rope positioning (for vertical lifting)

- Insert each end of the white rope into the grommets at the head end of the stretcher.
- Then pass it through the free grommets to the foot end, alternately inside and outside the stretcher. Adjust the length of the rope so that the two ends protruding from the stretcher are the same length (Fig. G - Appendix).
- Knot the ends of the rope so that it properly tightens the device, then pass the portions of rope that are still free through the last pair of foot handles and tie an additional knot (Fig. H - Appendix).

Extra handles

- The device is supplied with 4 additional handles whose use is recommended for horizontal lifting in especially heavy load conditions where more operators than those normally used are foreseen.
- To install these handles, align and join their edges, pass them through the free grommets and tie a knot so that they cannot come out of the grommets (Fig. I - Appendix).

Note: *It will not be possible to attach these handles if a rope has been installed, as all the grommets will be occupied.*

Extra drag handle

- The drag handle can be useful if 3 operators are needed for dragging.
- The drag handle should be attached to the belt closest to the patient's head by means of the spring catches (Fig. J - Appendix).

Applying the hoist straps

- Total can be hoisted using the two orange straps supplied with the device.
- The straps must be positioned underneath the stretcher, then passed through the dedicated slots (Fig. K - Appendix)
- After releasing the straps from the inside of the stretcher, adjust the length of all the ends so that the stretcher can be lifted horizontally, then attach the spring catch to the ends of the straps.
- Hoisting is only allowed from a fixed position.

IT
EN

DE

FR

ES

PT

EL

DA

NL

11.2 TRANSPORTING THE STRETCHER

TRANSPORT BY HAND

- Transport can only be performed if the patient is correctly positioned as described in 11.1 and all requirements specified in this user manual are met.
- Transporting the stretcher requires a minimum of 2 operators.
- A greater number of operators is necessary in the case of transport for long and/or difficult journeys. Moreover, in the case of particularly heavy patients, evaluate the weight that each rescuer can support in compliance with the legal requirements in terms of safety in the workplace and according to the physical conditions of the operator.
- Recovery stretchers are equipped with various gripping points along the entire perimeter of the device itself. This lets operators position themselves at the point considered most suitable by the coordinator of the handling operations.

STRETCHER SLIDING (Total Only)

- Some rescue situations may require sliding the stretcher on the ground. In these cases, pay attention to the presence of obstacles that may create danger to the patient, the operators and that may damage the device itself.

TRANSPORT BY OTHER MEANS

- If it is not possible to transport the stretcher by hand, but ropes, winches, ladders, lifting straps or other suitable means are used, this is considered high risk rescue and must be carried out only by personnel adequately trained in this area and with the necessary experience on the field.

11.3 STORING THE DEVICE

After carrying out cleaning and maintenance (par. 12), proceed as follows:

- Lay out the Total or Total Kon stretcher and place the straps with the metal grommets under the end of the stretcher.
- Starting at the head end, roll up the stretcher as tightly as possible.
- Continue rolling up the device, using your knees to prevent it from rolling out.
- Fasten the straps to the metal rings and place the stretcher in the backpack.

12. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the User Manual, identifying a reference person who meets the basic requirements set forth in the User Manual.
- **The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.**
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l. declines any responsibility for any direct or indirect damage which is a consequence of improper use of spare parts and/or any repair work carried out by unauthorised parties.
- Reconditioning, a process performed on the device to restore the technical and functional safety of the device used, for example re-registration, must be performed by the Manufacturer.
- Where foreseen, use only original or Spencer Italia S.r.l. approved components/spare parts and/or accessories in order to carry out all operations without causing alterations or modifications to the product.
- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the User Manual in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- If required, the product and all its components must be washed and left to dry completely before storage.

12.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.
The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Any metal parts exposed to external agents undergo surface treatments and/or coating in order to obtain better resistance. Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it creates the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

12.2 ROUTINE MAINTENANCE

You must establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged within this user manual.

All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.

Checks must be made before and after use and at least every 3 months to verify:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Holding and conditions of all straps and buckles
- Conditions of grommets
- Conditions of the spring catch and all its parts
- Fulfilment of the requirements of the user manual in section 5 Warnings and 6 Specific Warnings.
- Fulfilment of the requirements of the manual in section 11 Proper use.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

12.3 PERIODIC OVERHAUL

No periodic overhaul is foreseen for the device.

12.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

IT

EN

DE

FR

ES

PT

EL

DA

NL

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l..
The end-user can replace only the spare parts indicated in § 15.

12.5 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have exceeded the maximum allowable life span.

13. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
Damage on the support bed	Normal wear or improper use	Immediately remove the device from service and contact Customer Service
Damage to grommets	Improper use	Check the condition of the grommets after each use; if permanent deformation is evident, take the device out of service immediately and contact Spencer Customer Service
Damaged straps	Wear	Immediately remove the device from service and replace the damaged straps or contact Customer Service
The spring catch locking ring does not rotate	Possible dirt inside	Blow with compressed air until movement is released. Lubricate if necessary.
The spring catch does not close.	Possible oxidation or dirt	Try localised treatment with an anti-corrosive lubricant. If the problem persists, replace the spring catch.

14. ACCESSORIES AND SPARE PARTS

There are no accessories for this device. There are no spare parts for this device.

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

EL

DA

NL

Notice

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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