### Rigid splints with flexible internal structure Férulas rígidas con alma flexible

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Manually immobilise the injured limb in either the analgesic position or maintaining the patient's position. Prepare a Blue Splint/Blue Splint Pro which is the correct size for the limb (it must cover both the articulation above and below the fractured bone). Bend and adapt the metal bar inside to follow the morphology of the injured limb.

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With care apply the splint keeping the metal bar under the limb. Close the self-adherent strips straps or the buckles starting from the distal area and making sure that it is not too tight in the injury area.

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Follow the same procedure for the upper and the lower limbs.

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More than one Blue Splint/Blue Splint Pro can be applied on the same limb so as to increase the immobilized area.

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The Blue Splints/Blue Splints Pro are extremely light and can be overlapped with each other to create an infinity of configurations and solutions. They can be used for any size of adult or child.

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The arm and forearm splints can also be used as an armrest during intravenous infusions.

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Cleaning can be carried out using alcohol or other neutral disinfectants. Before machine washing, it is necessary to remove the metal core. The straps of the Blue Splint Pro model must be closed before machine washing.

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