**1** Remove the Davis traction splint from its case and place it next to the leg that has to be immobilized. Place the ischial padding part and align it with the patient's hip, loosen the extension rings and extend the device until you reach the patient's leg length. The best way to reach a perfect and comfortable length is to pull the ratchet of the splint 20cm beyond the heel of the other leg. Incline the device by using the stand and check that the stand is safely locked into position.

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**3** Whilst one operator (operator 1) carefully lifts the injured leg, the other (operator 2) applies the ankle cuff (not too tight but firm) around the patient's ankle. Adjust the traction strap so that it is as short as possible, in order to have bigger traction space.

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**5** Place the "Davis" underneath the leg, trying to avoid too many movements. Slightly pull the ankle cuff's traction strap and slide the Davis towards the patient's ischium, until the ischial pad touches the ischium. Attach the ankle cuff's strap to the hook of the ratchet, then pull gently.

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**7** After checking the ankle cuff fixation, turn the ratchet knob until the patient feels less pain (conscious patient), or when the rotation of the knob becomes hard and difficult (unconscious patient).

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## "DAVIS" TRACTION SPLINT

"Davis" Traction Splint is the best method for immobilizing lower limb fractures. The mechanism allows to reduce pain and limit risks of damage to the area surrounding the fracture, by keeping traction constant and allowing hemodynamic re-directioning of the limb. For device application at least two experienced rescuers must operate together, coordinating their movements and action.

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**2** Once you have tightened the extension rings, always check if they are effectively blocked. Open all the "Davis" straps so to be ready for inserting the patient's limb. Prepare the patient by removing clothes and shoes; manually keep the leg aligned with the splint.

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**4** Operator 2 pulls the ankle cuff's strap whilst operator 1 prepares to lift the leg. The two rescuers will then lift the foot and find enough space to apply the Davis splint.

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
**6** Fix the ischial strap by pulling the longer part through the ring connected to the shortest part, then close with self-adherent strips. Tighten the belt enough not to block the patient's blood circulation.

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**8** After checking the alignment of the leg, close all the Davis straps, starting from the top of the leg to the bottom (ankle). Pay attention to any blood circulation restrictions to the fracture area and do not tighten too much. Fasten the leg restraining straps by pulling the elastic part through the ring and close with self-adherent strips. Before proceeding to transport the patient, we suggest the use of a spine board and isothermic sheet.

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User's Manual

“DAVIS” traction splint user manual

This appliance conforms with the directive 93/42/CEE “Medical Devices” and following modifications.  
Guarantee of Quality system for the production and the final control of the products certified by the notifying body TÜV SÜD Product Service GmbH.

1 GENERAL INFORMATION

1.1 USE AND CONTENTS

This manual is intended to provide the health care operator with the all the necessary information for its safe and appropriate use as well as adequate maintenance of the device.  
**Note:** The manual is an integral part of the device. It must be kept for the duration of the device and must accompany the device in case of change of ownership or destination. If the operating instructions received relate to products not received, you must immediately contact the manufacturer before use.

1.2 SYMBOLS USED

- General or specific warning
- See instructions for use
- Lot number
- Product code
- The product is compliant with the specifications of the Directive 93/42/CEE

1.3 REQUEST FOR SERVICE

For any information regarding the correct interpretation of the instruction manual, the use, maintenance, installation and restore of the product, please contact the Spencer Customer Care Service ph. +39 0521 541111, fax +39 0521 541222, e-mail export@spencer.it for general informations about Spencer products service@spencer.it for technical assistance or write to Spencer Italia S.r.l. - Via Provinciale, 12 - 43038 Sala Baganza (Parma) - ITALY. In order to facilitate the assistance service, please always indicate the lot number (LOT) or serial number (SN) shown on the label applied on the box or on the device. **Conditions for warranty and assistance can be viewed on <http://support.spencer.it>.**

1.4 DEMOLITION

When the devices are no more suitable for being used, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations about demolition.

1.5 LABELLING

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the Manufacturer, the product, CE mark, lot number (LOT), lot number (LOT), lot number (LOT) or serial number (SN) shown on the label applied on the box or on the device. **Conditions for warranty and assistance can be viewed on <http://support.spencer.it>.**

2 WARNING

2.1 GENERAL WARNING

- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the Competent Authorities and/or Manufacturer if requested.
- Spencer Italia S.r.l. is always available for conducting training courses.
- Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- Onomatopoeic sounds may be emitted when the device is used as they may cause injury to the patient or themselves.
- Perform the required maintenance and to respect the life span of the device, as indicated by the Manufacturer in the User's Manual.
- Before each use of the device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the Manufacturer must be contacted.
- Use of the device in anyway other than described in this manual is forbidden.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the device itself, moreover CE certification and product warranty will be considered void.
- Handle with care.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments.
- When the device is being used, the assistance of qualified staff must be guaranteed.
- Never leave an unattended patient. The presence of at least one operator is essential at all times when the medical device is in use.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store in a cool, dry, dark place and do not expose to direct sun.
- Store and transport device in its original packaging.

2.2 SPECIFIC WARNING

- Establish a maintenance program and periodic testing, identifying a reference employee. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the Manufacturer in the user's manual.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the Competent Authorities and/or Manufacturer if requested.
- Use only accessories/spare parts that are original or approved by Spencer Italia S.r.l., in order to carry out any operation without causing any alteration or modification to the device, otherwise we assume the responsibility for the proper functioning or damage resulting from device to the patient or the operator and warranty will be considered void according to the compliance to the Medical Device Directive 93/42/CEE.
- Davis allows leg immobilization and traction.
- Device use must be carried out in compliance with the local Health Authority regulations.

CONTRA-INDICATIONS AND UNDESIRABLE EFFECTS

2.3 This device, if used by trained personnel, does not include contra-indications or undesirable effects.


3 DESCRIPTION OF THE DEVICE

3.1 INTENDED USE

Davis is a traction/immobilization system especially designed to immobilize femur fractures, relieve pain and prevent damage to tendons, muscle tissues and the vascular system. Its exclusive traction mechanism relaxes the fractured femur into the correct position. It restores the pressure of the tissues back to the traumatized leg, limiting haemorrhagic losses; it avoids damage to the nervous, vascular, muscular and bone tissues. The telescopic system, micrometrically adjustable, combined with the extendable standing arch and the unbreakable nylon transport bag, make the Davis a compact device, easy to use and to store. It is not intended for the patient to act on the device.

3.2 MAIN COMPONENTS

- The device is composed of the following parts:
1. traction gear mechanism and "S" shaped coupler
  2. control system
  3. heel stand and ring for blocking system release
  4. traction gear knob
  5. traction belt with hook
  6. telescopic poles extension locknuts
  7. four padded leg support straps
  8. strap and ischial padding
  9. frame telescopic part
  10. frame fixed part
  11. safety lock
  12. extensible strap for ankle
- Davis is provided with a carrying case.

3.3 MODELS

These basic models could be modified, with reference to codes and/or descriptions without any previous notification.

DESCRIPTION	CODE
DAVIS traction splint adult	SR01010A
DAVIS traction splint paediatric	SR02020D
DAVIS HP traction splint adult	SR01040B
DAVIS HP traction splint paediatric	SR02030B

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DESCRIPTION	ADULT	PAEDIATRIC
Width	205 mm	170 mm
Maximum Length	1350 ± 10 mm	1170 ± 10 mm
Height closed	230 mm	225 mm
Weight	1,78 ± 0,1 kg	1,58 ± 0,1 kg

3.5 ENVIRONMENTAL CONDITIONS

Acceptable conditions for operating, storage and transport: operating temperature from -20 to +50 °C Humidity: 10 - 90%  
storage and transport temperature from -20 to +60 °C

4 OPERATING INSTRUCTIONS

4.1 TRANSPORTATION AND STORAGE

Always keep the integrity of the stretcher in all its parts before use; before effecting transport ensure the device is packed adequately and that risk of bumping or breaking is minimised. Keep the original packing for possible future transportation. Any damages to the stretcher caused during transport and moving are not covered by the warranty. Repairing or replacing of damaged parts are at the Customer's charge. The device must be stocked in a dry and humidity-free environment.

4.2 PREPARATION

On receipt of the product remove packaging and check all pieces indicated on the included list. If the boxes or packaging are damaged, contact the freight forwarder and then Spencer Customer Service.

4.3 FUNCTIONING

See back

4.4 TROUBLESHOOTING CHART

PROBLEM	CAUSE	ACTION
The device does not hold leg traction	Shortening of the device	Lock the telescopic poles action locknuts
The device does not allow proper traction	Straps are loose	Check the ankle and ischial straps fastening
A peripheral ischemia of the immobilized leg is found	Ankle and ischial straps fastening	Release straps along the leg starting from ankle

5 MAINTENANCE

5.1 CLEANING

**During all checking, maintenance and cleaning procedures, the operator must wear adequate personal protection such as gloves, mask, glasses etc.**  
If cleaning operations are not regularly carried out there is a risk of cross infection due to the presence of residues on the device. For a correct device storage, clean it as follows:  
• Clean by means of a clean cloth and bactericidal or germicidal disinfectant.  
• Rinse with hot water and let it dry. Keep clean outward surface by means of a clean cloth.

5.2 MAINTENANCE

- The person who carries out the precautionary maintenance of the appliance has to guarantee the following basic requirements:
- Technical knowledge of the appliance, of the periodic maintenance procedures as described in these instructions;
- Use of technical personnel in possession of specific qualifications and training in the maintenance operations of the appliance in question;
- Use of components/replacement parts/accessories that are either original or approved by the supplier (if scheduled), in such a way as to carry out every operation without causing alterations or modifications to the appliance;
- Possession of the checklist of operations carried out on the appliance;
- Guarantee complete adherence to the instructions of the 93/42/CEE Directive also regarding the obligation towards the manufacturer to allow the aforementioned a post sales care and traceability of the appliance when requested.

5.2.1 PRECAUTIONARY AND PERIODIC MAINTENANCE

No programmed revision or maintenance by the manufacturer is needed. The product needs periodic maintenance and checking to protect its efficiency and to extend its life. It is important to check the Davis traction system at least once a month and more times if the product is often used or if it is exposed to intense mechanical stress. During the periodic check it is necessary to test the functionality of the traction system, in order to find out potential troubles which could damage the Davis traction system. The person responsible for every day maintenance can only substitute the spare parts indicated on the paragraph "Spare Parts". All other substitutions or repairs can be carried out only by the Manufacturer or by a centre authorised by the Manufacturer.

5.2.2 SPECIAL SERVICING

Only the Manufacturer or centers with written authorization are authorized to complete any special servicing operations. The device, if used as described above, has an average life span of 5 years.

6 ACCESSORIES AND SPARE PARTS

6.1 SPARE PARTS

DESCRIPTION	adult CODE	paediatric CODE	DESCRIPTION	adult CODE	paediatric CODE
Carrying case	SR01011A	SR02021A	Heel stand	SR01015A	SR02023A
Accessories kit	SR01012A	SR02022A	Control system (complete)	SR01013A	SR01013A
Ischial strap and pad	SR01014A	SR02024A			

Warning  
The information contained in this document can be modified without warning and is not to be intended as a commitment on the part of Spencer Italia S.r.l. Spencer products are exported to many countries in which the same identical regulations do not exist. It is for this reason that there might be any difference between the description and the actual product delivered. Spencer works continuously to reach the perfection of all the models sold. We therefore hope to have your understanding if we reserve the right to modify the shape, equipping, layout or technical aspects that are herein described, at any given time.

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