



LIENZ A.L.S. Alluminium case

Code CB06002A

Index:

1. GENERAL INFORMATION
2. DEVICE DESCRIPTION
3. SAFETY
4. OPERATIONS
5. MAINTENANCE

1. GENERAL INFORMATION

This User's Manual must be stored near the product, protected from anything which could compromise its integrity and legibility.

1.1. SYMBOLS USED

SYMBOL	MEANING	NOTES
	DANGER	It indicates a danger with mortal risk for the user.
	WARNING	It indicates a warning or a note about key functions or about useful information. Pay maximum attention to the blocks of text marked out by this symbol.

- **In case of doubts about the interpretation of the instructions, please contact the Manufacturer in order to obtain the necessary explanations.**
 - **Operators must read this manual carefully before undertaking any action**
 - **The device must be used by specialized medical personnel only**

2. DEVICE DESCRIPTION

LIENZ is a sturdy aluminium case for emergency care equipment. It's strong and easy to transport. Designed for Ambulance Services, Emergency Medical Teams, First Aid Teams.

2.1 DIMENSIONS AND TECHNICAL DATA

LIENZ characteristics:

- base rails for extra stability;
- 3 hooks for carrying straps or security straps;
- reinforced corners;
- large, comfortable handle, flexible core;
- anti-rip hinges to withstand extreme stress fracture panel;
- 2 locks with key
- atraumatic hygienic plastic movable divisions;
- Size: 22x52x38 cm;
- weight: 6.2 Kg.

3 – SAFETY

3.1 GENERAL WARNINGS



- **The operator must read with high attention the informations reported in the present Manual, with particular reference to the adequate precautions for safety and paying attention to the using methodologies.**
- **Make sure that all precautions were taken to avoid hazards deriving from with blood or body fluids.**

4. OPERATIONS



Damages of the device caused during the transport and the handling are **not** covered by guarantee. Reparings and replacements of damaged parts are to customer's charge.

If there are any damaged relating to the cartons or to any other packaging material, please contact the transporter and following The Customer Assistance Service Spencer tel. 0039 0521 541111, or write to Spencer Italia Srl.

5. MAINTENANCE

5.1 PRESERVATION

The product requires periodical maintenance and control in order to preserve the efficiency and to extend the duration.

The check must entail verifying the integrity of the LIENZ, researching of anomalies which can compromise the correct working

5.2 REQUEST FOR SERVICE

For any information concerning use, maintenance, etc., contact the manufacturer.

5.3 DEMOLITION

Follow laws in force.



Thank you for choosing

www.spencer.it

© Copyright Spencer Italia S.r.l.

All rights reserved

Strada Cavi, 7

43044 Collecchio (PR) Italy

Tel. +39.0521.541111 – Fax +39.0521.541222 – E-mail info@spencer.it