

## B-Life

### Reusable manual resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands  
 ADULT and REUSABLE marked on the device to minimize risks of improper applications  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply.  
 Low storage space required  
 Latex Free

#### Technical data

Materials	Silicone, PVC, PC
Bag volume	1600 ml
Stroke volume	770 ml with 1 hand 900 ml with 2 hands
Max. re expansion rhythm	45 bpm

#### Optional accessories

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
BG03091A	Reservoir bag
EV50011C	EVX 11 Overpressure non-rebreathing valve
RM20405A	B-LIFE Face mask size. 5

Class IIa med. dev. according to dir. 93/42/CEE

B-LIFE - BODY ADULT  
 RESUSCITATOR WITH R-VALVE  
**BG01003A**



Registration number 216104

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)