

B-Life

Reusable manual resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Caratteristiche specifiche

The material of construction and the surface finish increase grip to the operator's hands
 CHILDE and REUSABLE marked on the device to minimize risks of improper applications
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply.
 Low storage space required

Technical data

Materials	Silicone, PVC
Bag volume	500 ml
Stroke volume	300 ml with 1 hand 350 ml with 2 hands
Max. re expansion rhythm	100 bpm

Accessori opzionali

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
BG03091A	Reservoir bag
EV50007C	EVX 07 – Overpr.valve n/rebr. autoclav. ped/neon
RM20510B	ECO B-LIFE – RUBBER FACE MASK SZ.3

Class IIa med.dev. according to Dir. 93/42/CEE

B-LIFE - BODY PAED.
 RESUSCITATOR WITH R-VALVE

BG01004B



Registration number 216123

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it