

B-Life

Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use INFANT marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Low storage space required
 Latex Free
 Steam sterilizable

Optional accessories

EV50007C EVX 07 - Overpressure and non-rebreathing valve

Optionals

OX20131A OXX30 anti crush oxygen tubes 10pcs
 BG02001A Laxan transparent bag for manual resuscitator

Technical data

Materials	Silicone, PC
Bag Volume	280 ml
Max. re expansion rhythm	95 bpm

MD CE marked according to Dir. 93/42/CEE

B-LIFE - BODY INFANT
 RESUSCITATOR WITH R-VALVE
BG01005B



Registration number 216125

Rev.2 18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it