

Eco B-Life

Pediatric disposable manual resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The material of construction and the surface finish increase grip to the operator's hands
 Connection standard ISO 15/22 mm
 Low storage space required

Technical data

Materials	Rubber, PC
Bag volume	500 ml
Stroke volume	250 ml
Max. re expansion rhythm	100 bpm

Optional accessories

BG02001A	Laxan transparent bag for manual resuscitator
----------	---

Class IIa med. dev. according to dir. 93/42/CEE

ECO B-LIFE - BODY PAED.
 RESUSCITATOR W/O R-VALVE
BG01027B



Registration number	216135
---------------------	--------

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it