

## Eco B-Life

### Infant disposable manual resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands  
 Connection standard ISO 15/22 mm  
 Low storage space required

#### Technical data

Materials	Rubber, PC
Bag volume	240 ml
Stroke volume	130 ml
Max. re expansion rhythm	95 bpm

#### Optional accessories

BG02001A	Laxan transparent bag for manual resuscitator
----------	---

Class IIa med. dev. according to dir. 93/42/CEE

ECO B-LIFE - BODY INFANT  
 RESUSCITATOR W/O R-VALVE  
**BG01028B**



Registration number	216136
---------------------	--------

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)