

LAXAN

Carry bag for manual resuscitators



Semi-transparent bag for transport of devices used in manual resuscitation procedures.

Specific features

The non-permeable material offers great protection to the content
 The integrated handle allows easy carrying
 The wide top opening with zip closure, allows quick access to the content

Technical data

Bag dimensions	295x230x140 ± 5 mm
Empty bag weight	155 ± 10 g
Material	PVC

LAXAN - TRANSPARENT BAG FOR
 MANUAL RESUSCITATOR
 BG02001A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it