

# PRO B-Life

## PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

### Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply  
 Low storage space required  
 Latex Free

### Technical data

Materials	Medical grade PVC, PC
Bag Volume	1600 ml
Stroke volume	900 ml
Overpressure valve	60 cmH <sub>2</sub> O
Max. re expansion rhythm	45 bpm

### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
BG03091A	Reservoir bag
RM10894A	Cuffed mask Size 5

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC ADULT  
 RESUSCITATOR WITH R-VALVE  
**BG02006A**



Registration number 28803

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)