

## PRO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

### Specific features

The material of construction and the surface finish increase grip to the operator's hands  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply  
 Provided with tube for oxygen enrichment  
 Low storage space required  
 Latex Free

### Technical data

Materials	PVC, PC
Bag Volume	500 ml
Stroke volume	250 ml
Overpressure valve	45 cmH <sub>2</sub> O
Max. re expansion rhythm	100 bpm
Reservoir volume	2500 ml
Mask	Disposable sz. 3

### Optionals

BG02001A	Laxan transparent bag for manual resuscitator
----------	---

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE PVC PAED. RESUS BAG  
 OV.VALVE MASK SZ.3,  
 BG02015A



Rev.0 (21/02/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)