

**PRO B-LIFE**

Child resuscitator



The resuscitators PRO B-LIFE are made of PVC and for this reason they can be stored in a compressed form and return to their normal shape as soon as needed for use. The shape of the bag also ensures great performance of flexibility.

The connection diameters are ISO 15/22 mm.

The kit includes:

- R-valve
- Overpressure valve
- Reservoir
- Anti-crush oxygen tube
- Face Mask ECO B-LIFE size 3

**TECHNICAL DATA**

Bag Volume	500 ml
Reservoir Volume	2500 ml
Vol. max per cycle	250 ml
Max cycling rate	100 bpm
Overpressure valve	40±5 mbar
Valve dead space	7 ml
Bag	PVC
Reservoir	PVC
Mask	Rubber
Non-rebreathing valve	Polycarbonate
Patient valve	Polycarbonate
Anti-crush tube	2m PVC

BG02015A

PRO B-LIFE -PVC.PED.BAG  
W/OVERP. MASK ANTIC.TUBE

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>