

PRO B-Life

PVC Resuscitator



Disposable self-inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

- Connection standard ISO 15/22 mm
- Possibility to connect the bag to sources of oxygen supply
- Low storage space required
- Latex Free

Technical data

Materials	PVC, PC
Bag volume	500 ml
Maximum oxygen concentration <small>(With reservoir)</small>	98%
Max. re expansion rhythm	100 bpm

Standard equipment

EV50010C	EVX 10 DISPOSABLE - NON-REBREATHING VALVE
----------	---

Optional accessories

OX20131A	OXX30 anti-crush tube 10pz
BG03091A	Reservoir bag
BG02000A	B-HOLDER - manual resuscitator holder
BG02001A	LAXAN – transp. bag for man.res.

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PAEDIATRIC
RESUSCITATOR BAG IN PVC
BG02016A



Rev.0 (08/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it