

## PRO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

### Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply  
 Low storage space required

### Technical data

Materials	Medical grade PVC, PC
Bag Volume	280 ml
Max. re expansion rhythm	95 bpm

### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
BG03091A	Reservoir bag

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC INFANT  
 RESUSCIT. R-VALVE  
**BG02018A**



Registration number 58603

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)