

PRO B-Life

PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required

Technical data

Materials	PVC, PC
Bag Volume	280 ml
Overpressure valve	40 cmH ₂ O
Max. re expansion rhythm	95 bpm

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
BG03091A	Reservoir bag

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC INFANT
 RESUSCIT. W/OVERP. R-VALVE
BG02019A



Registration number 58623

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it