

PRO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use ADULT marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Latex Free

Standard equipment

BG03091A	Sacca Reservoir
EV50009C	EVX09 Overpressure/non-rebreathing valve
RM10894A	Cuffed mask Size 5

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02000A	B-HOLDER - manual resuscitator holder
BG02001A	LAXAN – transp. bag for man.res.

Technical data

Materials	Medical grade PVC, PC
Bag Volume	1600 ml
Reservoir volume	2500 ml
Overpressure valve	60 cmH ₂ O
Maximum oxygen concentration	95%
Max. re expansion rhythm	45 bpm

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC ADULT RES.
 W/OVERP.RESE.CUFF.MASK

BG02044A



Registration number 16983

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it