

PRO B-Life PVC Resuscitator



Disposable self-inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
Possibility to connect the bag to sources of oxygen supply
Low storage space required
Latex Free

Standard equipment

BG03091A	Reservoir bag
EV50008C	EVX 08 - PED/NEON DISPOSABLE OVERPR.VALVE N/REBR.
RM10892A	Cuffed Mask Sz.3

Optionals

OX20131A	OXX30 anti-crush tube 10pz
BG02000A	B-HOLDER - manual resuscitator holder
BG02001A	LAXAN – transp. bag for man.res.

Technical data

Materials	PVC, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration ^(With reservoir)	98%
Max. re expansion rhythm	100 bpm

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC PAED. RES.
W/OVERP.RESE.CUFF.MASK

BG02045A

Registration number



37003

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it