

PRO B-Life

PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The material of construction and the surface finish increase grip to the operator's hands
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Provided with tube for oxygen enrichment
 Low storage space required
 Latex Free

Technical data

Materials	PVC, Rubber, PC
Bag Volume	1600 ml
Stroke volume	900 ml
Overpressure valve	60 cmH ₂ O
Max. re expansion rhythm	45 bpm
Reservoir volume	2500 ml
Mask	Disposable sz. 5

Optionals

BG02001A	Laxan transparent bag for manual resuscitator
----------	---

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC ADULT RES.
 W/OVERP.MASK AND TUBE
BG02060A



Registration number 36963

Rev.1 (19/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it