

D - PRO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
Possibility to connect the bag to sources of oxygen supply.
Disposable
Suitable for patient from 33 kg approx.

Technical data

Materials	PVC
Bag Volume	1630 ml
Stroke volume	800 ml
Overpressure valve	60 cmH ₂ O
Inspiratory and expiratory resistance at 50l/min	< 5 cmH ₂ O
Reservoir volume	2000 ml
Weight	600 g

Included equipment

Reservoir
Overpressure and non-rebreathing valve
Oxygen tubing
Mask size 4

MD CE marked according to Dir. 93/42/CEE

D-PRO B-LIFE ADULT
RESUSCITATION SET PVC
BG02100



Distributed by Spencer Italia s.r.l.

Rev.0 (12/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it