

D - PRO B-Life

Infant PVC resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply.
 Disposable
 Suitable for patient up to 10Kg approx

Technical data

Materials	PVC
Bag Volume	340 ml
Stroke volume	150 ml
Overpressure valve	40 cmH ₂ O
Inspiratory and expiratory resistance at 50l/min	< 5 cmH ₂ O
Reservoir volume	1600 ml
Weight	350 g

Included equipment

Reservoir
 Overpressure and non-rebreathing valve
 Oxygen tubing
 Mask sz 4

MD CE marked according to Dir. 93/42/CEE

D - PRO B-LIFE PALLONE RIANIM.
 PVC NEONATALE
BG02102



Distributed by Spencer Italia s.r.l.

Rev.0 (19/04/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it