

## B-Life

### Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands  
 Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply  
 Low storage space required  
 Steam sterilizable

#### Standard equipment

EV50007C | EVX 07 - Overpressure and non-rebreathing valve

#### Optionals

BG03091A	Reservoir
RM20403B	B-LIFE – Silicone face Mask sz.3
OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

#### Technical data

Materials	Silicone, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Stroke volume range	300 ml with 1 hand 350 ml with 2 hands
Overpressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration <sup>(with reservoir)</sup>	98%
Max. re expansion rhythm	100 bpm

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE PEDIATR.  
 W/VALVE SOVR.R-VALVE

**BG03023A**

Registration number



21723

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)