

B-Life

Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use INFANT marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Comes with non-rebreathing valve
 Low storage space required
 Latex Free
 Steam sterilizable

Technical data

Materials	Silicone, PC
Bag Volume	280 ml
Maximum oxygen concentration	95%
Max. re expansion rhythm	95 bpm

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
FV50007C	EVX 07 - Overpressure and non-rebreathing

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE INFANT
 RESUSCITATOR
 BG03024A



Registration number	58647
---------------------	-------

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it