

B-Life

Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

- The material of construction and the surface finish increase grip to the operator's hands
- Intended use ADULT and REUSABLE marked on the device to limit risk of incorrect applications
- Connection standard ISO 15/22 mm
- Possibility to connect the bag to sources of oxygen supply
- Low storage space required
- Latex Free

Technical data

Materials	Silicone, medical grade PVC, PC
Bag volume	1600 ml
Reservoir	2500 ml
Stroke volume range	770 ml with 1 hand 900 ml with 2 hands
Overpressure valve	60 cmH ₂ O
Maximum oxygen concentration	95%
Max. re expansion rhythm	45 bpm
Overall weight	485 ±10 g

Standard equipment

BG03091A	Reservoir
EV50011C	EVX 11 Overpressure and non-rebreathing valve
RM20405B	B-LIFE face mask Sz. 5

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE ADULT
W/OVERP. RES. R-VALVE AND
MASK SIZE 5



BG03027A

Registration number 23763

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it