

B-Life

Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Features

The material of construction and the surface finish increase grip to the operator's hands
 Intended use ADULT and REUSABLE marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Latex Free

Technical data

Materials	Silicone, PC, PVC
Bag volume	1600 ml
Reservoir volume	2500 ml
Overpressure valve	60 cmH ₂ O
Maximum oxygen concentration ^(w/ reservoir)	95%
Max. re expansion rhythm	45 bpm

Optional accessories

RM10894A	Cuffed mask Mis. 5
OX20131A	OXX30 anti crush oxygen tubes 10pz
BG02001A	Laxan transparent bag for manual resuscitator

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE ADULT W/OVERP.
 RES. AND R-VALVE

BG03029A



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it