

B-Life

Pediatric manual resuscitation unit



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The material of construction and the surface finish increase grip to the operator's hands
Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications
Connection standard ISO 15/22 mm
Possibility to connect the bag to sources of oxygen supply Autoclavable <small>(except reservoir)</small>
Low storage space required

Standard equipment

BG03091A	Reservoir
EV50007C	EVX 07 - Overpressure and non-rebreathing valve
RM20403B	B-LIFE – Silicone face Mask sz.3

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Technical data

Materials	Silicone, Medical grade PVC, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Stroke volume range	300 ml with 1 hand 350 ml with 2 hands
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration <small>(with reservoir)</small>	98%
Max. re expansion rhythm	100 bpm
Overall weight	290 ±10 g

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE
PAED.W/OVERP.RESERV. R-VALVE
AND MASK SIZE 3



BG03033B

Registration number 21706

Rev.3 (29/09/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it