

## B-Life

### Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands  
 Intended use ADULT and REUSABLE marked on the device to limit risk of incorrect applications  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply  
 Latex Free

#### Technical data

Materials	Silicone, medical PVC, PC
Bag volume	1600 ml
Reservoir volume	2500 ml
Stroke volume range	770 ml with 1 hand 900 ml with 2 hands
Overpressure valve	60 cmH <sub>2</sub> O
Maximum oxygen concentration (with reservoir)	95%
Max. re expansion rhythm	45 bpm
Overall weight	510 ±10 g

#### Included equipment

BG03091A	Reservoir
EV50011C	EVX 11 Overpr.e and non-rebreathing valve
RM20405A	B-LIFE face mask Sz. 4

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILIC. ADULT  
 W/OVERP.RES. AND MASK SIZE 4 +  
 R-VALVE  
 BG03041A



#### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.0 (30/01/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)