

## B-Life

### Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands  
 Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply.  
 Low storage space required  
 Steam sterilizable

#### Included equipment

BG03091A	Reservoir
EV50007C	EVX 07 – Overp. and non-rebreathing valve
RM20403B	B-LIFE – Silicone face Mask sz.3
OX20130A	OXX30 anti crush oxygen tube

#### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

#### Technical data

Materials	Silicone, medical PVC, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Stroke volume range	300 ml with 1 hand 350 ml with 2 hands
Overpressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration <sup>(with reservoir)</sup>	98%
Max. re expansion rhythm	100 bpm
Overall weight	320 ±10 g

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SIL.RESUSCITATOR PAED.  
 W/OVERPR. RES. MASK OXYGEN  
 TUBE R-VALVE



**BG03043B**

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)