

## B-Life

### Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands  
 Intended use INFANT marked on the device to limit risk of incorrect applications  
 Connection standard ISO 15/22 mm  
 Low storage space required  
 Latex Free  
 Steam sterilizable

#### Standard equipment

EV50007C EVX 07 - Overpressure and non-rebreathing valve

#### Optionals

OX20131A OXX30 anti crush oxygen tubes 10pcs  
 BG02001A Laxan transparent bag for manual resuscitator

#### Technical data

Materials	Silicone, PC
Bag Volume	280 ml
Overpressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration	95%
Max. re expansion rhythm	95 bpm

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE INFANT  
 W/OVERP. VALVE R-VALVE  
 BG03025A



Registration number 34003

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)