

B-Life

Neonatal manual resuscitation unit



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The material of construction and the surface finish increase grip to the operator's hands .
 Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Latex free

Technical data

Materials	Silicone, medical grade PVC, rubber , PC
Bag volume	280 ml
Reservoir volume	600 ml
Over pressure valve	40 cmH ₂ O
Maximum oxygen concentration ^(with reservoir)	95%
Max. re expansion rhythm	95 bpm
Overall weight	280±30 g

Standard equipment

BG03092A	Reservoir
EV50007C	EVX 07 – Overpressure and non-rebreathing valve
RM20504A	ECO B-LIFE – Rubber face mask sz.0
RM20503A	ECO B-LIFE – Rubber face mask sz.1
RM20502A	ECO B-LIFE – Rubber face mask sz.2

MD CE marked according to Dir. 93/42/CEE

B-LIFE SILICON INF.
 RESUSCITATION BAG WITH
 OVERPRESSURE VALVE AND MASK
 0/1/2



BG03051A

Identifiant d'enregistrement 34023

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.1 (22/07/2017)

COPIE NON CONTROLÉE – toute révision sera disponible depuis le lien <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it