

## Reservoir bag



Reservoir bag to be used for oxygen enrichment with manual resuscitators.

### Specific Characteristics

Light and thin  
 Suitable for adult and pediatric uses  
 Disposable  
 Latex Free

### Technical data

Material	PVC
Reservoir volume	2500 ml

Class IIa med. dev. according to dir. 93/42/CEE

B-LIFE SILICON ADULT/PAEDIAT  
 RESERVOIR BAG

BG03091A



Registration number	16983
---------------------	-------

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)