

Reservoir bag



Reservoir bag to be used for oxygen enrichment with manual resuscitators.

Specific Characteristics

Light and thin
 Suitable for adult and pediatric uses
 Disposable
 Latex Free

Technical data

Material	PVC
Reservoir volume	2500 ml

Class IIa med. dev. according to dir. 93/42/CEE

B-LIFE SILICON ADULT/PAEDIAT
 RESERVOIR BAG

BG03091A



Registration number	16983
---------------------	-------

Rev.1 (06/10/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it