

## Reservoir bag



Reservoir bag to be used for oxygen enrichment with manual resuscitators.

### Specific Characteristics

Light and thin  
 Suitable for infant use  
 Disposable  
 Latex Free

### Technical data

Material	PVC
Reservoir volume	600 ml

Class IIa med. dev. according to dir. 93/42/CEE

RESERVOIR BAG R-VALVE INFANT  
 BG03092A



Registration number	308067
---------------------	--------

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)