

BAG-B

Resuscitation kit



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties, supplied with carry bag.

Specific features

The material of construction and the surface finish of the bag increase grip to the operator's hands
 ADULT and REUSABLE marked on the device to minimize risks of improper applications
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply.
 Low storage space required
 Latex Free

Technical data

Carry bag dimensions	295x230x140 ± 5 mm
Materials	Silicone, PVC, PC
Bag volume	1600 ml
Stroke volume	770 ml with 1 hand 900 ml with 2 hands
Max. re expansion rhythm	45 bpm
Overpressure valve	60cmH ₂ O
Mask	Size 5

Optional accessories

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG03091A	Reservoir bag

Ass. Med.dev. according to dir. 93/42/CEE

BAG B - REANIMATION KIT
BG03100A



Registration number 59063

Rev.0 (22/09/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it