

D-B-Life

Adult reusable resuscitator



The image can differ from the actual product

Reusable self inflating manual resuscitator set for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply.
 Steam sterilizable
 Suitable for patient starting from 33 kg approx.

Included equipment

Reservoir
 Overp. and non-rebreathing valve
 Adult mask (sz4)
 Anti crush oxygen tube

Technical data

Materials	Silicone, medical PVC, PC
Bag volume	1630 ml
Reservoir volume	2000 ml
Stroke volume range	800 ml
Overpressure valve	60 cmH ₂ O
Dead space	7ml
Inspirator and expiratory resistance at 50l/min	< 5 cmH ₂ O
Overall weight	600 g

Class IIa med. Dev. according to Dir. 93/42/CEE

D - B-LIFE ADULT MANUAL
 RESUSCITATOR SET
 BG03200



Distributed by Spencer Italia s.r.l.

Rev.0 (29/03/2021)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it