

D-B-Life

Pediatric reusable resuscitator



The image can differ from the actual product

Reusable self inflating manual resuscitator set for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply.
 Steam sterilizable
 Suitable for patient from 10 to 33 kg approx.

Included equipment

Reservoir
 Overp. and non-rebreathing valve
 Pediatric mask (sz2)
 Anti crush oxygen tube

Technical data

| | |
|---|------------------------|
| Materials | Silicone, PVC, PC |
| Bag volume | 680 ml |
| Reservoir volume | 1600 ml |
| Stroke volume range | 350 ml |
| Overpressure valve | 40 cmH ₂ O |
| Dead space | 7ml |
| Inspirator and expiratory resistance at 50l/min | < 5 cmH ₂ O |
| Overall weight | 410 g |

Class IIa med. Dev. according to Dir. 93/42/CEE

D - B-LIFE PEDIATRIC MANUAL
 RESUSCITATOR SET
 BG03201



Distributed by Spencer Italia s.r.l.

Rev.0 (22/02/2021)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it